# Intelligent Evaporation

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<b>About this Document</b>	<b>58</b>
Version and variants .....	58
About this manual .....	58
Symbols and keywords.....	59
<b>Safety Instructions</b>	<b>60</b>
General safety instructions.....	60
EU Declaration of Conformity.....	60
Intended use.....	60
Installation / electrical safety .....	61
Qualifications of employees.....	61
Obligations of the operator .....	61
Installation site .....	61
Modifications to the device.....	62
Employee safety .....	62
California Residents .....	62
Safety during use .....	62
Disposal .....	63
<b>Device Description</b>	<b>64</b>
Device overview .....	64
Hei-VAP Core: Overall view .....	64
Hei-VAP Core: Control panel (manual lift) .....	65
Hei-VAP Core: Control panel (motor lift) .....	65
<b>Start-Up</b>	<b>66</b>
Connecting/disconnecting the power cord .....	66
Transportsicherung entfernen.....	66
Switching the device on and off.....	66
Setting up the device .....	67
Heating bath and flask lift.....	67
Immersing flask lift.....	68
Filling the heating bath.....	68
Connecting the heating bath .....	69
Emptying the heating bath.....	69
Menu navigation.....	71
<b>Operation</b>	<b>72</b>
Adjusting the rotation speed .....	72
Adjusting the heating bath temperature.....	73
Protecting setpoints .....	74
Evaporating and receiving flasks.....	75
Assembling and removing the evaporating flask .....	75
Assembling and removing the receiving flask.....	76
Ventilation .....	76

<b>Cleaning and Maintenance</b>	<b>77</b>
Cleaning.....	77
Maintenance .....	78
<b>Troubleshooting</b>	<b>79</b>
Malfunctions and troubleshooting.....	79
Resetting temperature limiter.....	80
<b>Assembly</b>	<b>81</b>
Electrical connections.....	81
Installation .....	82
Assembly parts.....	82
Transportation lock .....	83
Heating bath .....	84
Vapor tube.....	85
Assembling the condenser .....	87
Evaporating flask: Angle and immersion depth.....	91
Control panel .....	93
Connecting peripheral devices .....	94
Connecting cooling water.....	94
Connecting vacuum.....	95
Optional accessories .....	97
Safety hood .....	97
Safety shield .....	98
Inlet valve .....	98
<b>Dismantling and Storage</b>	<b>100</b>
Dismantling, transportation and storage .....	100
Dismantling .....	100
Transportation and storage .....	101
<b>Accessories and Spare Parts</b>	<b>102</b>
Scope of delivery.....	102
Accessories.....	102
Solvent data .....	103
<b>Attachments</b>	<b>104</b>
Technical data.....	104
<b>Service</b>	<b>105</b>
Contact / Technical Service .....	105
Warranty.....	106
Certificate of decontamination.....	107
EU Declaration of Conformity .....	108
UKCA Declaration of Conformity .....	109
RoHS Declaration of Conformity.....	110
China RoHS Declaration of Conformity .....	111

### Copyright

This operating manual is protected by copyright and is exclusively intended for the internal use of the purchaser.

Handing over this operating manual to third parties, reproducing it in any type or form – in part or in full – and utilizing and/or communicating its contents is prohibited without prior written approval from the manufacturer.

Any violations will result in an obligation to provide compensation for damages. Other claims remain unaffected.

## ➤ Version and variants

### Version

This manual describes the function, operation and maintenance of the Hei-VAP Core evaporator.

Version	Alteration date
1.0	11/2022

### Variants

The device is available in different variants. Certain features or functions are only available in specific product variants. These variants are described in this operating manual.

## ➤ About this manual

This operating manual is an integral part of the device described here.

- Carefully read through this operating manual and observe all of the safety instructions and warnings notices.
- Ensure that every user has carefully read the operating manual before using the device for the first time.
- Store this operating manual so that it is accessible to all users at all times.
- Pass the operating manual on to any subsequent owner.

The current version of this operating manual can also be found at:  
[https://heidolph-instruments.com/en/service/downloads/operation-manuals#Rotary Evaporators \(Hei-VAP Core\)](https://heidolph-instruments.com/en/service/downloads/operation-manuals#Rotary-Evaporators-(Hei-VAP-Core))

In addition, also observe the local regulations at the installation site.

## ➤ Symbols and keywords

Standardized symbols and signal words are used in this operating manual to warn against any dangers and provide important instructions. These instructions must be strictly observed to avoid accidents and damage.

The following terms and basic symbols are used:

Symbol	Additional signal words / explanation
<p><b>Warning signs</b></p> 	<p>The yellow triangle indicates hazardous situations. It is used in combination with the following signal words:</p> <p><b>DANGER:</b> Indicates a hazardous situation which, if not avoided, will result in serious injury or death.</p> <p><b>WARNING:</b> Indicates a hazardous situation which, if not avoided, may result in serious injury or death.</p> <p><b>CAUTION:</b> Indicates a hazardous situation which, if not avoided, may result in property damage and minor or moderate injury.</p>
<p><b>Prohibitory signs</b></p> 	<p>Prohibited:</p> <p>The red circle indicates a situation that should be avoided under all circumstances and which, if not avoided, may result in serious injury or death.</p>
<p><b>Mandatory signs</b></p> 	<p>Must be observed:</p> <p>The blue circle indicates important information. Please observe this information to avoid any property damage.</p>

### Other symbols used:

Symbol	Description
→	Handling instruction, action required
✓	Result of action
▪	List of information
▪	
▪	
a.	List of variants
b.	
c.	

### ➤ General safety instructions

#### EU Declaration of Conformity

This device complies with the following EC-Directives\*:

- 2006/42/EC Machinery Directive
- 2014/30/EU Electromagnetic Compatibility Directive

This device was tested according to the following Directives:

- UL 61010-1 :2012/R:2016-04  
CAN/CSA-C22.2 No. 61010-1:2012/U2:2016-04
- UL 61010-2-010:2015  
CAN/CSA-C22.2 NO. 61010-2-010:2015

The device has been manufactured according to state-of-the-art technology and in compliance with recognized safety regulations. However, risks may still arise during installation, operation and maintenance.

→ Please ensure that the operating manual is available at all times.

#### The device may only be used under the following conditions:

- Only operate the device if it is in a technically perfect condition.
- Only operate the device if you are properly aware of the risks and required safety measures.
- Only operate the device in accordance with the instructions given in this operating manual.
- If any information is ambiguous or missing, ask your superior or contact the manufacturer.
- Do not operate anything on the device without authorization.
- Only use the device in accordance with its intended use.

#### Intended use

The device may only be operated by authorized personnel.

The device is suitable for the following uses:

- Distillation to separate substances
- Concentration
- Crystallization
- Powder drying

The device may be operated in research laboratories, other laboratories and production facilities in the following sectors:

- Chemistry
- Pharmacy
- Biology
- Environmental analytics
- Basic research
- Similar research laboratories

#### Unintended use/Reasonable foreseeable misuse

Any use which deviates from the device's intended use is considered to be improper. The manufacturer is not liable for damage that occurs as a result. The risk is borne by the operator alone.

When used in medical and food sectors the operator alone is responsible for compliance with the applicable regulations and standards.

### Installation / electrical safety

- The device may only be connected up if the available power supply voltage corresponds to the information stated on the rating plate for the device.
- The power supply connection must be easily accessible at all times.
- An RCD circuit breaker (residual current circuit breaker) must be fitted to the electrical system in the building so that the device is disconnected from the power circuit in the event of a fault.
- Repairs may only be performed by electricians authorised by Heidolph Instruments.
- Never operate the device with a damaged power cord.
- Always turn the device OFF and disconnect the power cord before carrying out any maintenance or repair work.

### Qualifications of employees

- The device may only be operated by trained personnel.
- The device may only be operated by persons who have been instructed and supervised in its proper use by trained specialist personnel.
- The device may only be operated by specialist personnel who are above the legal minimum age.
- Other persons may only work on the device under the constant supervision of experienced and trained specialist personnel.
- This operating manual must be read and understood by all persons working with the device.
- Personnel must receive safety training that ensures responsible and safe working practices.

### ➤ Obligations of the operator

#### Installation site

- The device must be installed in a suitable location.
- The device must be installed on a firm stable and temperature resistant surface.
- Ensure that the device and all of its components are easily accessible at all times.
- Maintain a sufficient safety clearance to, in particular, any moving and/or hot equipment components.
- It is not permitted to place or store any objects such as accessories, tools or chemicals within this clearance area during operation.
- All screw connections must be securely tightened.
- It is not permitted to operate the device near to highly flammable or explosive substances.
- Operate the device in conjunction with an extractor hood, if working with potential harmful media (see DIN EN 14175 and DIN 12924).

## Safety Instructions

- The device is designed for indoor use only and under the following ambient conditions:

Ambient temperature	5 - 31 °C at 80 % relative humidity 32 - 40 °C decreasing linearly to a maximum 50 % relative humidity
Installation altitude	0 - 2,000 m above sea level
Contamination level	2
Overvoltage category	II
Permissible supply deviations	± 10 %
- If the device is operated in corrosive atmospheres, the service life of the device will decrease based on the concentration, duration and frequency of the exposure to the corrosive atmosphere e.g. concentrated hydrochloric acid (HCl).

## Modifications to the device

- It is not permitted to make any unauthorised modifications or changes to the device.
- Do not attach or install any parts that have not been approved by the manufacturer.
- Unauthorised modifications or changes will void the EC Declaration of Conformity for the device and operation of the device will no longer be permitted.
- The manufacturer is not liable for any damage, dangers or injuries that result from unauthorised modifications and changes or due to the non-observance of the instructions in this manual.

## Employee safety

- Ensure that the device is only operated by qualified specialist personnel and trained employees.
- Observe the following instructions to avoid any personal injuries and property damage:
  - Laboratory regulations
  - Accident prevention regulations
  - Hazardous Substance Act
  - Other generally accepted occupational health and safety regulations
  - Local regulations

## California Residents

- Important information for California residents regarding Prop 65. Please visit [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov) for more information.

## ▶ Safety during use

- Wear the appropriate clothing when working on the device (protective clothing, safety glasses and, if necessary, safety gloves).
- Do not use the device in potentially explosive areas. The device is not protected against explosion. There is no explosion or ATEX protection available.
- Do not carry out work with exposed flames in the vicinity of the device (risk of explosion).
- Do not operate any devices in the vicinity which may exhibit electromagnetic fields in the frequency range  $9 \times 10^3$  Hz to  $3 \times 10^{11}$  Hz.

## Safety Instructions

- Do not operate or assemble devices in the vicinity which are emission or radiation sources (electromagnetic radiation) for the frequency range  $3 \cdot 10^{11}$  Hz to  $3 \cdot 10^{15}$  Hz (in the optical spectral range wavelengths from  $1,000 \mu\text{m}$  to  $0.1 \mu\text{m}$ ).
- Do not operate or assemble appliances in the vicinity of the device which constitute emission or radiation sources for ionizing radiation or in the ultrasonic range.
- Do not operate the device where adiabatic compression or shock waves might occur (shock wave combustion).
- Do not use materials which pose a potential risk of uncontrolled energy release which might result in pressure increase (exothermic reaction; spontaneous combustion of dusts).
- Do not work with hard, brittle materials (such as stones, soil samples, etc.) which could result in the destruction of the evaporating flask.
- Use only heating bath fluids which guarantee suitable heat transfer.
- Do not operate the device with excess pressure.
- Limit the pressure difference for glass components to a maximum of 2 bar.
- Do not increase the maximum excess pressure of the cooling fluid above 2 bar.
- Wear safety gloves and safety glasses when working with coolants in very low temperature ranges, such as with dry ice.
- Do not increase the flow speed of 1 m/s while siphoning liquids in with combustible materials (electrostatic charge; risk of ignition).
- Eliminate the production of gases classified in explosion group IIC during distillation by materials or due to chemical reactions, such as hydrogen.
- Do not use distillation material where the distillation residue might be explosive.
- Only connect clean ground joints.
- Attach ground joints securely in order to connect them leak-tight.
- Attach all tube connections securely to ensure safe operation.
- Connect all cables and tubes without kinks and locate them outside the operating and danger zone.
- Avoid putting pressure on the display when you are not operating the device.
- Do not spill liquids over the device or any parts of it.
- Remove any accidentally spilled liquids immediately.
- Eliminate errors immediately.
- Do not use abrasive material to clean the surface. Only wipe with damp cloth.
- Always switch the device OFF after use.

## Disposal

- Check the device components for hazardous substances and solvents.
- Clean all components before disposal.
- Dispose of the device in accordance with the relevant national regulations.
- Dispose of the packaging material in accordance with the appropriate national regulations.

## ➤ Device overview

### Hei-VAP Core: Overall view

## Hei-VAP Core: Control panel (manual lift)

## Hei-VAP Core: Control panel (motor lift)

### ➤ Connecting/disconnecting the power cord

#### Connecting the power cord

The power cord has a three-pole connector, and has a recess on its bottom side.


Power cord connector

Power cord

The device socket is located on the rear of the device.

- Insert the power cord connector into the device socket with the recess pointing downwards.
- ✓ The connector has been securely inserted when it has been pushed in up to the stop.
- Now insert the power cord plug into the power socket.

Device socket

#### Disconnecting the power cord

- First, remove the plug from the power socket.
- Then disconnect the power cord connector from the device.

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#### Caution:

#### Possible damage to the device (motorlift variants only)

Switching the device on before the heating bath and the control panel are connected can cause malfunctioning of the electronics the moment the transportation lock is removed.

Prior to switching on the device for the first time, remove the transportation lock exactly in the order given in the instructions (see chapter "Assembly", section "Transportation lock").

### ➤ Switching the device on and off

The ON/OFF switch is located at the lower left side of the device.

#### Switching the device on

- Press the ON/OFF switch on the left side.
- ✓ The currently set values will appear on the display.
- ✓ The device is switched on.

ON/OFF switch: OFF

### Switching the device off

- Press the ON/OFF switch on the left side.
- ✓ The display will go blank.
- ✓ The device is switched off.

## ➤ Setting up the device

### **Warning:** **Risk of slipping device!**

Due to an unsuitable surface or as a result of vibrations generated during the operation the device may tilt or start moving and slide off the table top.

The installation surface must be smooth, level, clean and temperature resistant.

Avoid shaking and unwanted impacts to the device as well as at the installation surface.

- Set up the device on a solid, smooth, level and temperature-resistant surface.
- Clean installation surface and feet of the device with a damp cloth and ethanol regularly.

## ➤ Heating bath and flask lift

The heating bath can be filled with different heat transfer fluids, e.g.,

20° C to 90° C	Tapwater
20° C to 170° C	Polyethylene glycol (z. B. Heating bath fluid Budde)*
150° C to 210° C	Silicon oil with low Viscosity (40 cP), Flammpunkt > 285 °C z. B. Wacker AP 200 Silikonöl)*

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### **Danger!** **Mortal danger! Risk of explosion!**

Heating bath fluid oil with a flash point < 285 °C could catch fire at high heating bath temperatures and lead to uncontrolled thermal reactions up to an explosion.

- Use only oil with a flash point > 285 °C.

The markings in the heating bath indicate the minimum and maximum liquid levels for the immersed evaporating flask.

The evaporating flask is moved up and down using a lift and is immersed in the heating bath in this manner.

When the lift moves downward, the evaporating flask displaces the already filled in water or oil.

To prevent the heating bath from overflowing:

- Prior to filling, move the lift with the mounted evaporating flask downward.
- Only then fill the water or oil up to at most the upper mark of the heating bath tank.

### Immersing flask lift

Immerse the piston in the heating bath only so far that it does not touch the edge or the bottom.

For adjustment of the immersion depth and adjustment of the flask angle, see Chapter "Assembly", Section "Evaporating Flask: Angle and immersion depth".

#### Immersing the flask Manual lift


- Press the hand lift lever downward.
- While pressed, move the manual lift lever to the left for lowering or to the right for raising.
- Release the hand lift lever at the desired position.
- ✓ The manual lift locks in the corresponding position.

Manual lift: Upper lift position

#### Immersing the flask Motor lift

- ✓ The device is switched on.
- Press the up ▲ or down ▼ lift button.
- Release the lift button at the desired flask position.
- ✓ The motor lift will stop.

Lift buttons: Up, Down

### Filling the heating bath

Over time deposits of lime may be observed. This may result in early development of rust.

We recommend to clean the bath regularly with a standard cleaning/polish agent. That removes lime and also prevents the development of new lime.

#### **Caution:** **Risk of corrosion**

Do not utilize deionized or distilled water. This may cause corrosion to your bath.

Utilize tap water as heating bath fluid only!

- ✓ The evaporating flask is immersed in the heating bath at the correct position.
- Fill the heating bath with the heating fluid at most up the upper marking.
- ✓ The liquid level is within the minimum and the maximum marking.

### Connecting the heating bath

**Warning:**  
**Risk of crushing! Risk of burns!**

If the device is switched on, the heating or rotation can be started unintentionally. This may injure and burn your hands.

Connect the heating bath to the device base only when the, device is switched off.

The heating bath socket is located on the right rear side of the heating bath.

- ✓ The device is switched off.
- ✓ The heating bath is locked in place on the baseplate.

Heating bath connector

Heating bath socket


- Insert the 7-pin heating bath connector with the recess pointing upward into the heating bath socket.
- Tighten the knurled collar nut clockwise to the right.

Collar nut

Heating bath socket

Recess

Key

### Emptying the heating bath

**Warning:**  
**Risk of scalding! Risk of slipping!**

When emptying the heating bath, you may burn yourself or slip on spilled fluid.

Before emptying heating bath:

- Please wear protective clothes.
- Switch the device OFF.
- Wait until the heating bath fluid has cooled down to <math>< 50^{\circ}\text{C}</math>.
- Disconnect base device and heating bath.
- Remove evaporating flask.
- Lift the heating bath only by the insulated handle.

## Start-Up

- ✓ The device is switched OFF.
- ✓ The evaporating flask is removed.
- Loosen the knurled collar nut anti clockwise until it is free.
- Pull the heating bath connector from the heating bath.
- Place your thumbs on the corrugated surface of the handles and then completely grasp them.
- Cautiously lift the heating bath at both handles and empty it.

Protective lip

Handles

- ✓ The evaporating flask has been removed.
- Turn the union nut at the heating bath connecting cable to the left anticlockwise direction, until it is released.
- Disconnect the connector socket contact from the heating bath.
- Carefully lift the heating bath at both handles and empty it.

### ➤ Menu navigation

The initial screen displays the current actual values. Use the rotary knobs to set the parameters Rotation and Temperature for distillation.

- ✓ The device is switched on.
- ➔ Turn the rotation knob or the temperature knob to set the desired value.
- ✓ The setpoints are shown in the display and designated as "set".
- ✓ If no further action is taken, the display changes after few seconds and shows the current actual values again.
- ✓ Pressing the respective rotary knobs will start the heating and the rotation.

- ✓ When rotation is switched on, the rotation knob lights up with a white illuminated ring.
- ✓ With heating switched on, the temperature knob lights up with an orange illuminated ring.
- ✓ The heating function is indicated by a heating symbol next to the temperature display.

Rotation switched on    Heating switched on

- ✓ Press the rotary knobs again to switch off the heating or the rotation.
- ➔ Press the setpoint button "show set" to show the setpoints in the display. The setpoint display is designated with "set".
- ➔ Press the Standby button at the bottom right to switch the standby mode on and off.

#### **Warning:** Possible unwanted parameter changes

As long as the display and the rotary knobs are active, parameters may be changed unintentionally.

If the device is to be operated without surveillance only leave when the rotary knobs and the display are locked (see chapter "Operation, Protecting setpoints").

### ➤ Adjusting the rotation speed

The distillation speed can be increased by increasing the rotation speed. The shortened distillation time contributes to the thermal protection of the distillation material.

However, depending on the size of the evaporating flask and the selected heating fluid, different rotation speeds are appropriate.

**Warning:**  
**Risk of injury! Risk of entanglement!**

When operating the device without a safety hood, you may get in contact with rotating parts at the evaporating flask.

Rotating items can injure your hands or other parts of the body by drawing in long hair, loose clothing or jewelry.

- Never touch any rotating parts.
- Do not wear loose clothing, belts or ties, which may be drawn in by rotating parts.
- Do not wear long necklaces or bracelets.
- Cover or tie back long hair when operating the stirrer.

### Rotation speed range and adjustment

The rotation speed can be adjusted between 10 rpm and 280 rpm.

**Warning:**  
**Risk of scalding! Risk of slipping!**

At high rotation speeds, the heating bath fluid can splash and/or spill over, resulting in slippery floors and severe scalding.

If you operate the device without a safety hood:

- Always test the fluid used before a process and document the spraying behaviour!
- Make this documentation available to all your employees.
- Wear protective clothing near the device.
- Make sure the floor is clean.

- Test the fluid selected by you
- at slowly rising temperatures
  - with your different evaporating flask sizes
  - with rising rotation speeds
- Note the results and make them available to all your employees.

The viscosity and volume of the oil change depending on the oil selected, the age of the oil and the setting of the bath temperature.

For this reason, document both the exact name/manufacturer specification and the period of use for oil and repeat the documentation at specified time intervals.

## Operation

- ✓ The base device is switched on.
- ✓ The evaporating flask is immersed in the heating bath.
- Adjust the desired rotation speed via the rotation knob.
- Press the rotation knob to start the rotation.
- Briefly press the rotation knob to end the rotation.

Rotation knob: Rotation switched on

## ➤ Adjusting the heating bath temperature

### Adjustable temperature range:

The temperature range can be adjusted between 20 °C and 210 °C.

- Permissible temperature range when using water: 20 °C - 100 °C
- Permissible temperature range when using oil: 20 °C - 210 °C

Starting at a desired bath temperature of 75 °C, we recommend using polyethylene glycol or silicon oil as the heating fluid (also see chapter "Accessories and Spare parts, Accessories").

### **Danger!** **Mortal danger! Risk of explosion!**

Heating bath fluid oil with a flash point < 285 °C could catch fire at high heating bath temperatures and lead to uncontrolled thermal reactions up to an explosion.

→ Use only oil with a flash point > 285 °C.

### **Warning:** **Risk of scalding!**

If you operate your device without a safety hood, severe scalding may occur at high bath temperatures.

Be sure to wear protective clothing (safety goggles, protective clothing, protective gloves) when you are near the device while a process is running.

### **Caution:** **Possible device damage**

The heating responds very fast. If the heating bath is empty and heating is switched on, the heating bath may turn dark.

Only switch heating on when the heating bath is sufficiently filled.

## Operation

### Adjusting the heating bath temperature

- ✓ The base device is switched on.
- ✓ The heating bath is filled with oil or water.
- ✓ The evaporating flask is immersed in the heating bath.
- Adjust the desired temperature via the temperature knob.
- ✓ As soon as you set the temperature to  $\geq 100\text{ }^{\circ}\text{C}$ , "oil" is shown as an advice in the display.


Temperature knob

Water is not permitted as heating bath fluid at temperatures  $\geq 100\text{ }^{\circ}\text{C}$ . For safety reasons, "oil" is shown as an advice in the display as soon as you set the temperature to  $\geq 100\text{ }^{\circ}\text{C}$ .

- Press the temperature knob to start heating.
- ✓ As long as heating is switched on, the temperature knob lights up orange.
- ✓ The heating function is indicated by a heating symbol next to the temperature display.
- Press the temperature knob again to switch off heating.

Heating symbol

The orange light ring flashes if residual heat is still present after switching off.. The flashing frequency depends on the level of the current residual heat. The flashing stops once the residual heat drops below  $50\text{ }^{\circ}\text{C}$ .

The display of the residual heat is only available when the device is switched on and when in standby mode. The residual heat is not displayed if the device has been switched off via the main switch (ON/OFF switch).

## ▶ Protecting setpoints

You can protect both the setpoint for the rotation speed and the setpoint for the heating bath temperature from unintentional adjustment.

- Press the respective rotary knob for approx. 2 sec. to lock the setpoint.
- ✓ The locked value is indicated in the display by a frame.
- ✓ Turning the rotary knobs has no effect on the setpoints.
- Press the rotary knob again for approx. 2 sec. to unlock the setpoint.
- ✓ Turning the rotary knobs will change the setpoints.

## ➤ Evaporating and receiving flasks

### Assembling and removing the evaporating flask

**Warning:**  
**Risk of injury!**

If the device is switched on during assembly of the evaporating flask, it is possible for the device to start the rotation due to unintentionally touching of the control panel and this could result in injuries.

Only assemble and remove the evaporating flask when the device is switched off.

**Warning:**  
**Risk of glass breakage!**

If one of the flasks is damaged, glass may break, and you may be seriously injured.

Only work with flawless glass. Inspect the glass regularly for damage.

#### Assembling the evaporating flask

- ➔ Loosen the coupling of the Easy-Clip and open the swivelling locking part of the Easy-Clip.
- ➔ Now guide the evaporating flask over the ground-joint of the vapor tube.
- ➔ Close the Easy-Clip by sliding the swivelling part of the Easy-Clip over the flange at the evaporating flask so that the Easy-Clip engages.
- ✓ The evaporating flask now hangs securely in the Easy-Clip lock.
- ➔ Tighten the Easy Clip on the vapor tube.

Easy-Clip opened

Easy-Clip closed

#### Removing the evaporating flask

To remove the evaporating flask:

- Heating must be switched off
- Rotation must be switched off
- The system must be ventilated, and ambient pressure must be established
- The lift with the evaporating flask must be in highest position (see chapter "Start-Up", section "Heating bath and flask lift")
- The device must be switched off

### **Warning:** **Risk of burns!**

Contact with the evaporating flask can cause severe burns if it has not cooled down sufficiently.

- Always allow the evaporating flask to cool down before removing it.
- The maximum temperature must not exceed 50 °C / 122 °F.
- Wear protective clothing.
- Wear protective gloves to ensure a secure grip.

- Loosen the coupling of the Easy-Clip flask clamp from the vapor tube using 1-2 anticlockwise turns.
- Grip the evaporating flask with one hand and open the Easy-Clip.
- Carefully remove the evaporating flask.

## Assembling and removing the receiving flask

### **Assembling the receiving flask**

- Hold the receiving flask with the opening facing up against the flange of the condenser.
- Slide the clamp over the flange on the condenser and the receiving flask.
- Carefully tighten the knurled screw of the clamp.
- ✓ The receiving flask is attached.

### **Removing the receiving flask**

- Grip the receiving flask with one hand.
- Loosen the knurled screw of the clamp.
- Remove the clamp.

Knurled screw

Clamp

## ➤ Ventilation

The system can be ventilated manually using the ventilation cap:

- For ventilation, turn the ventilation cap half a turn or a maximum of 1 turn anticlockwise to the left.
- ✓ The ventilation speed increases with the rotation.
- Slowly turn the ventilation cap clockwise to the right.
- ✓ The air supply will slowly close.

Ventilation cap

### ➤ Cleaning

The housing and surfaces of the device can be cleaned using a damp cloth and, if required, a mild soap solution.

**Warning!**  
**Risk of electric shock!**

Liquids penetrating into the device could cause the device to become live if it is still connected to the power supply.

Before cleaning the device:

- Switch the device off.
- Disconnect the power cord plug.

**Caution:**  
**Damage to the device**

Liquids penetrating into the device could damage the electronics.

Only wipe the surfaces with a damp cloth and not a wet cloth.

**Caution:**  
**Damage to surfaces**

The surfaces could be damaged by improper cleaning.

Do not use the following under any circumstances:

- Chlorine bleach or cleaning agents containing chlorine
- Substances containing solvents (e.g. acetone)
- Ammonia
- Abrasive cleaning agents such as cleaning wool, scouring agents or cleaning agents with metallic particles

### Cleaning the heating bath

Rust spots on the bottom of the heating bath after shorter or longer use are caused by the deposition of foreign substances (iron particles or rust particles).

Rust spots can be removed by cleaning with a commercial available stainless steel polish.

Lime deposits can be removed easily by cleaning with a commercially available lime solvent, e.g. on citric or acetic acid basis.

Please follow the instructions of the cleaning agent manufacturer.

Deposits of lime can lead to premature formation of rust as a result of crevice corrosion.

We therefore recommend to clean the heating bath at regular intervals with a stainless steel polish. At the same time, certain polishes can hinder depositions of lime.

### Cleaning of coated glassware sets

If the cleaning temperature is too high, the coating of the glassware sets may be damaged.

Rinse coated glass parts only up to a maximum temperature of 50 °C.

### ➤ Maintenance

---

**Caution:**  
**Damage to the device**

When operating with a damaged PTFE seal, the desired vacuum cannot be achieved, your sample can be contaminated, and your instrument can be damaged.

Never work with a damaged or worn PTFE seal. Examine the PTFE seal regularly and replace it if necessary.

---

The device is maintenance-free. Any repairs that may be required must be performed by a specialist authorised by Heidolph Instruments.

In this case, contact your authorised Heidolph dealer or a representative of Heidolph Instruments.

### ➤ Malfunctions and troubleshooting

Errors at the device	Possible cause	Remedy
<ul style="list-style-type: none"> <li>▪ The device cannot be switched on</li> <li>▪ Display remains dark after switching on</li> </ul>	No power supply	Check that the power cord plug is connected properly to the power socket and the device
	Fuse defective	Replace device fuses (see chapter " Assembly, Electrical Connection ")
	ON/OFF switch Off	Switch ON/OFF switch on
	ON/OFF switch defective	Contact local dealer
Device does not heat	Connecting cable to heating bath is not attached	Attach connecting cable
	Bath heating defective	Contact local dealer
	Overheat protection activated	<ul style="list-style-type: none"> <li>▪ Let the heating bath cool down</li> <li>▪ Emptying the heating bath</li> <li>▪ Reset overheat protection (see the following section "Resetting the temperature limiter")</li> </ul>
Drive does not turn	Rotation is switched off	Press rotation knob
	Drive defective	Contact local dealer
Motor lift does not run	Lift is at end stop	Move lift in the opposite direction
	Lift is at height stop	Change height stop
	Mechanics / motor defective	Contact local dealer
No vacuum	ON/OFF switch at the vacuum pump off	Switch ON/OFF switch on
	Vacuum valve defective	Replace vacuum valve
Vacuum insufficient	Leak in system	<ul style="list-style-type: none"> <li>▪ Check seals and connections</li> <li>▪ Check ground glass joints, grease if necessary</li> </ul>
	Vacuum pump defective	Follow information from the vacuum pump manufacturer

## Troubleshooting

Error messages	Possible cause	Remedy
1	<p><b>Only for motor lift version:</b></p> <p>The transportation lock is still activated.</p>	<ul style="list-style-type: none"> <li>Remove transportation lock</li> <li>Move lift upward</li> <li>Switch the device off and then on again</li> </ul>
2	<p>The cable connection between the control panel and the base device gas bee separated.</p>	<ul style="list-style-type: none"> <li>Re-establish the cable connection</li> <li>Switch the device off and then on again</li> </ul>
3	<ul style="list-style-type: none"> <li>Heating defective</li> <li>Overheat protection activated</li> </ul> <p>Overtemperature protection +5 °C exceeded</p> <p>Heating bath fluid completely evaporated</p> <p>Heating bath sensor difference &gt; 10 °C</p> <ul style="list-style-type: none"> <li>Temperature sensor heating bath defective</li> <li>Temperature of the heating bath fluid &gt; 217 °C</li> <li>Temperature of the heating bath fluid &lt; 0 °C</li> </ul>	<ul style="list-style-type: none"> <li>If necessary, allow the heating bath to cool down and empty, reset the overheat protection (see "Resetting temperature limiter" section).</li> </ul> <p>Contact local dealer</p> <ul style="list-style-type: none"> <li>Switch the device off</li> <li>Allow the heating bath to cool down, reset the overheat protection (see "Resetting temperature limiter" section).</li> <li>Filling with heat transfer liquid</li> </ul> <p>Contact local dealer</p> <p>Contact local dealer</p>
5	<p>Motor overloaded</p>	<ul style="list-style-type: none"> <li>Switching the device off</li> <li>Contact local dealer</li> </ul>

## ➤ Resetting temperature limiter

- ➔ Switch the device off Remove the evaporating flask
- ➔ Disconnect the heating bath plug remove the heating bath.
- ➔ If necessary, empty the heating bath and turn it over.
- ➔ Use a suitable object to briefly press the red button on the bottom until it clicks into place.
- ✓ Reset the temperature limiter

Temperature limiter

If you experience a malfunction that cannot be resolved using the information listed above, please inform your authorized Heidolph dealer immediately.

## Electrical connections

- Fuses must **only** be replaced by a **professional electrician**.
- **Repairs** to the device must **only be carried out by a qualified electrician approved by Heidolph Instruments**. Improper repairs could result in serious dangers. Please contact your specialist dealer if the device needs to be repaired.

### Warning: Risk of electric shock

The device must be connected to an earthed power socket.



The rated voltage for the device and the supply voltage must be the same. The rated voltage for the device can be found on the rating plate on the rear of the device.

The device is supplied with a grounded power cord. In case the original plug is replaced, the new plug must have a protective conductor!

If the device is still connected to power while changing the fuse, you might get in contact with live parts.

Prior to replacing the fuse ensure that the unit is switched off and pull the plug from the electric socket. Only use original spare fuses from Heidolph Instruments.

After changing a fuse check the device for safe operating condition according to IEC 61010-1.

### Colour code for electrical connections:

#### Europe

<b>GREEN/YEL-LOW</b>	PE: Protective conductor ( <b>Earth</b> )
<b>BLUE</b>	N: Neutral conductor
<b>BROWN</b>	P: Phase

#### North America

<b>GREEN</b>	PE: Protective conductor ( <b>Earth</b> )
<b>WHITE</b>	N: Neutral conductor
<b>BLACK</b>	P: Phase

The device is connected to the mains with the power cord supplied. The power connection is located on the rear of the housing.

For countries where a different plug than the standard plug supplied is used:

- The plug supplied may only be changed by a professional electrician.
- If you wish to use an adapter, ensure that it is grounded and approved in accordance with the local regulations.
- Wiring, adapter and plug must at least meet the electrical power stated on the rating plate.

## Installation

### Assembly parts

#### Basic device

Control panel support Base plate heating bath

#### Control panel

##### Control panel

##### Easy-Clip

#### Heating bath

##### Pouring rim

##### Min./Max. marks

Insulated handles

**Warning:**  
**Risk of crushing! Risk of burns!**

If the device is switched on, the heating or rotation can be started unintentionally. This may injure and burn your hands.

Install the device only when the device is switched off.

## Transportation lock

### (Motor lift only)

The transportation lock must be removed before start-up of the device. Retain the transportation lock for a possible return.

The transportation lock consists of three M5x8 screws and a connecting metal sheet with holes for securing the position.

### Remove transportation lock prior to start-up

- ✓ The basic device is positioned on a stable, level surface.
- Connect the heating bath.
- Connect the control panel.
- Switch the device on.
- Use the provided allen key to loosen the 3 screws and remove the metal sheet of the transportation lock.
- Move the flask lift to the highest position.
- Switch the device off.
- ✓ You can now switch the device back on and put it back into operation.
- Retain all parts of the transportation lock for later use.


Transportation lock

### Heating bath

#### Insert the heating bath

- ✓ A sliding base plate with a raised guide rail is located on the right side of the base device.

Guide rail

Heating bath

- Place the heating bath with the recess to the left on the guide rail.
- Move the heating bath slightly to the right and left, so that it locks securely.

Recess

#### Moving the heating bath

When using larger evaporating flasks or adapters between evaporating flask and vapor tube, it may be necessary to increase the distance between the heating bath and the drive.

- ✓ The basic device is positioned on a stable, level surface.
- ✓ There is about 20 cm space for movement on the same base to the right of the heating bath.
- Hold the heating bath by both handles and pull it to the right to the desired position.

## Vapor tube

### Assembly parts



- Vapor tube
- Flask clamp Easy-Clip
- Coupling
- Clamping sleeve

The Easy-Clip flask clamp is pre-mounted in the NS29 size. The Easy-Clip NS24 is included for smaller flasks.

When the device is delivered, the coupling for the vapor tube and the condenser are pre-assembled.

- Press and hold the locking knob below "lock" for a quarter turn.
- Unscrew the coupling of the vapor tube together with the Easy-Clip to the right of the drive.
- Unscrew the coupling for the condenser to the left of the drive.

- Remove the tension spring and the PTFE seal on the condenser side.

### Caution: Possible contamination, leakage, damage to your device

A damaged seal can prevent reaching the desired vacuum, can contaminate your sample or damage your device.

Never work with a damaged or worn PTFE seal. Inspect the PTFE seal regularly and replace it if necessary.

## Assembly

- ✓ Coupling and Easy-Clip are screwed together.
- Push the coupling of the vapor tube together with the Easy-Clip onto the vapor tube. (Direction and sequence as shown under "Installation parts")
- Push the clamping sleeve with the bead toward the glass onto the vapor tube until it clicks into place.
- Insert the vapor tube with the clamping sleeve first from below into the drive and screw it tight.
- Press the locking knob during the last quarter turn.

Vapor tube, Easy-Clip and coupling

### Assembling the condenser

**Warning:**  
**Risk of glass breakage!**

If a glass component is damaged, glass may break, and you may be seriously injured.

Only work with flawless glass. Inspect the glass regularly for damage.

The coated glassware sets from Heidolph Instruments reduce the risk of injury in the event of glass breakage. However, they are only suitable for a temperature range up to max. 80 °C.

#### PTFE seal

- Slide the PTFE seal with the mark "Motor Side" first forward onto the vapor tube.
- ✓ The PTFE seal has been pushed to the stop.

PTFE seal

#### Condenser:

##### Glassware set G1, G3 and G3 XL

- Slide the coupling for the condenser over the cooler flange.
- Slide the tension spring over the cooler flange.
- Hand-tighten the cooler with the coupling on the drive.

Tension spring      Coupling

- ✓ The cooler is mounted on the drive.

### Condenser:

#### Glassware set G5 Dry ice cooler

- ✓ The cooler is mounted on the drive.
- Screw the ventilation cap onto the opening opposite the vapor tube.
- If you are working with inert gas, screw the inlet tube with inert gas valve onto the opening above the ventilation cap.

Inert gas valve

Ventilation cap


- Place the red sealing ring in the groove at the top of the cooler flange.
- Place the white PTFE ring above it with the collar facing upwards in the cooler flange.
- Now insert the jacket trap into the condenser.
- ✓ The flange is drawn in by the vacuum and the system is closed vacuum-tight by the installed seal.
- Place the black lid on top to protect your dry ice.

**Condenser: Glassware set G6  
(Cooler for reflux distillation)  
Additional assembly parts at the  
vapor tube:**


- PTFE seal ring
- Spacer sleeve
- Additional coupling


- ✓ The vapor. tube is installed together with the PTFE seal.
- Install the PTFE seal ring in the spacer sleeve.
- Install the spacer sleeve together with the PTFE seal ring with the narrow side first in the vapor tube on the condenser side.
- Screw the additional coupling over it with the male thread pointing left.
- Tighten the distributor with the coupling including the tension spring onto the male thread.
- ✓ The distributor is installed.

**Distributor**


- Insert the condenser into the upper ground glass joint on the distributor.
- Screw the ventilation cap onto the opening opposite the vapor tube.

## Assembly

- Turn reflux valve clockwise until the seal ring is clearly visible and white plug blocks the output to the receiving flask.
- ✓ The valve is closed, and the sample is guided back to the evaporating flask.
- Turn reflux valve anticlockwise until the seal ring is clearly visible and white plug unblocks the output to the receiving flask.
- ✓ The valve is open, and the sample is guided into the receiving flask.

The glassware set is no longer vacuum-tight if the reflux valve is opened too far. Open the valve only so far that the white seal ring underneath the knurled screw remains visible about 0.5 cm wide.


### Support system for vertical condenser. Glassware set G1, G3 XL and G6

#### Assembly parts:

- Support rod
- Cross clamp with condenser clamp

A hole is located on the left side of the device for receiving the support rod.

- ✓ First mount the condenser and then slide the condenser clamp over the condenser.
- Insert the support rod with the flattened side downward into the hole.
- Secure the support rod with the knurled screw.

## Assembly

- Slide the cross clamp over the support rod.
- Use the knurled screw to secure the cross clamp on the support rod.
- Secure the condenser clamp with the knurled screw around the condenser.

Knurled screw cross clamp

### Support system for vertical condenser. Glassware set G5

#### Assembly parts:

- Support rod
  - Cross clamp with tension band
- Mount the support system analog to the support for G3, G3 XL and G6.
  - Use the tension band to secure the condenser.

Tension band

Knurled screw

Cross clamp

## Evaporating flask: Angle and immersion depth

### Flask angle

The evaporating flask angle can be adjusted by turning the drive.

- ✓ The rotation is switched off.
  - ✓ The evaporating flask is mounted (see chapter "Operation", section "Mounting the evaporating flask").
- Slide the locking mechanism on lift column to the right and hold it against the spring force.
  - Pivot the condenser carefully until the connected evaporating flask has the correct angle. Ensure that the flask can rotate freely and does not touch anything.
  - Release the locking mechanism.
  - If necessary, pivot the condenser slightly until it locks in place.

Lock

Lift column

### Flask immersion depth (Height stop)

The stop band of the height stop is located on the right side of the lift column that can be used to adjust a constant immersion depth for the evaporating flask.

- ✓ The rotation is switched off.
- ✓ The evaporating flask is mounted (see chapter "Operation", section "Mounting the evaporating flask").
- ✓ The height stop is at the top position.

#### Adjusting the lift stop without the safety hood:

- Lower the lift until the desired immersion depth is reached (see chapter "Start-Up" section "Heating bath and flask lift").
- Turn the locking screw on the stop band of the height stop to the right of the lift column by about 180° anticlockwise
- Slide the stop band down as far as it will go using the locking screw.
- Tighten the locking screw clockwise.
- Return the lift to the top.
- ✓ The locking screw holds the stop in position.
- ✓ The lift moves down no further than the stop and ensures a constant immersion depth of the evaporating flask

#### Adjusting the lift stop with the safety hood:

- Hold the safety hood with one hand.
- Guide the safety hood by hand when moving the safety hood (see chapter "Assembly, Safety hood").

### Control panel

The control panel support is located in front of the flask drive.

Control panel support

**Warning:**  
**Risk of poisoning!**


If you operate the device inside a fume hood and have to reach the control panel during a running process, you may be endangering your health when you open the fume hood.

Remove the control panel and position it outside the fume hood.

**Mounting the control panel on the support.**

- Connect the cable from the back of the control panel to the connection cable in the support.
- Place the control panel on the support so that the buttons on the back of the panel are seated in the recesses of the support.
- Then slide the control panel down to the stop.
- ✓ The control panel is now ready for use.

Connection cable

## ➤ Connecting peripheral devices

### Connecting cooling water


**Warning:**  
**Risk of explosion**

A closed valve in the recirculation system will build up high pressure in the glass coils and may result in glass bursting.

The cooling coils of the condenser may be exposed to a maximum cooling water pressure of 2 bar. To prevent a higher pressure from ever occurring, no valve or throttle must be installed in the cooling water return line.

**Assembly parts:**


- 2x tube nipple  
(included in scope of delivery)
- 1x tube, Ø inside 7-8 mm  
(not included in scope of delivery)

Coolant inlet and outlet are located on the back and bottom of the condenser.

- Remove the screw caps from the cooling water inlet and outlet.
- Push the screw caps over the tube nipples so that the black seal is located inside the cap.
- Push the respective inlet and outlet hoses for the cooling water through the screw cap over a tube nipple.
- Tighten the screw caps with the seal first on one of the respective openings of the condenser.
- ✓ Cooling water inlet and return line are connected.

Cooling water inlet and outlet

## Connecting vacuum

**Warning:**  
**Risk of implosion, risk of poisoning**

Glass flasks with even the finest damage may burst and, depending on the sample, the ambient air may be poisoned.

Do not increase the vacuum pressure unless necessary.

Never exceed the permissible pressure of your distillate.

To connect the vacuum, you need:


- a vacuum source
- a vacuum tube ( $\varnothing$  7-8 mm inside)

### Glassware sets G1, G3, G3 XL and G6

The opening for the vacuum connection is located on the back and top of the condenser.

- Unscrew a screw cap with the tube nipple.
- Push the vacuum tube onto the tube nipple.
- Tighten the screw cap with the tube nipple and the tube.
- Connect the other end of the vacuum tube with your vacuum source.

### Vacuum connection G1, G3, G3 XL and G6

Screw cap vacuum connection

### Glassware set G5

The opening for the vacuum connection is located on the back and top of the condenser.


### Vacuum connection G5

Opening vacuum connection

## Assembly

### Vacuum connecting scheme (demonstration)

## Optional accessories

### Safety hood


**Caution:**  
**Possible damage of the protective material**

At to high temperatures and/or unsuitable heating bath liquids the material becomes soft and may deform and possibly crack.


The safety hood is only suitable for heating bath temperatures < 130° C. (For suitable heating bath liquids please refer to Chapter "Start Up, Heating bath and flask lift" and chapter "Accessories and Spare parts, Accessories").

Before you can use the safety hood, you need to screw the support of the safety hood to the height stop.

- Unscrew the locking screw from the height stop band.
- Slide the safety hood off the support.
- Tighten the safety hood carrier on the height stop band using the two supplied Allen keys.
- ✓ The knurled screw on the safety hood carrier takes over the function of the locking screw for the height stop band.


- Loosen the clamping screw and slide the safety hood horizontally onto the carrier until the hood sits above the heating bath.
- Tighten the clamping screw.
- ✓ The safety hood is installed. The safety hood can be lifted open at the handle until it clicks into place.

### Safety shield

- Move the flask lift upward
- Hang the safety shield to the edge of the heating bath and position it so that the (potentially present) safety hood can be lifted open easily.
- Turn the eccentric screw 180° and use to secure the safety shield on the heating bath.


### Inlet valve

#### (Optional accessories) Mounting parts:

- 2x inlet tube
- 1x inlet valve

During a running process, the sample flow for the distillation is metered via the inlet valve.

- ✓ The glassware set is under vacuum.
- Turn the inlet valve anticlockwise.
- ✓ The in-feed is opened, and the distillate is drawn into the evaporating flask.
- ✓ The more you open the inlet valve, the faster your sample will be supplied.
- Turn the inlet valve clockwise.
- ✓ The in-feed is reduced, or the inflow of the sample is stopped.


#### Mounting the inlet valve with the inlet tubes

All glassware sets (with the exception of glassware set G1) have a screw opening on the lower left side of the condenser to accommodate the inlet tube with the inlet valve.

## Assembly


### Inlet tube 1

- Shorten the supplied PTFE tube so that the end of the tube ends in the evaporating flask.
  - Slide the drain disk onto the PTFE tube so that it is positioned before the vapor tube of the evaporator flask.
  - ✓ This prevents condensed fluid from flowing back into the evaporating flask via the inlet tube.
  - Unscrew the screw cap from the condenser.
  - Slide the seal ring and screw cap about 3 cm onto the PTFE tube.
- Insert the inlet valve with the inlet tube into the inlet opening of the condenser and tighten the screw cap.

### Inlet tube 2

- Slide the PTFE tube onto the inlet opening.
- ✓ The inlet hoses are installed.

### Inlet hoses Glassware set G1

The opening for the inlet tubes is located at the top of the condenser.

# ➤ Dismantling, transportation and storage

## Dismantling

---

### **Warning** **Risk of scalding!**

Hot glass surfaces and/or heating bath fluids may result in scalding. Allow all components to cool down before dismantling the device.

### **Risk of poisoning!**

Leakage of hazardous substances and solvents at leak points can cause hazardous substances to be released into the environment.

- Collect samples from the device in suitable containers.
- Before dismantling, ensure that no hazardous substances and solvents adhere to the device components.
- If necessary, run a distillation process with harmless and solvent-free samples before dismantling.

### **Risk of electric shock!**

If the device is still connected to the power supply, you may come in contact with live parts when you dismantle the device.

Before dismantling the device, switch it off and disconnect the mains plug from the socket. Do not reconnect the device to power under any circumstances.

---

- ✓ The device is switched off.
- ✓ The heating bath is emptied and cleaned.
- ➔ Remove all tube and cable connections to cooling liquid, vacuum and other devices.
- ➔ Remove all glassware and remove the heating bath.

## ➤ Scope of delivery

Component	Quantity	Product number
Hei-VAP Core: Hand lift*	1	571-01000-00
Hei-VAP Core: Motor lift*	1	572-01000-00
Glassware set	1	depending on variant
Vacuum seal PTFE/FKM	1	23-30-01-01-30
Vapor tube	1	514-00000-01
Clamping sleeve	1	23-30-01-05-31
Transportation lock (motor lift)	1	11-300-006-28
Coupling	1	23-09-03-01-03
Tension ring	1	22-03-02-01-05
Flask clamp Easy-Clip NS 29	1	23-30-01-05-29
Flask clamp Easy-Clip NS 24	1	23-30-01-05-57
Joint clamp	1	515-42000-00
Operating manual English / German	1	01-005-005-85
Warranty registration / Confirmation of condition	1	01-006-002-78
EU Declaration of Conformity	1	01-001-025-03
Power cord	1	depending on country

\* The product numbers are valid for 230 V devices in the European Community. Please contact Heidolph Instruments for further product numbers on country specific varieties.

## ➤ Accessories

Component	Quantity	Product number
Device fuses 230 V	2	14-002-028-22
Device fuses 115 V	2	14-002-028-29
Safety hood	1	569-00010-00
Safety shield	1	569-00020-00
5 L Heatingbath fluid Budde (20° C bis 170°C)	1	515-31000-00
5 L Heatingbath fluid Wacker AP 200 Siliconöl (150° C bis 210° C)	1	569-00600-00
Tube (suitable for vacuum and water)	1	591-35000-00
Vacuum seal PTFE/FKM	1	23-30-01-06-70
Extension cord for control panel 1.5 m	1	14-300-020-12
Glassware set	-	see homepage
Vacuum pumps and chillers (type and size vary)	-	see homepage

Please observe the following **permissible temperature range** when using **coated glassware sets**:

- G5 (dry Ice glassware set): -76° C bis 150° C
- any other glassware sets: -40° C to 80° C

Additional accessories are shown on our website at [www.heidolph-instruments.com](http://www.heidolph-instruments.com)

## Solvent data

The graph shows the relationship between the pressure and vapor temperature of a selection of solvents.

The temperature difference between the vapor temperature and the cooling medium should be at 20 K to result in sufficient condensation.

The temperature difference between the heating bath and vapor temperature should be at 20 K to reach a sufficient Distillation rate (dT).

i.e.: Set a vacuum for a vapor point at 40 °C, set the heating bath temperature at 60 °C.


Figure 11-1: Graph

Torr to mbar conversion: [mmHg]  $\approx$  3/4 [mbar]

## ➤ Technical data

### Hei-VAP Core

Standard supply voltage	230 V 50/60 Hz or 115 V 50/60 Hz	
Power consumption(W)	1400	
Device fuses	T8 H (230 V) or T15A H (115 V)	
Protection class (IEC 61140)	1 ⊕	
Protection class (IEC 60529)	<ul style="list-style-type: none"> <li>• Device: IP 20</li> <li>• Cable heating bath: IP 67</li> </ul>	
Noise pressure level (dB(A)) (based on IEC 61010)	< 85 dB(A)	
Lift distance (mm)	155	
Drive	EC motor	
Rotation speed (rpm)	10-280	
Speed control	electronic	
Heating bath	<ul style="list-style-type: none"> <li>• inner Ø: 253 mm, outer Ø: 291 mm</li> <li>• Material: V4A-steel (1.4404) stainless steel X2CrNiMo 17-12-2</li> <li>• Volume 4.5 L</li> </ul>	
Heating capacity (W)	1,300	
Temperature range (°C)	20 - 100 H <sub>2</sub> O / 20 - 210 oil	
Temperature control	electronic	
Temperature accuracy heating bath (°C)	±1	
Overtemperature protection heating bath	cut off at 5 °C difference to set temperature	
Overheat protection (°C)	250	
Condensing surface (cm <sup>2</sup> )	1,400	
Evaporation rates (L/h) ΔT* = 40°C (continuous operation)	<ul style="list-style-type: none"> <li>• Toluol 8.5</li> <li>• Acetone 5.8</li> <li>• Ethanol 3.5</li> <li>• Water 1.2</li> </ul>	
Display bath temperature (°C)	digital	
Display rotation speed	digital	
<b>Lift</b>	<b>Hand lift</b>	<b>Motor lift</b>
Weight (without glassware) (kg)	14,7	15
Lifting speed (mm/s)	-	25
Dimensions with glassware set G3 (w x d x h) (mm)	739 x 487 x 845	739 x 432 x 845

\* ΔT = Difference between heating bath temperature and boiling temperature

## ➤ Contact / Technical Service

### Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please get in touch with the following address.

For repairs please call Heidolph Instruments or your local authorized Heidolph Instruments distributor.

### **Warning!** **Danger of poisoning!**

Contaminated devices can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact substance
  - take adequate protective action towards our parts receiving and service personnel
  - mark the pack in accordance with Ordinance on Hazardous Substances
- 

A "Confirmation of Condition" form can be found at the end of this operating manual.

→ Prior to shipping the device for repair, complete a copy of this form and submit it in advance.

### Contact details:

---

#### **Heidolph Instruments Germany**

Heidolph Instruments GmbH & Co. KG  
Technical Service  
Walpersdorfer Str. 12  
D-91126 Schwabach / Germany

Tel.: +49 - 9122 - 9920-74

Fax: +49 - 9122 - 9920-84

E-Mail: [service@heidolph.de](mailto:service@heidolph.de)  
[www.heidolph.com](http://www.heidolph.com)

#### **Heidolph Instruments North America**

Phone: 1-866-650-9604  
E-mail: [service@heidolph.com](mailto:service@heidolph.com)  
[www.heidolphNA.com](http://www.heidolphNA.com)

#### **Heidolph Instruments United Kingdom**

Phone: 01799 - 5133-20  
E-mail: [service@radleys.co.uk](mailto:service@radleys.co.uk)  
[www.heidolph-instruments.co.uk](http://www.heidolph-instruments.co.uk)

#### **All other countries**

You will find contact details of your local Heidolph distributor at [www.heidolph.com](http://www.heidolph.com)

---

## ➤ Warranty

---

Heidolph Instruments provides a 3 year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online ([www.heidolph.com](http://www.heidolph.com)). The warranty period begins with the date of purchase. The serial number is also valid without registration.

This warranty covers material and manufacturing defects.

---

In the event of a material or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damage incurred as a result of improper handling or transport.

Warranty claim?

- Please notify Heidolph Instruments or your local Heidolph distributor if you wish to make a warranty claim.

## ➤ Certificate of decontamination

Enclose the Certificate of decontamination, duly completed, with your device return. Submissions without a Certificate of decontamination cannot be processed!

### Certificate of decontamination

#### IN CASE OF RETURNS

Please fill in the required fields.

**Note: The sender must package the goods properly and appropriately for transport.**

Heidolph Instruments GmbH & Co. KG  
 Walpersdorfer Straße 12  
 91126 Schwabach, Germany  
 Phone: +49 (0) 9122 9920-380  
**Fax: +49 (0) 9122 9920-19**  
 Email: service@heidolph.de

#### SENDER

Name _____	First name _____
Company/institution _____	Department _____
_____	Workgroup _____
Address _____	
ZC/City _____	
Country _____	Phone _____
Email _____	

#### DEVICE DETAILS

Article number _____	Serial no. _____
Ticket number _____	
Reason for sending in _____	
_____	

**Has the device been cleaned, decontaminated/disinfected?**      **Yes**      **No**      (Please mark as applicable)

If yes, which measures were carried out?

\_\_\_\_\_

**Yes**      **No**      (Please mark as applicable)

**Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous?**

If yes, with which substances did the device come into contact?

\_\_\_\_\_

#### LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

_____	_____	_____
Date	Signature	Company stamp

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hereby declare, that the product of the brand "Heidolph" (from serial number on 200067818) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Hei-VAP Core	517-61000-00, 571-00000-00, 572-00000-00 569-50000-00, 569-51000-00
Hei-VAP Expert	517-61000-00, 573-00000-00, 574-00000-00 569-60000-00, 569-61000-00
Hei-VAP Ultimate	517-61000-00, 573-00000-00, 574-00000-00 569-70000-00, 569-71000-00
Hei VAP Expert Control	517-61000-00, 573-00000-00, 574-00000-00 569-60000-00, 569-61000-00, 569-00200-00
Hei-VAP Ultimate Control	517-61000-00, 573-00000-00, 574-00000-00 569-70000-00, 569-71000-00, 569-00200-00

Maschinenrichtlinie / Machinery Directive 2006/42/EG  
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU  
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU  
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:  
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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## UKCA Declaration of Conformity

# UK CA

### Declaration of Conformity In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,  
Heidolph Instruments GmbH & Co. KG  
Walpersdorfer Straße 12  
91126 Schwabach / Germany

Product: Laboratory benchtop evaporator

Model:

Hei-VAP Core	517-61000-00, 571-00000-00, 572-00000-00, 569-50000-00, 569-51000-00
Hei-VAP Expert	517-61000-00, 573-00000-00, 574-00000-00, 569-60000-00, 569-61000-00
Hei-VAP Ultimate	517-61000-00, 573-00000-00, 574-00000-00, 569-70000-00, 569-71000-00
Hei VAP Expert Control	517-61000-00, 573-00000-00, 574-00000-00, 569-60000-00, 569-61000-00, 569-00200-00
Hei-VAP Ultimate Control	517-61000-00, 573-00000-00, 574-00000-00, 569-70000-00, 569-71000-00, 569-00200-00

Description:

Rotary Evaporators with heating bath and human machine interface

The object of the declaration described above is in conformity with the relevant UK Statutory Instruments (and their amendments):

2008 No. 1597

*The Supply of Machinery (Safety) Regulations*

2008 2016 No. 1091

*The Electromagnetic Compatibility Regulations 2016*

2012 No. 3032

*The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012*

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP ( UK ) Ltd.  
8. Northumberland Av.  
London WC2N 5BY

Signed for and on behalf of Heidolph Instruments GmbH & Co. KG  
Walpersdorfer Straße 12, 91126 Schwabach / Germany

Schwabach, 01.09.2021

Wolfgang Jaenicke  
Managing Director

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Instruments GmbH & Co. KG fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

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## China RoHS Declaration of Conformity

### China RoHS DECLARATION OF CONFORMITY

Heidolph Instruments GmbH & Co.KG has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

#### Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Instruments GmbH & Co.KG products is 25 years.

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This table is created according to SJ/T 11364-2014

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

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O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

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X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Instruments GmbH & Co.KG may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Instruments GmbH & Co.KG will not take responsibility for the EFUP of those products and appliances.

Schwabach, 06.08.2021

Wolfgang Jaenicke  
Chief Executive Officer CEO

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Instruments GmbH & Co. KG PRODUCTS
© Heidolph Instruments GmbH & Co. KG