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Datum \_\_\_\_\_

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Firmenstempel \_\_\_\_\_

Translation of the original instructions  
Page 58 – 105

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## Contents

### Introduction

About this document.....	58
Typographic conventions.....	58
Copyright protection.....	58
Safety instructions in the official languages of the European Union .....	58

### General notes

Directives applied, product certification.....	59
Copyright notice .....	59
California Residents .....	59
Residual risk.....	59
Intended use.....	59
Compliant use.....	60
Reasonably foreseeable misuse .....	60
Transportation.....	60
Storage.....	60
Acclimatization .....	60
Permissible ambient conditions .....	60
Energy saving potential and cost efficiency .....	60

### Safety

General safety instructions.....	61
Electrical safety.....	61
Data security .....	61
Operational safety .....	61
Overpressure operation .....	62
Environmental protection .....	62
Biohazard.....	62
Other regulations.....	62

### Assembly

Assembly .....	63
----------------	----

### Device description

Principle of operation .....	64
------------------------------	----

Mechanical design .....	65
Hei-VOLUME Distimatic Pro & Hei-VAP Industrial.....	66
Hei-VOLUME Distimatic Pro & Hei-VAP Expert/Ultimate .....	66
Heating tape (optional) .....	67
Signaling of the operating states (optional) .....	67
Structure of the user interface .....	68
Control and display elements .....	68
System time, timer, stopwatch.....	70
Color coding.....	70
Symbols.....	71
Input keys.....	71
User input .....	72

## Commissioning

Commissioning .....	72
Switching the device on/off .....	72
Calibrate the flask level sensor .....	73
Initial filling of the evaporation flask .....	74
Device configuration .....	75
Main Menu.....	75
[Applications] menu .....	76
[Settings] Menu .....	76
Screen lock.....	77
Display brightness .....	78
Calibrating the flask level sensor.....	78
Maximum pressure.....	78
Emergency stop sensor S3 sensitivity .....	78
Factory settings.....	78
Selecting the user language .....	78
Select formats .....	79
Automatic drainage of residue .....	79
Select evaporator type/glassware.....	80
Data recording .....	80
Connectivity.....	81
System check .....	82
Testing valves and signaling .....	82
Perform overpressure test.....	83
Testing the safety switch .....	84

## Operation

Select operating mode .....	85
Sensor mode .....	85
Parameter sensor mode .....	88
Hei-VAP parameters .....	90
Data logging .....	90
Recipes .....	91
Recipe list .....	91
Time mode .....	92
Time mode parameters .....	93
Rinse mode .....	96
Rinse mode parameters .....	97
Manual mode .....	98
Opening the error list .....	99
Resetting the device to factory settings .....	99
Opening system information .....	100
Updating the firmware and software .....	100

## Troubleshooting

Troubleshooting .....	101
-----------------------	-----

## Appendix

Technical specifications .....	102
Scope of delivery .....	102
Device service .....	103
General cleaning instructions .....	103
Repairs - Return of equipment .....	103
Maintenance .....	104
Disposal .....	104
Warranty statement .....	104
Contact details .....	104
Certificate of decontamination .....	105

### About this document

These operating instructions describe the functions and operation of the automatic Hei-VOLUME Distimatic Pro module. The operating instructions are an integral part of the described device.

### Typographic conventions

This document uses the following symbols, signal words, and highlights:

Description
Warning symbols in combination with a signal word indicate dangers: <b>DANGER</b>
 Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury. <b>WARNING</b>
Indicates a potential danger. Failure to respect the indications will result in serious injuries. <b>CAUTION</b>
Indicates a potential hazard. If not avoided, damage to property and minor to moderate injuries can occur.
Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.
[ ... ] Square brackets identify labels of control elements on the device as well as labels and entries in software masks and user interfaces.
→ The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.

### Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

### Safety instructions in the official languages of the European Union

A summary of all safety instructions in the official languages of the European Union can be found in our Safety Guide for the product group rotary evaporators (Ref. 01-005-006-53). This document is available for download on our homepage in the most up-to-date version.

## Directives applied, product certification

### CE marking

The device meets all requirements of the following directives:

- Low voltage Directive 2014/35/EU
- EMC Directive 2014/30/EU

### NRTL Certification

The device has been tested in accordance with the following standards:

- UL 61010-1:2012/R:2018-11  
CAN/CSA-C22.2 No. 61010-1:2012/A1:2018-11
- UL 61010-2-010:2019  
CSA-C22.2 No. 61010-2-010:2019
- Supplemented by UL 61010-2-081:2019
- Supplemented by CSA-C22.2 No. 61010-2-081:2019

## Copyright notice

The software implemented in this product is protected by copyright laws. The rights holder is Heidolph Scientific Products GmbH, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website [www.heidolph.com](http://www.heidolph.com).

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

## California Residents

Important information for California residents regarding Prop 65:  
Please visit [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov) for more information.

## Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

## Intended use

The automatic Hei-VOLUME Distimatic Pro module has been specially developed for the automatic filling and draining of Heidolph rotary evaporators. It is partly completed machinery as defined by the Machinery Directive 2006/42/EC.

The following processes can be controlled with the automatic module:

- Dosed addition of a starting material
- Drainage of distillate
- Drainage of residue (optional)

## Compliant use

The user is generally responsible for evaluating the conformity of his application and, if necessary, for taking additional measures.

## Reasonably foreseeable misuse

Additional measures may be necessary, and/or specific directives and safety regulations may have to be observed for use under conditions or for purposes deviating from the intended use (see, for example, section "Other regulations" on page 62). Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

## Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

## Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

## Acclimatization

After each transportation and after storage under critical climatic conditions (e.g., high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

## Permissible ambient conditions

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

## Energy saving potential and cost efficiency

Even in standby mode, small amounts of energy are used to supply individual device functions. However, especially over longer periods of nonuse, this resource consumption adds up and causes avoidable costs.

Therefore, always switch off the device when not in use, especially for long periods of time (devices with main switch: Main switch in position [O]; devices without main switch: Unplug) to effectively avoid unnecessary energy consumption and the resulting costs.

### General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

### Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device.

Prior to use, check that the device and the power supply cord are free of visible damage.

Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

### Data security

The user is responsible for ensuring data security when transferring data between the described device and other devices.

Only use secure networks for the data transfer and avoid use of critical infrastructure.

Only use high-quality shielded data cables for the data transfer.

For data transfer via a USB B connection, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

### Operational safety

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately. Shut down the device and disconnect it from the power supply if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all relevant general and safety instructions for the connected peripheral devices (observe the supplied documentation!).

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

## Overpressure operation

The distillate and residue are drained applying overpressure, i.e. the used distillate and residue containers are pressurized accordingly!

When using unsuitable distillate or residue containers, there is a risk that the containers will rupture or that liquid will escape uncontrolled and contaminate the environment.

There is also a high risk of injury, especially when processing toxic or corrosive substances!

Only use the distillate and residue containers recommended by Heidolph or others that are sufficiently pressure resistant and equipped with an adequately dimensioned pressure relief valve.

## Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

## Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

## Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operator is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unauthorized or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

# Assembly

### CAUTION

Due to improper system tubing and wiring as well as unauthorized changes to the system, there is a risk of direct and indirect damage to property!

For system tubing and wiring, follow all notes and instructions in the assembly instructions provided (part number 01-001-009-12-0).

Do not make any unauthorized changes to the device, system, tubing, nor the existing cable connections.

If necessary, contact our technical service (see "Contact details" on page 104).

If required, Heidolph Scientific Products GmbH offers a professional full assembly and commissioning service, including the device and all supplied accessories.

The majority of the sensors and valves of the automatic module are plugged into the coded connectors on the back of the Distimatic Pro (sensor cards).

Thanks to the coding of the connectors, all units are automatically detected and assigned in the device software!

All other units must be plugged into the correspondingly marked connection sockets. For a detailed description of all connections, please refer to the assembly instructions supplied. This document describes all steps for the correct connection of all tubes and cables between the Distimatic Pro, the components supplied (containers, sensors, valves, glassware) and the evaporation system in detail!

Faulty or disconnected units are detected and displayed when performing a system test (integrated software function).

The professional mounting and installation of the device including the supplied accessories is within the sole responsibility of the operator.

### Principle of operation

The automatic Hei-VOLUME Distimatic Pro module has been specially developed for use in combination with rotary evaporators. The automatic module enables the automatic filling and draining of the system in continuous operation.

The following processes are controlled with the automatic Hei-VOLUME Distimatic Pro module:

- Filling the evaporation flask
- Drainage of distillate
- Drainage of residue

The following operating modes are available to the user:

- Sensor mode: Standard mode for automatic filling and draining in continuous operation (see section "Sensor mode" on page 85).
- Time mode: The time mode is used specifically for processing media that form sticky residues or that tend to foam (see section "Time mode" on page 92).
- Rinse mode: The rinse mode is used to flush the supply lines and drains with a solvent for cleaning purposes (see section "Rinse mode" on page 96).
- Manual mode: In manual mode, the connected rotary evaporator can be operated for processing small quantities (see section "Manual mode" on page 98).

In automatic operation (sensor and time mode), the distillate and – depending on the variant (with or without automatic drainage of residue) – the residue is drained from the collector or the evaporation flask by applying overpressure. Distillate and residue can also be drained manually by the push of a button. The required overpressure is generated by a compressor integrated in the automatic module.

The evaporation flask is filled from the educt container using the vacuum in the rotary evaporator system. The required vacuum is generated by an external vacuum pump.

A pressure mode valve is used to switch between the compressor (overpressure operation) and the vacuum pump (underpressure operation).

### Mechanical design

The Distimatic Pro Base is housed in a sheet steel housing with protection class IP 42 with removable side covers and a removable control panel. It contains the main circuit board, all function modules and the compressor (without illustration) for generating the overpressure to drain the residue and the distillate.

#### Front panel controls

#### CAUTION

The control elements of the central display are touch-sensitive (touch panel) and can be operated by finger or with a suitable stylus.

Pointed or sharp objects can damage the surface of the display!

To operate the device functions, touch the display only with your fingertips or a special stylus with a soft-touch tip.

- 1 On/off switch
- 2 Removable control panel
- 3 Display
- 4 Push-and-rotary controls
- 5 Standby button

#### Connection sockets, rear side

The sensors and valves of the system are plugged into the coded connection sockets of the sensor cards. Thanks to the coding of the sockets, all units are automatically detected and assigned in the device software. Observe the marking of the existing connection sockets and all information on installing optional expansion cards in the supplied assembly documentation!

(Example illustration, the actual assignment depends on the scope of the order!)

- 1 Connection sockets for valves and sensors
- 2 Free slots for expansion cards
- 3 SUB-D socket recirculation chiller
- 4 Connection socket Hei-VAP
- 5 Connection external emergency off switch
- 6 Overpressure connection (HVI sensor)
- 7 Aspiration connection: filter or inert gas
- 8 Measuring amplifier terminal level sensor S4
- 9 IEC appliance inlet
- 10 Switchable built-in power socket

### Hei-VOLUME Distimatic Pro & Hei-VAP Industrial

Use of the automatic Hei-VOLUME Distimatic Pro module with a rotary evaporator Hei-VAP Industrial requires a special set of glassware consisting of a Woulff bottle [1], a collector vessel [2] and an educt sensor vessel [3].


The figure shows the maximum configuration with Woulff bottle, collector vessel and educt sensor vessel, including all valves and sensors, completely mounted on the housing of a rotary evaporator:

### Hei-VOLUME Distimatic Pro & Hei-VAP Expert/Ultimate

A Hei-VOLUME Distimatic Pro & Hei-VAP Expert/Ultimate system consists of pre-assembled assemblies (condenser, valve matrix).

The figures show the pre-assembled G9 XL condenser with level sensor S4 and emergency off sensor S3 and the valve matrix in standard configuration:

See also section "Assembly" on page 63.

### Heating tape (optional)

-if necessary, the automatic Hei-VOLUME Distimatic Pro module can be operated with a heating tape (not included in the standard scope of delivery, see section "Commissioning" on page 72, Heat control expansion card).

#### Principle of operation

The heating tape is wound around the tubing and preheats it during operation so that viscous media are passed reliably out of the system.

The heating tape is connected to the Distimatic Pro and is integrated in the system software via the corresponding Heat Control expansion card.

As soon as the heating tape has been detected correctly, a corresponding icon appears on the software interface. Touch the associated icon to switch on/off the heating tape:

The operating behavior is defined by the process parameters of the selected mode (automatic switching on, on-period/pause time in percent). Example: Period 100 %, defined heating tape on-period 30 %, pause time 70 %, for a 100-second cycle, equates to an on-period of 30 seconds per cycle.

### Signaling of the operating states (optional)


The different operating states of the automatic module can also be made visible over large distances by an optional signal tower and the associated signal light expansion card (signal tower accessories set, article number 591-66000-00). The installation of the additional card is described in detail in the enclosed sheet associated.

The following operating states are indicated:

flashing green	Process is finished
continuously lit green	Process is in progress
flashing yellow	Educt is empty or collection container is full
continuously lit yellow	Note the warning on the display
flashing red	Emergency stop has been triggered
continuously lit red	Device in fault condition, note the error message on the display

## Structure of the user interface

All device functions are controlled via the integrated control panel and the graphical user interface.

After switching on the device and initializing the device control, the Start window of the last active operating mode appears on the display (in the example: Sensor mode).

From the start window, you can access the main menu, select an operating mode directly, and select the time display.

The function keys in the footer can be used to drain the distillate, to drain the residue, and to start the process directly.

- 1** Main menu
- 2** Operating mode
- 3** System time/Timer
- 4** Parameter window with progress bar
- 5** Process Visualization
- 6** Heating tape indicator (only with heating tape connected)
- 7** Function buttons

## Control and display elements

Control element	Function
Function button [Main menu]	<p>Touch this function button to open the main menu with the following options:</p> <ul style="list-style-type: none"> <li>• [Applications] → Operating mode</li> <li>• [Settings] → System parameters</li> <li>• [Languages] → Menu language</li> <li>• [Formats] → System configuration and display formats</li> <li>• [Datalogging] → Data recording</li> <li>• [System check] → complete system test</li> <li>• [Error list] → Open error list</li> <li>• [Info] → Open system info</li> </ul>
Function button [Mode]	<p>Touch this function button to open the [Applications] menu directly. Operating modes:</p> <ul style="list-style-type: none"> <li>• [Sensor mode]</li> <li>• [Time mode]</li> <li>• [Rinse mode]</li> <li>• [Manual mode]</li> </ul>


Control element	Function
[System time/Timer] function button	<p>Touch this function button to toggle between the following functions:</p> <ul style="list-style-type: none"> <li>• System time display</li> <li>• Timer</li> <li>• Stopwatch</li> </ul> <p>See the following section "System time, timer, stopwatch" on page 70 for detailed information.</p>
[Drain distillate] button	<p>Touch this button to drain the distillate from the collector → only when no process is active!</p> <p>By touching this button, valves V2 and V4 are opened, see section "Principle of operation" on page 64.</p>
[Drain residue] button	<p>Touch this button to drain the residue from the evaporation flask → only when no process is active!</p> <p>By touching this button, valve V3 is opened, see section "Principle of operation" on page 64.</p>
[Start process] button	Touch this button to start the selected process.
Main display [Process Visualization]	The current process step is visualized in the main display.
Left push-and-rotary control	<p>In open lists, you can switch between entries by turning the left push-and-rotary control. A selected value (highlighted in orange) can be accepted by pressing the push-and-rotary control (similar to touching the entry on the display).</p> <p>The color of the LED ring indicates process status:</p> <ul style="list-style-type: none"> <li>• green → overpressure</li> <li>• blue → vacuum</li> <li>• white → ventilate &amp; atmospheric pressure</li> </ul>
Right push-and-rotary control	<p>Turn the right push-and-rotary control to switch between the display screens (in the case of multi-page displays, similar to the wiping movement on the display).</p> <p>When the on-screen keyboard is open, the display value can be increased or decreased by turning the right push-and-rotary control. The step width (1, 25 or 50) varies depending on the rotation speed of the push-and-rotary control.</p> <p>Color code LED ring:</p> <ul style="list-style-type: none"> <li>• LED ring lights up white when the on-screen keyboard is open</li> <li>• LED ring flashes red when the system is in a safe operating condition</li> </ul>
Standby button	<p>Press and hold both push-and-rotary controls simultaneously for approximately three seconds to lock the display.</p> <p>Press the standby button to put the device in standby mode. Press the button again to put the device into the previously active operating mode.</p> <p>Pressing the standby button will stop the current process!</p>

## System time, timer, stopwatch

Briefly touch the [System time/Timer] function button in the header of the start page to open the selection menu with the following options:

- System time display: current time (12 or 24-hour display, depending on the selection in the [Settings] menu)
- Stopwatch: In this mode, the actual time elapsed since the start of the process is displayed.
- Timer: In this mode, the time remaining until the scheduled end of the process is displayed.

This setting is applied in all operating modes and can be changed at any time!

### Set time

Switch to [Time] mode: Touch the [System time/Timer] function button for approx. two seconds to open the on-screen keyboard and correct the system time (date, time).

### Programming the Timer

Switch to [Timer] mode: Touch and hold the [System time/Timer] function button for approx. two seconds to open the on-screen keyboard and to program the timer.

### Color coding

In the process visualization area, all process lines and function elements are color-coded to illustrate the operating states:

Color	Meaning
Gray	[passive] or off state (e.g. non-active/closed line, non-switched valve, sensor [low])
Orange	[active] or on state (e.g. active/open line, switched valve, sensor [high])
White	[ventilated] condition, atmospheric pressure present
Light blue	[vacuum] condition, vacuum pump active, evacuation
Green	[overpressure] condition, compressor on, draining

### Symbols

In the process visualization, the following symbols are used to illustrate the process sequences and the integrated function elements:

Symbol	Meaning
	Compressor
	Vacuum pump
	Pressure mode valve
	2-way valve
	Sensor
	Heating tape

### Input keys

The following input keys are available in any menu:

Symbol	Meaning
	[delete] button. Press this button to correct a numeric or alphanumeric value entered.
	[ok/confirm] button. Press this button to confirm and accept a selection or input value.
	[discard] button. Press this button to discard a selection or input value and return to the previous view.

## User input

User input for configuration, parameterization and programming is carried out via the numerical and the alphanumeric on-screen keyboard. The on-screen keyboard appears automatically as soon as an input field is touched:

## Commissioning

### Switching the device on/off

To switch on and off, use the main switch on the front of the device. Observe the following sequence when turning the automatic module and attached peripherals on and off:

- Switching on the system: First, turn on all peripheral devices so that they are recognized when initializing the automatic module software.
- Turning off the system: First, turn off all peripherals to put the entire system in a safe operating condition.

#### **CAUTION**

If the Distimatic Pro is switched on and off incorrectly, there is a risk that peripheral devices will not be detected correctly or will switch in an uncontrolled manner! Observe the order given when switching the Distimatic Pro and the connected peripheral devices on and off.

In particular, before turning off the Distimatic Pro, make sure that all processes have been completed properly.

## Calibrate the flask level sensor

The flask level sensor (S4) must be calibrated during the initial commissioning and every time the solvent is changed.

### **WARNING**

In case of inadequate calibration or failure of a level sensor (S4), overflow/leaking of the flask contents is prevented by an emergency stop sensor (S3)!

As soon as this sensor responds in the cooler, the Distimatic Pro and the connected peripheral devices are switched to a safe operating state and no further medium is fed in.

In extreme cases, there is a danger of the emergency stop sensor (S3) not clearly detecting the overflow of very low polar media such as hexane or heptane. If you change to such a medium, contact our application specialists if necessary (see "Contact details" on page 104).

### **CAUTION**

The function of the level sensor in the evaporation flask (S4) may be affected by electromagnetic fields. This can lead to a loss of the operational quality of the entire system.

During calibration, commissioning and operation, ensure that there are no sources of interference in the immediate vicinity of the evaporation flask to ensure proper operation of the level sensor at all times.

- Make sure that the sensor is connected to the Distimatic Pro correctly.
- Switch on the Distimatic Pro and wait approximately 10 minutes for the electronics to warm up to operating temperature. In this way, measurement inaccuracies are avoided during calibration!

#### **Calibration to maximum sensitivity, low-polar solvents:**

- Remove the sensor from the evaporation flask.
- Hold the sensor in the open space.
- Make sure there are no objects or bodies in a hemispherical space one (1) meter in front of the sensor head. These would be detected as a signal and falsify the calibration result!
- Start the calibration process in the device software:
  - Open the main menu using the function button [Main menu], [1], in the header of the start window:
  - Open the [Settings] menu, select [Calibrate level sensor] and follow the instructions displayed.


#### **OR**

#### **Calibration on high-polar solvents and mixtures:**

- Make sure the sensor is in the flask.
- Fill the evaporation flask with at least 500 ml of the product.
- Slide the sensor into the stopped flask until the sensor head is just above the fluid level or – for example, in the case of mixtures with water – until the sensor head touches the fluid level slightly.
- Start the calibration process in the device software:
  - Open the main menu using the [Main menu] function button in the header of the start window.

- Open the [Settings] menu, select [Calibrate filling level] and follow the instructions displayed.
- Start the calibration process with [ok].
  - The calibration process can be canceled at any time with [discard].
  - The progress is visualized by the traffic light display. When the calibration is complete, the green dot on the right lights up.
- Confirm the completion of the calibration process with [ok].

### PRACTICAL TIP

When processing solvents with a tendency toward boiling retardation (e.g. ethanol), it is recommended to calibrate the sensor slightly less sensitively. This prevents the boiling retardation from being detected as a fluid and the fluid level in the rotary flask from decreasing excessively before adding new medium.

The sensor should only respond to liquids, never to boiling retardation or foam!

## Initial filling of the evaporation flask

Complete the following steps for filling the evaporation flask to ensure continuous supply of the educt during operation:

- Turn on the rotary evaporator and wait for the software to initialize (start screen appears on the control panel).
- Switch on the heating bath heater and the recirculation chiller.
- Turn on the Distimatic Pro and wait for the software to initialize (start screen appears on the control panel).
- Start any given process by touching the [Start process] button on the control panel of the Distimatic Pro.
- Open the inlet valve on the educt sensor vessel.
- Close the inlet valve when the displayed level is reached.

## Device configuration

The basic settings for operation are made in the main menu. After switching on the Distimatic Pro and initializing the software, the start page of the last active mode appears.

Open the main menu using the [Main menu] function button [1] in the header of the start window.

## Main Menu

### [Home] button

Open start window from any view

### [Applications]

Open [Applications] menu, select operating mode, see "[Applications] menu" on page 76.

### [Settings]


Open [Settings] menu: various basic settings can be adjusted here and the automatic module can be reset to the factory settings, see section "[Settings] Menu" on page 76.

### [Languages]

Open [Languages] menu and load user language, see "Selecting the user language" on page 78.

### [Formats]

Open [Formats] menu: here you can select the desired format for the display of numerical values, time, date and temperature. At the same time, the connected rotary evaporator / glassware is assigned and operation with/without drainage of residue is selected, see "Select formats" on page 79.

### [Data logging]

Start datalogging function, see "Data recording" on page 80.

### [Hei-CONTROL Pro]

Open [Network parameters] menu, see "Connectivity" on page 81.

### [System check]

Touch this button to perform a full system check, see "System check" on page 82.

### [Error list]

Touch this button to open the list of all error messages, see "Opening system information" on page 100.

### [Info]


Touch this button to open the information window of the system, including firmware and hardware data, see "Opening system information" on page 100.

## [Applications] menu

The desired operating mode can be selected in the [Applications] menu.

- Complete all ongoing processes before changing the operating mode.
- Touch the entry of the desired operating mode. The background color of the selected entry changes to orange, as shown.
- Press [ok] to confirm the selection and activate the selected mode.
  - The start window of the activated mode appears on the display.
- Press [discard] to discard the selection and to return to the previous view.

The [Applications] menu can also be opened directly by touching the [Mode] function button in the header of the start window.

## [Settings] Menu

Different device settings can be adjusted in the [Settings] menu.

These are described in detail in the following sections.

### Screen lock

The [Automatic screen lock] menu item can be used to activate the automatic screen lock for the display (active) or to deactivate it (inactive). When the screen lock is active, the display locks after a certain period of inactivity (no operator actions on the control panel).

The process progress window is displayed in the locked screen. As soon as the display is touched, the unlock slider appears.

The screen lock prevents accidental operation during operation, but does not provide protection against unauthorized access.

To protect the device from unauthorized access, the [PIN for Screen lock] option must be enabled!

#### Unlocking the screen

To unlock the screen, touch any point on the screen and swipe the slider from left to right:

#### Time until automatic screen lock is activated

Using the option [Look screen after ... min] you can specify after which period of time the automatic screen lock should be activated.

- Touch the function button to open the numeric on-screen keyboard and set the desired time interval in minutes.
- Confirm the input with [ok].
- Touch the [discard] button to discard the entry and return to the previous menu without making any changes.

From a continuous on-time of approx. 60 minutes, there is a risk of the display becoming burnt in and static display elements remaining visible as a visible pattern temporarily, even after switching off. It is therefore advisable to set a maximum time of 60 minutes before the screen lock is activated.

#### Screen lock PIN

The [PIN for Screen lock] option allows you to set a security code that must be entered to unlock the screen. This feature can be used to protect the device from unauthorized access.

- Touch the [PIN for Screen lock] function button to open the on-screen numeric keypad.
- Enter a four-digit PIN and press [ok] to confirm the entry.
- Re-enter the PIN to confirm (confirmation prompt!) and end with [ok] to activate the security code → [PIN for screen lock is active].
- Touch the [PIN for Screen lock] function button again to disable the security code → [PIN for screen lock] inactive.

The security code is automatically disabled when the screen lock is disabled.

### Display brightness

The brightness of the display can be adjusted via the menu item [Brightness].

- Touch the [Brightness] function button to open the on-screen numeric keyboard.
- Enter a percentage value between 10 and 100 (in increments of ten) and confirm the input with [ok].
- Touch the [discard] button to discard the selection and return to the [Settings] menu without making any changes.

### Calibrating the flask level sensor

The flask level sensor S4 can be calibrated via the menu item [Calibrate filling level], see section "Calibrate the flask level sensor" on page 73.

### Maximum pressure

The menu item [Max. pressure] can be used to set the maximum permissible absolute target pressure (corresponding to atmospheric pressure + compressor pressure) in a range from 1,100 – 1,400 mbar.

The higher the target pressure, the more efficiently the system is drained. An adjustment may also be necessary when greater distances or heights (e.g. tanks) need to be overcome.

Please note that exceeding a target pressure of 1,400 mbar during operation activates the evaporator's overpressure cut-out and the system will be ventilated.

### Emergency stop sensor S3 sensitivity

The menu item [Emergency stop sensor S3 sensitivity] can be used to adjust the sensitivity of the emergency stop sensor within a range of 0 – 100% in increments of [1%].

The parameter [Emergency stop sensor S3 sensitivity] may only be changed by authorized technical personnel during installation or after consultation with the application engineering department at Heidolph.

The default setting covers most applications; adjustment is only necessary with extremely low- or high-polar media.

### Factory settings

The device can be reset to the factory settings via the menu item [Factory settings], see section "Resetting the device to factory settings" on page 99

### Selecting the user language

The user language can be changed at any given time in the [Languages] menu. Available languages are English, German, Mandarin, French, Spanish, Portuguese, Italian, Russian and Polish.

- Touch the function button [Languages]: The list of all available menu languages opens.
- Touch the entry for the desired language: the background color of the selected entry changes to orange.
- Press [ok] to confirm the selection and to return to the start window of the active mode.

Touch the [discard] button to discard the entry and to return to the start window of the active mode without making any changes.

## Select formats

The following system settings can be made in the [Formats] menu:

- Set date format (US or EU)
- Set time format (12/24-hour)
- Set number separator (US or EU)
- Set temperature unit (degrees Celsius, Fahrenheit, Kelvin)

The active selection is highlighted in orange. Touch the desired function button to switch between options. Several changes can be made at the same time!

Confirm the new settings with [ok]. All changes are immediately accepted without a confirmation prompt!

To discard any changes and to return to the previous menu, touch the [discard] button.

X	Formats	✓
	DD.MM.YYYY	MM/DD/YYYY
	15:00	03:00 pm
	with autom. drainage of residue	without auto drainage of residue
	Hel-VAP Industrial	
,	,	
10,0	10.0	Hel-VAP Ultimate
°C	°F	K

## Automatic drainage of residue

In the [Formats] menu, you can define whether the residue is automatically drained from the evaporation flask (blowing out in overpressure operation), or whether the flask is drained by hand.

A system with automatic drainage of residue, in sensor or time mode can be operated automatically in continuous operation (condition: continuous supply with educt!).

### Operation with automatic drainage of residue

In operation with automatic drainage of residue, the residue is automatically drained under the following conditions:

- Stagnation reached in sensor mode.
- Distillate drain counter expired in sensor mode.
- Refill counter expired in sensor mode.
- Additional evaporation time for distillate drain or refill counter has expired in sensor mode.
- Number of defined steps (filling/evaporation) reached in time mode.

### Operation without automatic drainage of residue

In operation without automatic drainage of residue, the evaporation flask must be drained by hand.

### Drainage of residue selection

- The active selection is highlighted in orange when the [Formats] menu is opened. To change, touch the desired function button [with / without automatic drainage of residue].
- Confirm the new setting with [ok]. The change is accepted immediately without a confirmation prompt!
- Touch the [discard] button to discard the change.

## Select evaporator type/glassware

In the [Formats] menu, you can define the connected evaporator type including the glassware used. The selection stored in the system appears in the form of function buttons:

- The active selection (in the example: Hei-VAP Industrial with set of glassware R) is highlighted in orange. Touch the desired function button to switch between options.
- Confirm the new setting with [ok]. The change is accepted immediately without a confirmation prompt!
- Touch the [discard] button to discard the change.

The selection evaporator type/glassware is necessary for display purposes only.

## Data recording

Use the [Data logging] function to record the individual steps of a process. The function records each individual step with start and end time and stores this information in the Data logging memory. The function must be activated before the process starts:

- To activate data recording, switch from the start/overview page of the active mode to the parameters page.
- Touch the [Data logging] button in the footer.

When the data logging function is activated, the red recording dot appears in the function button and in the process progress display.

The recorded values can be read out, saved and deleted using the [Data logging] option.

- Open the [Main menu] and touch the [Data logging] function button.
- The [Data logging] list opens.
- The recorded values show when a process was started and stopped again.

Datalogging	
Start time: 06.02.2021 09:02	
End Time: 06.02.2021 09:02 2 Datapoints	

Datalogging	
Start time: 12.12.2020 14:36	
End Time: 12.12.2020 14:36 1 Datapoints	

In the event of a power interruption or when the device is switched off during the current data recording, the recording process is aborted and the corresponding entry marked with "(!)". The displayed data set contains the recorded process data that was recorded up to the power interruption.

In the Datalog export file, the corresponding file names are also given the addition "(!)" before the file extension: "(!).csv", or "(!).json".

## Deleting entries

To delete an entry from the list, touch the associated Delete button.

### Saving entries on a USB stick

Recorded entries can be saved on an USB-stick. To do this, connect a VFAT-formatted USB stick to the USB port of the display and open the data list as described.

As soon as a USB flash drive is connected, a [Save] button (disk icon) is assigned to each list entry. Touch the associated [Save] button to save an entry to the attached USB stick.

Large amounts of data are generated during the long-term recording of process data (several days). It is therefore recommended that you save the current data recording memory on a USB flash drive and remove the existing entries from the internal memory before starting the process.

## Connectivity

All necessary network parameters and basic settings for time synchronization can be defined in the [Hei-CONTROL Pro] menu:

- The connection status of the Hei-CONTROL Pro application is displayed in the [Status] field.
- All associated parameters can be defined manually in [static IP configuration] mode (checkbox deselected, figure on the left). Select the checkbox (figure on the right) to obtain the IP address of the device automatically.
- Select the [NTP server] checkbox to synchronize the system time of the device with a time server. This function requires a valid IP configuration!
- Place the cursor in the corresponding input field and use the on-screen keyboard to enter the required server address.

If the time is synchronized via a time server, the [System time/Timer] function button in the display header can only be used to select a time zone. If an ethernet connection exists, the system time is updated automatically via the defined time server. To define the system time manually, the [NTP server] option must be deselected!

The full functionality of the Hei-CONTROL Pro application is not available in the current device version.

## System check

All connected sensors and actuators or peripheral devices are displayed in the [System check] menu. Various tests can be carried out via the corresponding [System test] sub-menu to ensure the functional and operational safety of the system:

- Open the [Main menu] and touch the [System check] function button in the footer. The [System check] menu appears:
  - Error-free components are marked by a green checkmark
  - Faulty components are marked by a red X. Please note the corresponding error messages.
- Swipe up or down on the touch panel to scroll through the list.
- Touch the [TEST] function button in the header of the [System check] menu to open the [System test] sub-menu.

System check		✓
S4 flask level sensor (evaporation flask)	⚠	
S1 educt sensor (educt sensor vessel)	✓	
S2 collector sensor (collector vessel)	✓	
S3 Cut-out sensor	✓	
S8 overflow sensor (distillate container)	✓	
S9 overflow sensor (residue container)	✓	
BV1 venting valve	✓	
UV1 switchover vacuum-overpressure	✓	

### CAUTION

Faulty connections or signaling can lead to malfunctions or uncontrolled switching operations! Leaks reduce the performance of the entire system.  
Rectify displayed errors immediately.

## Testing valves and signaling

[All actuators] can be used to test the function of all connected valves and signal lights.

- Open the [System check] menu.

### Test all valves and signal lights

- In the footer of the [System check] overview, touch the [All actuators] function button to activate the test function (the button icon changes from gray to orange).
- Start the test with [Start test]:

System test		✓
BV1 venting valve		
UV1 switchover vacuum-overpressure		
Compressor		
V1 from educt container		
V2 to distillate container		
V3 to residue container		
V4 intermediate valve		
V5 to condenser		
V6 to collector		
Signal lights		
All actuators	TEST	Start test
Overpressure	Safety switch	

## Perform overpressure test

With this test function, the system can be tested for leaks.

### WARNING

Due to faulty components, faulty assembly of components or faulty connections, there is a risk that pressurized components will suddenly loosen or rupture during an overpressure test!


Before performing an overpressure test, ensure that all components are properly installed. In particular, make sure that all screw connections are properly tightened.

Before starting the overpressure test, close all safety covers, hoods and doors of the entire system.

Ensure that third parties are clear of the danger zone.

Use the required personal protective equipment (eye protection, hand protection, etc.).

- Open the [System check] menu.
- Touch the [Overpressure] function button in the footer to activate the test function (button changes to orange).
- Touch the [Start test] function button.
- Confirm the subsequent confirmation prompt with [ok] to carry out the overpressure test:
  - During the test, the system pressure is built up to the defined maximum pressure ([Settings] menu). This pressure level must be reached within a maximum of 30 seconds!
  - When the maximum pressure is reached, a measurement is made for a period of 120 seconds to determine whether the system pressure drops below a defined limit value.
  - After this time has elapsed, the system message [Overpressure test successfully completed] is shown on the display, provided that all required values have been achieved.

- The [Overpressure test failed] error message appears if the defined minimum pressure is not reached during the overpressure test or the system pressure falls below a defined limit value. In this case, follow the instructions on the display.
- The test can be stopped at any time with [Stop test].
- Close the menu with [ok].

### Testing the safety switch

The function of the integrated safety switch can be tested with the test function [Safety switch].

The safety switch is activated when a defined system overpressure is reached during compressor operation. When the safety switch is activated, the compressor is shut down so that the system pressure can be released.

- Open the [System check] menu.
- In the footer, touch the function button [Safety switch] to activate the test function (button changes to orange).
- Touch the [Start test] function button.
- Confirm the subsequent confirmation prompt with [ok] to perform the safety switch test:
  - During the test, the compressor is used to build up pressure in the system. After a certain time has elapsed, the pressure mode valve is automatically switched to vacuum mode by software and the compressor line is blocked abruptly.
- A corresponding error message is displayed if the safety switch does not trigger.
- If the safety switch triggers correctly, [The safety switch works properly] appears.

If the safety switch triggers correctly, the standard error message for switching off the system after the safety switch has been triggered also appears!

Follow the instructions on the display.

- The test can be stopped at any time with [Stop test].
- Close the menu with [ok].

## Select operating mode

After switching on the device and initializing the device control, the start window of the last active operating mode appears on the display. Touch the [Mode] function button to open the [Applications] menu and to activate the desired mode (see section “[Applications] menu” on page 76).

### Sensor mode

#### Activating sensor mode

- Open the menu [Applications] and touch the function button [Sensor mode]. The entry is highlighted in orange.
- Confirm the selection with [ok] to activate the mode and to open the sensor mode start window:

In sensor mode, the filling quantity in the evaporation flask is regulated by the flask level sensor S4. The following sensor statuses are possible:

- The filling process is started and runs as long as the sensor does not have contact with the medium → S4 in the status [[low]].
- The filling process is stopped when the sensor comes into contact with the medium → S4 in the status [[high]].

#### Preconditions



The following conditions must be met for operation in sensor mode:

- Heating bath heater of the rotary evaporator is switched on
- Recirculation chiller is switched on
- Vacuum pump is switched on
- Rotary evaporator is under filling vacuum
- Level sensor S4 calibrated
- Educt sensor vessel is filled

#### Filling the evaporation flask

Also observe information given in section “Initial filling of the evaporation flask” on page 74. After the initial filling of the evaporation flask, the filled educt is brought to process temperature.

- Touch the [[Start process]] button in the footer of the start window. The following message appears on the display:
- The filled educt is saved within a preset time (factory setting: 300 s) to the heating bath temperature.
- Observe the progress bar in the message box.
- Press [ok] to skip the heating process.

Note that skipping the heating process will automatically start the evaporation process. This activates the refill, the emergency stop and the educt sensor!

- After the pre-heating time has elapsed or after the step has been skipped, the filling pressure is automatically regulated to the set process pressure. The following message appears on the display:
- Start automatic mode with [ok] as soon as the evaporation process begins (distillate precipitates in the condenser).

### Automatic mode

The process time starts when automatic mode is started.

- As the evaporation starts, the filling level in the evaporation flask decreases:
  - The flask level sensor S4 loses contact with the medium and changes from [high] to [low]. With this status change, the filling time is started ( $\rightarrow$  [low] = start filling time).
- As soon as the defined setpoint for the filling time has been reached (filling time expired), the evaporation flask is automatically filled with new educt.
  - The stagnation time starts when the level sensor S4 changes from [low] to [high] during filling ( $\rightarrow$  [high] = start stagnation time).
- The filling level in the evaporation flask rises again. As soon as the level sensor S4 detects media contact, the status changes again from [low] to [high] ( $\rightarrow$  [high] = start stagnation time).
- The continuous evaporation process reduces the filling level in the evaporation flask again:
  - The level sensor S4 loses contact with the medium and changes from [high] to [low] ( $\rightarrow$  [low] = start filling time).
  - When the status changes again, the filling time is started and the stagnation time is set to zero.
  - As soon as the defined filling time has been reached again, the evaporation flask is filled with new educt.
- The filling level in the evaporation flask rises again and the described process steps are repeated until the end of the process.

The filling time defines the time until the next filling of the evaporation flask and thus the amount of educt that is evaporated until the next refilling.

The stagnation time runs as soon as the level sensor S4 changes to the [high] status after filling. If the stagnation time expires completely, the automatic residue emptying is initiated (depending on the set of glassware and the parameterization) or the flask must be emptied manually.

### End of process in sensor mode

In sensor mode, the process sequence is stopped under one of the following conditions:

- Educt container runs empty
  - After the filling time has elapsed, the evaporation flask is no longer filled → status change of the feed sensor S1 ([high] to [low]).
  - Unplanned end is initiated: Additional evaporation time is started (setting via process parameters, see section "Parameter sensor mode" on page 88).
  - After the additional evaporation time has elapsed, the system switches to a safe operating state: Peripheral shutdown (recirculation chiller and vacuum pump off, evaporator in standby).
- Maximum filling rate reached
  - Using the process parameter [Max. filling: additional evaporation time] a maximum filling rate can be set in order to achieve a defined concentration of the residue (setting via process parameters, see section "Parameter sensor mode" on page 88).
- Maximum number of distillate drainage steps reached
  - Using the process parameter [Max. [Drain distillate]: additional evaporation time], a maximum number of distillate drainages can be specified in order to achieve a defined concentration of the residue (setting via process parameters, see section "Parameter sensor mode" on page 88).
- Residue level in the evaporation flask reached, flask level sensor S4 remains permanently in the status [high] (medium contact), nominal stagnation time is reached.
  - In this case, the sensor control is stopped and the residue is drained automatically or manually. The draining process depends on the glassware (with/without drainage of residue) and the [Drain residue on/off] default in the parameters screen form (setting via process parameters, see section "Parameter sensor mode" on page 88).
- Maximum filling level in distillate/residue container → level sensor S8 (distillate container) or S9 (residue container) indicate by a status change from [low] to [high] that the maximum level has been reached (option).


Both of the level sensors S8 (distillate container) and S9 (residue container) are available as optional accessories and are not included in the standard scope of delivery!

The level sensors prevent overflow of the monitored vessel. When a level sensor is triggered, the system is switched to a safe operating state.

## Parameter sensor mode

### General parameters

- To adjust the general parameters, switch from the start page of sensor mode to the corresponding parameters screen form (tap the parameters overview or swipe from right to left).
- Use the two [Default] navigation keys to return to the start window or to continue to the [Recipes] view.

The following parameters can be adjusted:

Parameters	Function
[Cycles]	<p><b>Maximum number of evaporation cycles</b></p> <ul style="list-style-type: none"> <li>• Selection switch in position [infinite] → no limit.</li> <li>• Selection switch in position [n] → define a certain number of cycles in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> <p>Touch the [infinite] / [n] icon to switch between the two options.</p>
[Fill]	<p><b>Maximum filling rate</b></p> <ul style="list-style-type: none"> <li>• Selection switch in position [infinite] → no limit.</li> <li>• Selection switch in position [n] → define a maximum number of filling operations in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> <p>Touch the [infinite] / [n] icon to switch between the two options.</p>
[Drain residue]	<p><b>Residue draining time</b></p> <ul style="list-style-type: none"> <li>• Selection switch in position [off] → no defined draining time.</li> <li>• Selection switch in position [on] → define a residue draining time in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> <p>Touch the [on / off] icon to switch between the two options.</p>
[Drain distillate]	<p><b>Distillate draining time</b></p> <ul style="list-style-type: none"> <li>• In the selection field, define a distillate draining time (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul>

	<b>Maximum number of distillate drainage steps</b>
[Distillate]	<ul style="list-style-type: none"> <li>Selection switch in position [infinite] → no limit.</li> <li>Selection switch in position[n] → define a maximum number of distillate drainage steps in the selection window (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> <p>Touch the [infinite] / [n] icon to switch between the two options.</p>
[End of process]	Specify here whether the distillate and/or residue are drained at the end of the process. Touch the corresponding checkbox to activate (checkmark set) or deactivate the selection (x).

The [Fill] and [Distillate] parameters are mutually exclusive!  
For an evaporation cycle, you can specify either a maximum filling rate or a maximum number of distillate drainage steps.

- Confirm or discard all changes with [Accept] or [Cancel].

#### Process parameters

- Touch the [Process parameters] button in the footer to open the menu.
- Tap the button of the desired parameter in the menu and adjust the value using the on-screen keyboard, see also "User input" on page 72.

Parameters	Function
[Heating tape]	On-period of the heating tape in operation. Range: 10 – 100% in increments of ten. (see "Heating tape (optional)" on page 67)
[Heating tape automatic]	Automatically switch on heating tape yes/no. Range: on / off.
[Evacuation after draining the collector]	Evacuation period of the bypass line after draining the collector vessel. Range: 1 – 9999 s.
[Warm-up time, filling in sensor mode]	Warm-up time when filling the evaporation flask in sensor mode. Range: 1 – 9999 s.
[Educt empty: additional evaporation time]	When the educt container runs empty: Additional evaporation time for evaporation of the remaining educt in the evaporation flask. Range: 1 – 9999 s.


Maximum filling rate: Selectable when a specific number of filling operations is specified in the start window:

[Max. filling: additional evaporation time]

In the example: five fillings, then optionally an additional evaporation time runs. This allows a desired concentration to be achieved!

Selectable if a specific number of distillation drainage steps is specified in the start window:

[Max. drain distillate: additional evaporation time]

In the example: [Drain distillate] six times, then the optional additional evaporation time runs. This allows a desired concentration to be achieved!

The parameters [Max. filling: additional evaporation time] and [Max. [Drain distillate]: additional evaporation time] are mutually exclusive and cannot be activated at the same time!

Forced draining of the collector active: Collector will be drained after 60 – 9999 seconds, regardless of the filling level.

Range: 60 – 9999 s.

Activate/deactivate forced draining of the collector (determines previous parameter!).

Range: active / inactive

[Forced draining of the collector]

Time until the standby mode is activated after the educt container runs empty.

Range: 1 – 9999 s.

[Standby timer after educt container runs empty]

## Hei-VAP parameters

In the sub-menu [Hei-VAP parameters] various parameters of the connected Hei-VAP can be stored (vacuum, hysteresis, rotational speed, heating bath temperature, recirculation chiller temperature).

These values are for informational purposes only and may be included in a new recipe.

Touch the [Hei-VAP parameters] button in the footer of the Parameter screen form to open the submenu and to enter the required values for the connected rotary evaporator by hand.

## Data logging

See section "Data recording" on page 80.

## Recipes

By touching the [Save] button, the actual settings and values can be stored in the system as a recipe.

- Touch the [Save] button to open the on-screen keyboard (see also section "User input" on page 72):
- Use the on-screen keyboard to enter any recipe name.

It is not possible to assign recipe names twice!  
If a desired recipe name is already stored in the system, a corresponding warning message appears.

In this case, you can enter a different name or update the existing recipe with the new values.

- Confirm your entry with [ok] to store the new recipe or the new settings in the recipe list.

## Recipe list

To open the recipe list, switch from the start window of sensor mode to the corresponding [Recipes] view (two swipes from right to left, or use the navigation button from the parameters window).

Here, new recipes can be created ([New] button) and existing recipes can be loaded ([Load] button) or deleted ([Delete] button).

The entire recipe list can be saved on a USB stick. Recipe lists can also be loaded from a USB stick to the device memory.

To do this, connect a USB stick to the USB port on the control panel and follow the instructions on the display.

## Time mode

### Activating time mode


- Open the menu [Applications] and touch the function button [Time mode]. The entry is highlighted in orange.
- Confirm the selection with [ok] to activate the mode and open the time mode start window:

In the time mode, the filling of the evaporation flask, the distillation and the draining of the distillate and residue are carried out in a purely time-controlled manner.

- Within a cycle, a desired amount of concentrate is built up as a residue in the evaporation flask in the time mode.
- A cycle can consist of 1 – 9999 steps.
- One step consists of the filling and evaporation processes.
- In each step, the desired end concentration is reached.

### Preconditions

The following conditions must be met for operation in time mode:

- Heating bath heater of the rotary evaporator is switched on
- Recirculation chiller is switched on
- Vacuum pump is switched on
- Rotary evaporator is under filling vacuum
- Process parameters are set correctly on the rotary evaporator
- Level sensor S4 calibrated

Before the initial filling of the evaporation flask or before the first process step, it must be ensured in time mode that all process parameters (vacuum, temperature, speed, cooling) are set correctly on the rotary evaporator. Otherwise, there may be inaccuracies in the recipe processing!

### Procedure

- Initial filling in time mode: This step is only necessary to prevent drying or crystallization of the educt. Activate with the [Initial filling] parameter in the parameters screen form, see following section "Time mode parameters" on page 93. See also "Initial filling of the evaporation flask" on page 74.
- In time mode, a defined quantity is filled into the evaporation flask and immediately vaporized within a specified time:
  - The fill quantity is defined via [Filling time].
  - The [Evaporation time] is defined by the user.
- Due to the pre-evacuation of the evaporation flask, the evaporation process starts immediately while filling. After processing a step, the evaporation flask is filled again and the following step is processed:
  - The number of steps in a cycle is parameterized by the user. This value depends on the desired target quantity, which is to form in layers in the evaporation flask.

- After processing the specified number of steps, a cycle is run through and the residue is pressed out of the evaporation flask by means of overpressure.
  - The number of cycles is parameterized by the user.
- At the end of the last cycle, all connected devices are switched off (operating safe state, end of process).


In time mode, the system containers are monitored by sensors! The end process can be initiated and/or an ongoing process can be terminated and the evaporator can be put into a safe operating state when a corresponding signal is emitted by a sensor.

## Time mode parameters

### General parameters

- To adjust the general parameters, switch from the start page of time mode to the corresponding parameters screen form (tap the parameters overview or swipe from right to left).
- Use the two [Default] navigation keys to return to the start window or to continue to the [Recipes] view.

The following parameters can be adjusted:

Parameters	Function
[Pre-evacuation]	<p><b>Adjustable time for the pre-evacuation phase</b></p> <p>Before starting the process loops evaporation – filling – draining, the system is pre-evacuated for the time set here.</p>
[Initial filling]	<p><b>Adjustable time for the initial filling</b></p> <p>Before starting the process loops consisting of evaporation – filling – draining, the system can be filled initially.</p> <p>An initial fill ensures that the process is started with minimum quantity, thus avoiding possible recrystallization.</p> <ul style="list-style-type: none"> <li>• Selection switch in position [off] → no initial filling.</li> <li>• Selection switch in position [on] → define a time period for the initial filling in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> <p>Touch the [on/off] icon to switch between the two options.</p>

[Fill]	<b>Filling time setting</b> Adjust the filling time as required using the plus/minus buttons.
[Evaporate]	<b>Evaporation time setting</b> Adjust the evaporation time as required using the plus/minus buttons.
[End concentration]	<b>Additional evaporation step to reach a specific end concentration</b> Activate this option to pass an additional evaporation step after passing through the process loops to achieve a desired end concentration or to remove any residual moisture from the residue in solid-liquid mixtures. <ul style="list-style-type: none"> <li>Selection switch in position [off] → no additional evaporation step.</li> <li>Selection switch in position [on] → define a time period for one additional evaporation step (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> Touch the [on/off] icon to switch between the two options.
[Drain residue]	<b>Residue draining time</b> <ul style="list-style-type: none"> <li>Selection switch in position [off] → no defined draining time.</li> <li>Selection switch in position [on] → define a residue draining time in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> Touch the [on/off] icon to switch between the two options.
[Drain distillate]	<b>Distillate draining time</b> In the selection field, define a distillate draining time (adjustment via plus/minus buttons or tap/on-screen keyboard).
[End of process]	Specify here whether the distillate and/or residue are drained at the end of the process. Touch the corresponding checkbox to activate (checkmark set) or deactivate the selection (x).

- Confirm or discard all changes with [Accept] or [Cancel].

**Process parameters**

- Tap the [Process parameters] button in the footer.
  - The process parameters menu opens.
- Touch the function button of the desired parameter and adjust the value using the screen keyboard, see also "User input" on page 72.

Parameters	Function
[Heating tape]	Performance of the heating tape during operation. Range: 10 – 100% in increments of ten. (see "Heating tape (optional)" on page 67)
[Heating tape automatic]	Automatically switch on heating tape yes/no. Range: on / off.
[Evacuation after draining the collector]	Evacuation period of the bypass line after draining the collector vessel. Range: 1 – 9999 s.
[Educt empty: additional evaporation time]	When the educt container runs empty: Additional evaporation time for evaporation of the remaining educt in the evaporation flask. Range: 1 – 9999 s.
[Forced draining of the collector]	Forced draining of the collector active: Collector will be drained after 1 – 9999 seconds, regardless of the filling level. Range: 1 – 9999 s. Activate/deactivate forced draining of the collector (determines previous parameter!). Range: active / inactive
[Standby timer after educt container runs empty]	Time until the standby mode is activated after the educt container runs empty. Range: 1 – 9999 s.

- Hei-VAP parameters, section "Hei-VAP parameters" on page 90.
- Data logging, section "Data recording" on page 80.
- Save, section "Recipes" on page 91.

## Rinse mode


### Activating rinse mode

- Open the [Applications] menu and touch the [Rinse mode] function button. The entry is highlighted in orange.
- Confirm the selection with [ok] to activate the mode and to open the rinse mode start window.

The rinse mode allows the fill and concentrate line to be flushed with a suitable solvent. The flushing process removes residues from previous processes and impurities. When rotation is activated, the evaporation flask is also cleaned.

Multiple cycles should be carried out, depending on the degree of contamination!

To clean the complete system including the distillation area and the condenser, an evaporation process with a suitable solvent is necessary. Adjusting the parameters appropriately, the distillate produced this way will moisten all glass parts, thereby cleaning them.

### PRACTICAL TIP

- Select a vacuum value adequate for the solvent and set the hysteresis to 50 mbar.
- A large pressure range is passed through by the large hysteresis.
- All glass parts are first wetted, after which drops are formed, which then run off.
- As a result, the glass parts in the distillate area are cleaned within a short time.

### Preconditions

The following conditions must be met for rinsing the distillate tube:

- Heating bath heater of the rotary evaporator is switched on
- Recirculation chiller is switched on
- Vacuum pump is switched on

### Procedure

- Depending on the selection, the various lines of the system are flushed.
- To fill the evaporation flask, see also "Commissioning" on page 72.

## Rinse mode parameters

### General parameters

- To adjust the general parameters, switch from the start page of rinse mode to the corresponding parameters screen form (tap the parameters overview or swipe from right to left).
- Use the two [Default] navigation keys to return to the start window or to continue to the [Recipes] view.

The following parameters can be adjusted:

Parameters	Function
[Pre-evacuation]	<b>Adjustable time for the pre-evacuation phase</b> Before starting the process loops evaporation – filling – draining, the system is pre-evacuated for the time set here.
[Fill]	<b>Filling time setting</b> Adjust the filling time as required using the plus/minus buttons.
[Evaporate]	<b>Evaporation time setting</b> Adjust the evaporation time as required using the plus/minus buttons.
[Drain residue]	<b>Residue draining time</b> <ul style="list-style-type: none"> <li>• Selection switch in position [off] → no defined draining time.</li> <li>• Selection switch in position [on] → define a residue container draining time in the selection field (adjustment via plus/minus buttons or tap/on-screen keyboard).</li> </ul> Touch the [on/off] icon to switch between the two options
[Drain distillate]	<b>Distillate draining time</b> In the selection field, define a distillate draining time (adjustment via plus/minus buttons or tap/on-screen keyboard).
[End of process]	Specify here whether the distillate and/or residue are drained at the end of the process. Touch the corresponding checkbox to activate (checkmark set) or deactivate the selection (x).

## Process parameters

Tap the [Process parameters] button in the footer.

- The process parameters menu opens.
- Touch the function button of the desired parameter and adjust the value using the screen keyboard, see also "User input" on page 72.

Parameters	Function
[Heating tape]	Performance of the heating tape during operation. Range: 10 – 100% in increments of ten (see "Heating tape (optional)" on page 67)
[Heating tape automatic]	Automatically switch on heating tape yes/no. Range: on / off.
[Standby timer after educt container runs empty]	Time until the standby mode is activated after the educt container runs empty. Range: 1 – 9999 s.

- Hei-VAP parameters, section "Hei-VAP parameters" on page 90.
- Save, section "Recipes" on page 91.

## Manual mode

### Activating manual mode

- Open the [Applications] menu and touch the [Manual mode] function button. The entry is highlighted in orange.



Confirm the selection with [ok] to activate the mode and to open the manual mode start window.

In manual mode, single batches can be processed. For automatic filling of the evaporation flask and draining of the distillate, the [Drain distillate] and [Start process] buttons are available on the operating panel.

Operating parameters must be individually defined on the connected rotary evaporator.

### Procedure

- Touch the [[Start process]] button
  - The vacuum pump is switched on, the pressure in the evaporation flask is lowered to process pressure.
  - In the footer, the [Fill educt] and [Stop process] buttons are displayed

- Touch the [Fill educt] button to start the automatic filling of the evaporation flask.

By briefly tapping the [Fill educt] button, a small amount can be injected into the evaporation flask.

Touch the [Fill educt] button for approximately 2 seconds (observe button progress) to fill the evaporation flask continuously. The filling has to be stopped manually.

- Touch the [Fill educt] button again to stop filling.
- Touch the [Stop process] button to stop the process (end process).
- When the process is ended, the [Drain distillate] button appears on the control panel.
- Touch the [Drain distillate] button → distillate is drained (default: approx. 30 s).
- Touch the [Drain distillate] button again → drainage of the distillate is stopped.

The evaporation flask must be drained by hand.

Please observe the relevant information and instructions in the operating instructions for the rotary evaporator!

### Opening the error list

- To display the error list, open the Main menu and touch the [Error list] function button in the footer.
- To view an individual error message, tap the corresponding entry:
  - The error message appears in plain text on the display.
  - The entries cannot be deleted for service reasons!
- Close the error list with [ok]. The display returns to the Main menu.

### Resetting the device to factory settings

Use the following procedure to reset the device to the factory settings:

- From the Main menu, open the [Settings] menu.
- Touch the [Factory settings] function button.
- The adjacent prompt appears on the display.
- Confirm the confirmation prompt with [ok] to reset the device to the factory settings. The device must be restarted!

- Touch the [discard] button to discard the reset and return to the previous menu without making any changes.

### Opening system information

The [Information] view displays all version and license information for your automatic module.

- To view system information, open the Main menu and touch the [Info] function button in the footer.
- Use the [Versions] and [Licenses] buttons to access the desired information about the software and firmware for your product.
- Scan the stored QR code with a suitable device to access additional product information and relevant documents via our homepage.
- Close the view with [ok]. The display returns to the Main menu.

### Updating the firmware and software

The system is updated using an update file stored on a USB stick. The upgrade includes the system's software and firmware by default.

- Open the safety cover of the USB port on the control panel and insert the USB stick with the update file.
- Open the main menu.
- Tap the [System update] button at the bottom of the main menu (this button appears only if a USB stick with a valid update file is detected on the USB port!):  

The image shows a row of four dark grey rectangular buttons with white icons and text. From left to right: 1. 'System check' with a gear icon. 2. 'System update' with a circular arrow icon. 3. 'Error list' with a red 'X' icon. 4. 'Info' with an info circle icon.
- Confirm the confirmation prompt with [ok] to update the device software.
- Follow the instructions on the display
- Touch the [discard] button to discard the operation and to return to the previous menu without making any changes.
- Restart the device after the software update.
- After the restart, the firmware of the device can be updated. To do this, also confirm the corresponding confirmation prompt with [ok],
- After the update, unplug the USB stick and put the security cover back into the USB port.

## Troubleshooting

Error message	Possible remedy
Cable breakage or sensor defect	Check the cable of the displayed component.
Short circuit	Check the cable connections of the displayed component.
Values out of range!	Correct the values of the displayed parameter.
Compressor has been stopped	Excessive pressure from compressor, reduce pressure.
Pressure sensor failed!	Check pressure sensor for function, system check.
Error due to barometer!	Check the function of the barometer.
Temperature sensor failed!	Check temperature sensor for function.
Emergency stop has been triggered!	Eliminate the cause of the emergency stop from the external emergency stop switch.
Educt container S1 is empty!	Refill educt.
Overflow sensor S8 to the distillate container has triggered!	Drain the distillate container.
Overflow sensor S9 to the residue container has triggered!	Drain the residue container.
S3 cut out sensor has triggered!	Eliminate the cause of the emergency stop on the rotary evaporator.
Failure message from recirculation chiller	Observe the indicator on the recirculation chiller.
Invalid sensor plugged in!	Check the coding of the sensor, check that the sensor is in perfect condition.
Invalid valve plugged in!	Check the coding of the valve, check that the valve is in good condition.
Communication error!	Check the cable connection and connectors of the displayed components.
Operation mode cannot be started!	Open the error list, remove the indicated error cause(s).
... has been plugged in twice	Identical sensors or actuators plugged in, remove one unit.
Set switchpoint failed!	Switching error, contact technical service!
Calibration failed!	Repeat calibration procedure.
Update failed!	Check USB stick for correct connection, restart update.

All error messages and warning messages are shown in plain text on the display.

Follow the instructions on the display.

In case of recurring faults, please contact the responsible sales department or our technical service, see "Contact details" on page 104.

## Technical specifications

### General device data

Model	Automatic <b>Hei-VOLUME Distimatic Pro</b> module
Dimensions (W × H × D)	186 × 429 × 521 mm, including control panel
Weight	approx. 17 kg
Housing degree of protection	IP42
Control panel degree of protection	IP42
Control panel	detachable, 7 inch touchpanel, two control knobs

### Electrical data

Rated voltage	100 – 240 V AC; 50/60 Hz
Connection	L+N+PE
Protection class	I
Oversupply category	II
Degree of pollution	2
Power input	1,500 W max.
Fuse	15 A

### Permissible ambient conditions

Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl

## Scope of delivery

Item	Quantity	Product no.
Hei-VOLUME Distimatic Pro Basis	1	591-60001-00
Safety guide	1	01-005-006-53
Warranty registration	1	01-006-002-78

Further information, in particular about the available expansion modules and the available accessories, can be found on our website at [www.heidolph.com!](http://www.heidolph.com)

### Device service

#### DANGER

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

When cleaning, avoid the penetration of liquids.

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses. Further information on available accessories can be found on our website at [www.heidolph.com](http://www.heidolph.com).

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

### General cleaning instructions

#### CAUTION

Improper cleaning can damage the surfaces of the device. Penetrating liquid can damage the electronic components inside the device.

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

### Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
  - [service@heidolph.de](mailto:service@heidolph.de).
- In your message, please provide us with the following information in addition to an error description:
  - Item number
  - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 105.

### Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example) contact our technical service.

### Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

### Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG)

Heidolph Scientific Products GmbH, headquartered in 91126 Schwabach/Germany, address Walpersdorfer Str. 12, is registered under the number DE 50705753 in the electronic waste equipment register of the foundation ear ([www.ear-system.de](http://www.ear-system.de)).

### Warranty statement

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage [www.heidolph.com](http://www.heidolph.com).

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

### Contact details

#### Heidolph Scientific Products GmbH

Technical Service  
Walpersdorfer Str. 12  
D-91126 Schwabach/Germany  
Email: [service@heidolph.de](mailto:service@heidolph.de)

#### Representations

To find your local Heidolph distributor please visit [www.heidolph.com](http://www.heidolph.com)

### Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

### CERTIFICATE OF DECONTAMINATION IN CASE OF RETURNS

Please fill in the required fields.

**Note:** The sender must package the goods  
properly and appropriately for transport.

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12  
91126 Schwabach

Phone: +49 (0) 9122 9920-380  
**Fax:** +49 (0) 9122 9920-19  
E-Mail: service@heidolph.de

#### SENDER

Name \_\_\_\_\_

First name \_\_\_\_\_

Company/institution \_\_\_\_\_

Department \_\_\_\_\_

Address \_\_\_\_\_

Workgroup \_\_\_\_\_

ZC/City \_\_\_\_\_

Phone \_\_\_\_\_

Country \_\_\_\_\_

Email \_\_\_\_\_

#### DEVICE DETAILS

Article number \_\_\_\_\_

Serial no. \_\_\_\_\_

Ticket number \_\_\_\_\_

Reason for sending in \_\_\_\_\_

Has the device been cleaned, decontaminated/disinfected?

Yes

No

(Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due  
to the processing of substances that are hazardous to health,  
the environment and/or are biohazardous?

Yes

No

(Please mark as applicable)

If yes, with which substances did the device come into contact?

#### LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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Jörg Ziel  
Quality Manager

## UKCA Declaration of Conformity

In accordance with UK Government guidance

# UK DECLARATION OF INCORPORATION

This declaration is issued under the sole responsibility of the manufacturer,

**Heidolph Scientific Products GmbH**

Walpersdorfer Straße 12

91126 Schwabach / Germany

**Product:** Laboratory automatic device

*Model:*

Distimatic base (HVI/HVB)	591-60001-xx-x
Distimatic Pro Industrial	591-61xxx-xx-x
Distimatic Pro Benchtop	591-62xxx-xx-x
Distimatic Bench-top w/o drainage	591-51000-xx-x
Distimatic Industrial w/o drainage	591-52000-xx-x
Distimatic Bench-top w drainage	591-53000-xx-x
Distimatic Industrial w. drainage	591-54000-xx-x
Distimatic Industrial w. drainage	591-55000-xx-x

**Description:**

Device for filling and emptying of evaporation equipment.

The object of the declaration described above is classified as Partly Completed Machinery in accordance with UK Statutory Instrument 2008 No. 1597 Supply of Machinery (Safety) Regulations 2008.

The object of the declaration described above also is in conformity with the relevant UK Statutory Instruments (and their amendments):

- 2008/2016/EC The Electromagnetic Compatibility Regulations 2016  
2012/30/EU The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

This product is supplied as an optional accessory for incorporation to a Machine Tool.

The device has to be installed in accordance with the instructions contained in the Installation Guide. The partly completed machinery must not be put into service until it has been established that the machinery into which the partly completed machinery is to be incorporated complies with the provisions of a relevant product directive as a whole.

And complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP ( UK ) Ltd.

8. Northumberland Av. London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863.  
All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04.September.2024

Andreas Hahn  
Managing Director

Jörg Ziel  
Quality Manager

## China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in

electrical and electronic products will not leak or mutate under normal operating conditions.

During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

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O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

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X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

Schwabach, 04.September.2024

Andreas Hahn  
Managing Director

Jörg Ziel  
Quality Manager

© Heidolph Scientific Products GmbH

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Subject to change without notice. The latest version of this document can be downloaded at any time from [www.heidolph.com](http://www.heidolph.com).

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MATERIAL CONTENT DECLARATION FOR Heidolph Scientific Products GmbH PRODUCTS