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Operating instructions
Notice d'instructions

Hei-PLATE Mix 20 l
Hei-PLATE Mix'n'Heat Core

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Seite 4 – 21

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Datum _____

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Firmenstempel _____

Translation of the original instructions

Page 26 – 43

Traduction de la notice originale

Page 48 – 65

Zertifikate/Certifications

Contents

Introduction

About this document.....	26
Typographic conventions.....	26
Copyright protection.....	26
Safety instructions in the official languages of the European Union	26

General notes

Directives applied, product certification.....	27
Copyright notice.....	27
California Residents.....	27
Residual risk.....	27
Intended use.....	27
Compliant use.....	27
Reasonably foreseeable misuse.....	28
Transportation.....	28
Storage.....	28
Acclimatization	28
Permissible ambient conditions	28
Processing liquids at extreme cold temperatures	28
Energy saving potential and cost efficiency	29

Safety

General safety instructions.....	29
Electrical safety.....	29
Specific safety instructions for device interfaces	29
Operational safety	30
Occupational safety	30
Personal protective equipment (PPE).....	30
Environmental protection.....	30
Biohazard.....	31
Other regulations.....	31

Device description

Mechanical design.....	32
Hei-PLATE Mix 20 L.....	32
Hei-PLATE Mix'n'Heat Core	33

Commissioning

Positioning the device	34
Power supply	34
Switching the device on/off	34
Automatic restart	34
Activating the [Automatic restart] function	34
Deactivating the [Automatic restart] function	35

Operation

General notes	36
Setting the rotation speed.....	36
Heating mode (Hei-PLATE Mix'n'Heat Core only)	37
Set the heating temperature	37
Switching the heating function on/off	37
Operation with external temperature controller	37

Troubleshooting

Troubleshooting.....	38
----------------------	----

Appendix

Technical specifications	39
Scope of delivery.....	40
Accessories	40
Device service.....	40
General cleaning instructions.....	40
Repairs - Return of equipment.....	41
Maintenance.....	41
Disposal.....	41
Warranty statement.....	42
Contact details.....	42
Certificate of decontamination.....	43

About this document

These operating instructions describe all the functions and the operation of the magnetic stirrers, type Hei-PLATE Mix 20 I and Hei-PLATE Mix'n'Heat Core. The operating instructions are an integral part of the described device!

Typographic conventions

This document uses the following symbols, signal words, and highlights:

	Description
	Warning symbols in combination with a signal word indicate dangers: DANGER Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury. WARNING Indicates a potential danger. Failure to respect the indications will result in serious injuries. CAUTION Indicates a potential hazard. If not avoided, damage to property and minor to moderate injuries can occur.
	Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.
[...]	Square brackets identify labels of control elements on the device as well as labels and entries in software masks and user interfaces.
→	The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

Safety instructions in the official languages of the European Union

A summary of all safety instructions in this document, in all official languages of the European Union, can be found in our Safety Guide for the Hei-PLATE product group (Ref. 01-005-006-68). This document is available for download on our homepage in the most up-to-date version.

Directives applied, product certification

CE marking

The device meets all requirements of the following directives:

- Machinery Directive 2006/42/EC
- EMC Directive 2004/108/EC
- RoHS 2 directive 2011/65/EU, including amendment 2015/863

Copyright notice

The software implemented in this product is protected by copyright laws. The rights holder is Heidolph Scientific Products GmbH, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website www.heidolph.com.

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

California Residents

Important information for California residents regarding Prop 65:
Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The described devices have been developed and designed by the manufacturer for the following laboratory activities:

- Heating (only Hei-PLATE Mix'n'Heat Core)
- Mixing
- Stirring
- Titrating

Any other use of this device is not considered as intended!

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Compliant use

The compliance of each individual application must always be evaluated by the user. Any additional measures necessary to ensure compliance are always the responsibility of the user.

Reasonably foreseeable misuse

Additional measures may be necessary, and/or specific directives and safety regulations may have to be observed for use under conditions or for purposes deviating from the intended use (see, for example, section "Other regulations" on page 31). Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g., high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

Processing liquids at extreme cold temperatures

The device is suitable for the processing of liquids at extreme cold temperatures within the stated permissible ambient conditions and observing the regulations for proper use.

Always use a suitable insulating vessel (e.g., a Deware flask) to process liquids above a temperature of $< -15\text{ }^{\circ}\text{C}$ to avoid damages to the device due to condensation and frost.

Energy saving potential and cost efficiency

Even in standby mode, small amounts of energy are used to supply individual device functions. However, especially over longer periods of nonuse, this resource consumption adds up and causes avoidable costs.

Therefore, always switch off the device when not in use, especially for long periods of time (devices with main switch: Main switch in position [O]; devices without main switch: Unplug) to effectively avoid unnecessary energy consumption and the resulting costs.

General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device.

Prior to use, check that the device and the power supply cord are free of visible damage.

Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

Specific safety instructions for device interfaces

If the interface inlets and outlets are assigned too high voltages and are insufficiently insulated, metallic parts such as the housing can become live in the event of a fault.

- Safely isolate extra low voltage inputs and outputs via 25 V AC or 60 V DC according to EN 61140, or by double or reinforced insulation according to EN 60730-1 or DIN 60950-1.
- Use only shielded connection cables. Connect the shield to the connector housing. Liquid can penetrate open interfaces and come into contact with live parts. This causes a risk of malfunctions through to short-circuit.
- Make sure that unused interfaces are always covered (covers on the front and back of the device).

Operational safety

Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately.

Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Occupational safety

Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.

Persons with pacemakers or other susceptible implants must maintain a minimum distance of 10 cm from the device to prevent possible malfunctions.

Do not operate any other devices in the immediate vicinity of the device ...

- which can generate electromagnetic fields in the frequency range between 9×10^3 Hz to 3×10^{11} Hz,
- which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0.1 μm),
- which generate ultrasonic or ionizing waves.

Do not process any substances that could release energy in an uncontrolled manner (e.g. self-ignition).

Do not process substances in which the energy input through mixing poses a danger.

Wipe off any liquid that may have spilled on the device immediately.

Do not work with an open flame in the immediate vicinity of the device (risk of explosion).

Only use heat transfer agents and media that are explicitly approved and suitable for laboratory operation and for the respective working environment.

Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The required instruction of the personnel is solely within the operating company's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operating company is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Mechanical design

Hei-PLATE Mix 20 l

- 1 Stirring platform
- 2 On/off switch
- 3 Rotary control [rotation]

- 4 Appliance inlet

The appliance inlet at the rear of the device is fitted with a protective cover. No tools are needed to open the protective cover.

Hei-PLATE Mix'n'Heat Core

- 1 Heating plate
- 2 On/off switch
- 3 Rotary control [temperature]
- 4 Button + heating function on/off indicator
- 5 Rotary control [rotation]

- 6 Appliance inlet
- 7 Six-pole connection socket for external temperature controller EKT Hei-Con

The connections at the rear of the device are fitted with a protective cover. No tools are needed to open the protective cover.

Positioning the device

CAUTION

- To use the device, place it on a clean, stable, level, horizontal, and heat-resistant surface.
- Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.
- Before switching on, make sure that the device is sufficiently stable.
- Before switching on, always make sure that all setups are centered as far as possible and are fixed on the base unit with sufficient stability.
- Keep all support and contact surfaces clean and dry.
- During operation, make sure that the required minimum distance of 15 cm between the magnetic stirrer and other devices and structures is maintained.

Power supply

DANGER

Observe the notes and instructions given in section “Electrical safety” on page 29.

Switching the device on/off

To turn on the device, press the on/off button. When the device is switched on, the on/off button lights up green.

To turn off the device, press the on/off button again. When the device is turned off, the on/off button lights up white.

Automatic restart

If the [automatic restart] function is activated, the device switches on again automatically with the set parameters after a power interruption when the power is restored.

In the as-delivered condition, the [Automatic restart] function is disabled and has to be enabled manually.

Activating the [Automatic restart] function

Hei-PLATE Mix 20 I

- Make sure that the device is disconnected from the mains (mains plug removed!)
- Set the [Rotation] controller to 750 rpm.
- Press and hold the on/off button and connect the power supply cord of the appliance to a properly grounded mains socket.
- Release the on/off button as soon as it lights up green:
 - On/off button flashes green five times, function active, device ready for operation.

Hei-PLATE Mix'n'Heat Core

- Connect the power supply cord to a properly grounded power socket.
 - Device in standby mode, on/off button is lit white.
- Press the [heating on/off] button and keep it pressed.
- Briefly press the on/off button (on/off button lights up green) and hold down the [heating on/off] button until it flashes orange five times:
 - The LED indicator of the heating on/off] button goes out and does not light up until the heating function is switched on! Function active, device ready for operation.

Deactivating the [Automatic restart] function

Hei-PLATE Mix 20 l

- Make sure that the device is disconnected from the mains (mains plug removed!)
- Set the [Rotation] controller to 750 rpm.
- Press and hold the on/off button and connect the power supply cord of the appliance to a properly grounded mains socket.
- Release the on/off button as soon as it lights up green:
 - On/off button flashes green twice, function deactivated.

Hei-PLATE Mix'n'Heat Core

- Switch off all device functions (heating, rotation).
 - Device in standby mode, on/off button is lit white.
- Let the appliance cool down (heating plate temperature < 50 C° - observe the residual heat indicator)
- Press the [heating on/off] button and keep it pressed.
- Briefly press the on/off button and press the [heating on/off] button until it flashes orange twice:
 - The LED indicator of the heating on/off] button goes out and does not light up until the heating function is switched on!
 - On/off button lights up continuously, function deactivated, device ready for operation.

General notes

CAUTION

At high speeds in combination with a large load and/or a tall setup on the heating plate, there is a risk of the device swinging up and falling.

- Only increase the speed to the required level gradually under the named conditions, and pay attention to the stability of the setup.
- Lower the speed or reduce the total load on the heating plate if the device starts to move independently during operation.

If the distribution of the load on the heating plate is uneven, there is a risk of the device swinging up and falling.

- Always ensure that the vessels are evenly distributed on the heating plate.
- Pay particular attention to even distribution of the weight load when vessels of different size and/or differently filled vessels are placed on the heating plate at the same time!

High speeds carry the risk of improperly fastened vessels falling off the heating plate and contaminating the environment by leaking chemicals.

- Before switching on the device, make sure that all vessels are properly fixed on the heating plate.
- Clean contaminated surfaces immediately.

Always place individual vessels in the center of the heating plate.

Use suitable attachments for processing multiple samples/vessels simultaneously.

Do not switch on the appliance until all the vessels are securely positioned on the heating plate.

Always wear the specified personal protective equipment (PPE).

Setting the rotation speed

WARNING

When using open vessels, there is a risk that fluid will spray out.

- Whenever possible, use closed vessels for processing corrosive, toxic or biohazardous substances and seal them safely.
- Especially using open vessels, increase the rotation speed only gradually and observe the fluid movements.
- **Hei-PLATE Mix'n'Heat Core:** Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.

- Turn on the device as described in section "Switching the device on/off" on page 34
- Turn the [rotation] control clockwise or counter-clockwise to set the desired speed (setting range: 100 – 1,400 rpm)
- The rotation speed can be adjusted at any time.

Heating mode (Hei-PLATE Mix'n'Heat Core only)

Set the heating temperature

WARNING

The heating plate is designed for a maximum operating temperature of 300°C. Contact with surfaces above 50 °C can lead to serious injury.

- Avoid direct skin contact with the heating plate during operation and also observe the residual heat indicator after switching off!
- Do not place heat-sensitive objects on the heating plate.
- Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.

- Turn on the device as described in section "Switching the device on/off" on page 34
- Turn the [temperature] control clockwise or counter-clockwise to set the heating temperature (setting range: 20 – 300 °C).

Switching the heating function on/off

- Make sure the device is switched on and the heating temperature is set.
- Press the [heating on/off] button to activate the heating function.
 - When the heating function is activated, the LED ring of the [heating on/off] button and the corresponding [heating] indicator LED will light orange.
- The heating temperature can be adjusted at any time during operation by turning the [temperature] control.
- Press the [heating on/off] button again to switch off the heating function.

Operation with external temperature controller

WARNING

As long as the temperature controller is not immersed in the sample, the temperature of the ambient air is displayed. The heating plate can heat up to 300 °C unnoticed!

- Before switching on the device and for calibration, immerse the temperature controller in the sample!
- Do not place heat-sensitive objects on the heating plate.
- Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.

- Make sure the device is switched on and the heating temperature is preselected (see section "Set the heating temperature" on page 37).
- Immerse the temperature controller at least 20 mm deep in the sample to obtain a stable reading.
- When making measurements, follow the instructions in the EKT Hei-con manual!

For detailed instructions on setting up, connecting and operating with a temperature controller, refer to the operating instructions for the EKT Hei-Con (article number 01-005-004-61-61).

Troubleshooting

The following table includes possible failures and corresponding corrective measures:

Malfunction	Possible cause/ remedy
On/off button illumination remains off	→ Mains voltage not present: Check the power supply cable for damage, check the connection plug for correct seating, check the fuse of the house installation. → LED defective, contact technical service.
LED heating on/off does not light up	→ Heating plate does not heat: Device switched off or mains voltage not present: Switch on device, check the power supply cable for damage, check the connection plug for correct seating, check the fuse of the house installation. → LED defective, contact technical service.
No mixing function	→ No magnetic stirring bar in the vessel, insert stirring bars. → Rotation speed set to zero, increase rotation speed.
No heating function	→ Heating plate defective, contact technical service!
Deviation setpoint/ current value of sample	→ Setpoint too low, take into account heat losses. → Heating plate defective, contact technical service!

If a fault cannot be rectified with the described suggestions, please contact an authorized sales representative or our technical service (see section "Contact details" on page 42).

Technical specifications

General device data		
Hei-PLATE model	Mix 20 l	Mix'n'Heat Core
Dimensions (W × H × D)	168 × 101 × 299 mm	
Usable surface heating plate	Ø 145 mm	Ø 135 mm Ø 145 mm
Weight	approx. 3 kg	
Maximum permissible load	25 kg	
Drive	EC motor, left-turning	
Speed range	100 – 1,400 rpm	
Protection class (EN 60529)	IP42	
Acoustic pressure	< 50 db(A)	
Heating		
Heating power	n.a.	800 W at 230 V (EU) 600 W at 115 V (US)
Heating temperature range		20 – 300 °C
Indication/measurement range		-80 – 350°C (with temperature control unit EKT Hei-Con)
Heating control		PID
Electrical data		
Rated voltage	1/N/PE 230V AC 50/60 Hz or 1/N/PE 115V AC 50/60 Hz	
Protection class	I	
Overvoltage category	II	
Degree of pollution	2	
Power input	30 W	230 V: 830 W 115 V: 630 W
EMC class	B, Group 1	
Permissible ambient conditions		
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)	
Installation altitude	up to 2,000 m asl	

Scope of delivery

Item	Quantity	Product no.
Hei-PLATE Mix'n'Heat Core	1	506-11100-00
OR		
Hei-PLATE Mix 20 l	1	506-54100-00
Safety Guide	1	01-005-006-68
Warranty registration	1	01-006-002-78

Accessories

Further information on available accessories can be found on our website at www.heidolph.com.

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

DANGER

Electric shock

Live components are installed inside the device. When opening the device, there is a risk of touching live components.

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Penetrating liquid poses the risk of an electric shock.

When cleaning, avoid the penetration of liquids.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION

Property damages

Improper cleaning can damage the surfaces of the device. Penetrating liquid can damage the electronic components inside the device.

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
 - service@heidolph.de.
- In your message, please provide us with the following information in addition to an error description:
 - Item number
 - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 43.

Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example) contact our technical service.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG)

Heidolph Scientific Products GmbH, headquartered in 91126 Schwabach/Germany, address Walpersdorfer Str. 12, is registered under the number DE 50705753 in the electronic waste equipment register of the foundation ear (www.ear-system.de).

Warranty statement

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Contact details

Heidolph Scientific Products GmbH

Technical service

Walpersdorfer Str. 12

D-91126 Schwabach/Germany

Email: service@heidolph.de

Representations

To find your local Heidolph distributor please visit www.heidolph.com

Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

CERTIFICATE OF DECONTAMINATION IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12
91126 Schwabach

Phone: +49 (0) 9122 9920-380

Fax: +49 (0) 9122 9920-19

E-Mail: service@heidolph.de

SENDER

Name _____

First name _____

Company/institution _____

Department _____

Workgroup _____

Address _____

ZC/City _____

Country _____

Phone _____

Email _____

DEVICE DETAILS

Article number _____

Serial no. _____

Ticket number _____

Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected?

Yes

No

(Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous?

Yes

No

(Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200456156) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Hei-Mix S	503-02000-xx
Hei-PLATE Mix'n'Heat Core	506-1xxxx-xx
Hei-PLATE Mix'n'Heat Core+	506-2xxxx-xx
Hei-PLATE Mix'n'Heat Expert	506-3xxxx-xx
Hei-PLATE Mix'n'Heat Ultimate	506-4xxxx-xx
Hei-PLATE Mix 20 I	506-5xxxx-xx

Maschinenrichtlinie / Machinery Directive 2006/42/EG
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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Jörg Ziel
Quality Manager

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,
Heidolph Scientific Products GmbH
Walpersdorfer Straße 12
91126 Schwabach / Germany

Product: Laboratory magnetic stirrer

Model:

Hei-PLATE Mix'n'Heat Core	506-1xxxx-xx
Hei-PLATE Mix'n'Heat Core+	506-2xxxx-xx
Hei-PLATE Mix'n'Heat Expert	506-3xxxx-xx
Hei-PLATE Mix'n'Heat Ultimate	506-4xxxx-xx
Hei-PLATE Mix 20 I	506-5xxxx-xx
Radleys HP Core+	RR 91213-xx

Description:

Magnetid stirrer with or without heating plate and human machine interface

The object of the declaration described above is in conformity with the relevant
UK Statutory Instruments (and their amendments):

2008 No. 1597

2008 2016 No. 1091

2012 No. 3032

The Supply of Machinery (Safety) Regulations

The Electromagnetic Compatibility Regulations 2016

*The Restriction of the Use of Certain Hazardous Substances in
Electrical and Electronic Equipment Regulations 2012*

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP (UK) Ltd.
8. Northumberland Av.
London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH
Walpersdorfer Straße 12, 91126 Schwabach / Germany

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

© Heidolph Scientific Products GmbH

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Scientific Products GmbH PRODUCTS