**Safe**  
Heating and Mixing

# Betriebsanleitung Operating instructions Notice d'instructions

Hei-PLATE Mix'n'Heat Core<sup>+</sup>

 **heidolph**  
research made easy

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## Contents

### Introduction

About this document .....	28
Typographic conventions .....	28
Copyright protection .....	28
Safety instructions in the official languages of the European Union.....	28

### General notes

Directives applied, product certification .....	29
Copyright notice .....	29
California Residents .....	29
Residual risk .....	29
Intended use .....	29
Compliant use.....	29
Reasonably foreseeable misuse .....	30
Transportation.....	30
Storage.....	30
Acclimatization.....	30
Permissible ambient conditions.....	30
Processing liquids at extreme cold temperatures.....	30
Energy saving potential and cost efficiency .....	30

### Safety

General safety instructions.....	31
Electrical safety .....	31
Specific safety instructions for device interfaces.....	31
Operational safety .....	31
Occupational safety .....	32
Personal protective equipment (PPE) .....	32
Environmental protection.....	32
Biohazard .....	32
Other regulations .....	32

### Device description

Mechanical design .....	33
-------------------------	----

### Commissioning

Positioning the device .....	34
Power supply .....	34

Switching the device on/off.....	34
Settings .....	34
Open the settings menu.....	34
Automatic restart .....	35
Soft start.....	35
Auto standby function .....	36
Safety function.....	36
Sensor calibration.....	36

## Operation

General notes .....	37
Push-and-rotary control.....	37
Adjust values.....	37
Switch heating/mix function on/off.....	37
Display set values .....	37
Lock function.....	38
Set the heating temperature .....	38
Residual heat indicator .....	38
Setting the rotation speed .....	39
Operation with external temperature sensor.....	39
Troubleshooting.....	41
Error codes .....	41

## Appendix

Technical specifications.....	42
Scope of delivery .....	43
Accessories.....	43
Device service.....	43
General cleaning instructions .....	43
Repairs - Return of equipment .....	44
Maintenance.....	44
Disposal.....	44
Warranty statement .....	45
Contact details .....	45
Certificate of decontamination .....	46

### About this document

These operating instructions describe all the functions and the operation of the magnetic stirrers, type Hei-PLATE Mix'n'Heat Core. The operating instructions are an integral part of the described device!

### Typographic conventions

This document uses the following symbols, signal words, and highlights:

Description
Warning symbols in combination with a signal word indicate dangers: <b>DANGER</b>
 Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury.
<b>WARNING</b>
Indicates a potential danger. Failure to respect the indications will result in serious injuries.
<b>CAUTION</b>
Indicates a potential hazard. If not avoided, damage to property and minor to moderate injuries can occur.
Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.
[ ... ]
Square brackets identify labels of control elements on the device as well as labels and entries in software masks and user interfaces.
→ The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.

### Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

### Safety instructions in the official languages of the European Union

A summary of all safety instructions in this document, in all official languages of the European Union, can be found in our Safety Guide for the Hei-PLATE product group (Ref. 01-005-006-68). This document is available for download on our homepage in the most up-to-date version.

## Directives applied, product certification

### CE marking

The device meets all requirements of the following directives:

- Machinery Directive 2006/42/EC
- EMC Directive 2004/108/EC
- RoHS 2 directive 2011/65/EU, including amendment 2015/863

## Copyright notice

The software implemented in this product is protected by copyright laws. The rights holder is Heidolph Scientific Products GmbH, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website [www.heidolph.com](http://www.heidolph.com).

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

## California Residents

Important information for California residents regarding Prop 65:  
Please visit [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov) for more information.

## Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

## Intended use

The device has been developed and designed by the manufacturer for the following laboratory activities:

- Heating
- Mixing
- Stirring
- Titrating

Any other use of this device is not considered as intended!

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

## Compliant use

The compliance of each individual application must always be evaluated by the user. Any additional measures necessary to ensure compliance are always the responsibility of the user.

## **Reasonably foreseeable misuse**

Additional measures may be necessary, and/or specific directives and safety regulations may have to be observed for use under conditions or for purposes deviating from the intended use. Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

## **Transportation**

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

## **Storage**

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

## **Acclimatization**

After each transportation and after storage under critical climatic conditions (e.g., high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

## **Permissible ambient conditions**

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

## **Processing liquids at extreme cold temperatures**

The device is suitable for the processing of liquids at extreme cold temperatures within the stated permissible ambient conditions and observing the regulations for proper use.

Always use a suitable insulating vessel (e.g., a Deware flask) to process liquids above a temperature of < -15 °C to avoid damages to the device due to condensation and frost.

## **Energy saving potential and cost efficiency**

Even in standby mode, small amounts of energy are used to supply individual device functions. However, especially over longer periods of nonuse, this resource consumption adds up and causes avoidable costs.

Therefore, always switch off the device when not in use, especially for long periods of time (devices with main switch: Main switch in position [O]; devices without main switch: Unplug) to effectively avoid unnecessary energy consumption and the resulting costs.

### General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

### Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device.

Prior to use, check that the device and the power supply cord are free of visible damage.

Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

### Specific safety instructions for device interfaces

If the interface inlets and outlets are assigned too high voltages and are insufficiently insulated, metallic parts such as the housing can become live in the event of a fault.

- Safely isolate extra low voltage inputs and outputs via 25 V AC or 60 V DC according to EN 61140, or by double or reinforced insulation according to EN 60730-1 or DIN 60950-1.
  - Use only shielded connection cables. Connect the shield to the connector housing.
- Liquid can penetrate open interfaces and come into contact with live parts. This causes a risk of malfunctions through to short-circuit.
- Make sure that unused interfaces are always covers (covers on the front and back of the device).

### Operational safety

Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately.

Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

### Occupational safety

Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.

Persons with pacemakers or other susceptible implants must maintain a minimum distance of 10 cm from the device to prevent possible malfunctions.

Do not operate any other devices in the immediate vicinity of the device ...

- which can generate electromagnetic fields in the frequency range between  $9 \times 10^3$  Hz to  $3 \times 10^{11}$  Hz,
- which generate emission or radiation sources in the frequency range  $3 \times 10^{11}$  Hz to  $3 \times 10^{15}$  Hz (in the optical spectral range wavelengths from 1,000 µm to 0.1 µm),
- which generate ultrasonic or ionizing waves.

Do not process any substances that could release energy in an uncontrolled manner (e.g. self-ignition).

Do not process substances in which the energy input through mixing poses a danger.

Wipe off any liquid that may have spilled on the device immediately.

Do not work with an open flame in the immediate vicinity of the device (risk of explosion).

Only use heat transfer agents and media that are explicitly approved and suitable for laboratory operation and for the respective working environment.

### Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The required instruction of the personnel is solely within the operating company's responsibility.

### Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

### Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

### Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operating company is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

## Mechanical design

- 1 Heating plate
- 2 On/off switch
- 3 Rotary control [temperature]
- 4 Display
- 5 Rotary control [rotation]

- 6 Appliance inlet
- 7 Six-pole connection socket for external temperature controller EKT Hei-Con

The connections at the rear of the device are fitted with a protective cover.  
No tools are needed to open the protective cover.

### Positioning the device

#### CAUTION

- To use the device, place it on a clean, stable, level, horizontal, and heat-resistant surface.
- Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.
- Before switching on, make sure that the device is sufficiently stable.
- Before switching on, always make sure that all setups are centered as far as possible and are fixed on the base unit with sufficient stability.
- Keep all support and contact surfaces clean and dry.
- During operation, make sure that the required minimum distance of 15 cm between the magnetic stirrer and other devices and structures is maintained.

### Power supply

#### DANGER

Observe the notes and instructions given in section "Electrical safety" on page 31.

### Switching the device on/off

To turn on the device, press the on/off button. When the device is switched on, the on/off button lights up green.

To turn off the device, press the on/off button again. When the device is turned off, the on/off button lights up white.

### Settings

The following adjustments can be made in the device settings:

- Enable/disable auto standby
- Activating/deactivating the [automatic restart] function
- Activating/deactivating the [soft start] function
- Calibration of external temperature sensor (PT1000)
- Activating/deactivating the safety function

### Open the settings menu

- Turn off the device (on/off button lights up white).
- Press simultaneously and hold both the push-and-rotary controls [temperature] and [rotation].
- Press the on/off button until AR off (first setting option) appears on the display.
- Release the push-and-rotary controls and the on/off button.
- Use the [temperature] control to select one after another the following options:
  - [AR]: Activate or deactivate automatic restart
  - [SoS]: Activate or deactivate soft start
  - [ASd]: Activate or deactivate auto standby
  - [SAF]: Activate or deactivate the security function
  - [CAL]: Calibrate the connected temperature sensor

- Use the [rotation] control to select the desired mode:
  - [on] or [off] for automatic restart, soft start, auto standby and safety function
  - Temperature deviation in [°C] for sensor calibration

If no selection is made within about five seconds, the menu is automatically closed and the device is in the on state!

### Automatic restart

If the [automatic restart] function is activated, the device switches on again automatically with the set parameters after a power interruption when the power is restored.

By default, the [automatic restart] function is deactivated. The function can be activated or deactivated at any time via the device settings:

- Open device settings (see section "Open the settings menu" on page 34)
- Set [AR] to [on] to activate the function

#### Operation behavior when the mains is restored after a power interruption

- On/off button flashes white/green.
- Heating plate reheats to the set temperature value.
- Speed is ramped up to the set value.
- Once the set operating values have been reached, the on/off button will light up green again.

- Set [AR] to [off] to disable the function

### Soft start

When the [soft start] function is activated, the speed is gradually increased to the set value after the rotation is switched on.

By default, the [soft start] function is activated. The function can be activated or deactivated at any time via the device settings:

- Open device settings (see section "Open the settings menu" on page 34)
- Set [SoS] to [off] to disable the function
- Set [SoS] to [on] to enable the function

### Auto standby function

When the [auto standby] function is activated, the device is automatically switched off from idle status (no heating function, no agitation function) if no operation is performed for a period of three minutes.

The device must then be switched on again using the on/off button.

By default, the [auto standby] function is deactivated. The function can be activated or deactivated at any time via the device settings:

- Open device settings (see section "Open the settings menu" on page 34)
- Set [ASd] to [on] to enable the function
- Set [ASd] to [off] to disable the function

### Safety function

(Only relevant and available operating the device with an external temperature sensor, see also section "Operation with external temperature sensor" on page 39.)

If the safety function is activated, the heating function is automatically switched off when the sensor detects a sudden drop in temperature or when the measured temperature rises too slow.

By default, the safety function is activated. The function can be activated or deactivated via the device settings:

- Open device settings (see section "Open the settings menu" on page 34)
- Set [SAF] to [off] to disable the function
- Set [SAF] to [on] to activate the function

### Sensor calibration

(Only relevant and available operating the device with an external temperature sensor, see also section "Operation with external temperature sensor" on page 39.)

The calibration function can be used to set the connected temperature sensor according to the specific ambient conditions in your laboratory.

The sensor calibration function is only available when a compatible temperature sensor is connected. To calibrate the connected temperature sensor, follow these steps:

- Immerse the sensor in the reference sample
- Open device settings (see section "Open the settings menu" on page 34)
- Select parameter [CAL]
- Adjust the displayed temperature value as desired

The value is retained until a new calibration is performed.

## General notes

### CAUTION

At high speeds in combination with a large load and/or a tall setup on the heating plate, there is a risk of the device swinging up and falling.

- Only increase the speed to the required level gradually under the named conditions, and pay attention to the stability of the setup.
- Lower the speed or reduce the total load on the heating plate if the device starts to move independently during operation.

If the distribution of the load on the heating plate is uneven, there is a risk of the device swinging up and falling.

- Always ensure that the vessels are evenly distributed on the heating plate.
- Pay particular attention to even distribution of the weight load when vessels of different size and/or differently filled vessels are placed on the heating plate at the same time!

High speeds carry the risk of improperly fastened vessels falling off the heating plate and contaminating the environment by leaking chemicals.

- Before switching on the device, make sure that all vessels are properly fixed on the heating plate.
- Clean contaminated surfaces immediately.

Always place individual vessels in the center of the heating plate.

Use suitable attachments for processing multiple samples/vessels simultaneously.

Do not switch on the appliance until all the vessels are securely positioned on the heating plate.

Always wear the specified personal protective equipment (PPE).

## Push-and-rotary control

The device is equipped with two push-and-rotary controls for setting [temperature] and [rotation], for adjusting various device parameters, and for switching the heating and mix functions on and off (see section "Mechanical design" on page 33).

The operation of the push-and-rotary controls is described in detail in the following sections:

### Adjust values

Quickly turn the two knobs clockwise or counter-clockwise to increase or decrease the displayed value (temperature, speed, setting value).

### Switch heating/mix function on/off

Press the [temperature] control to turn the heating function on and off.

Press the [rotation] control to turn the mixing function on and off.

### Display set values

During normal operation, move the two controls one detent position clockwise or counter-clockwise to display the set point (temperature, speed).

### Lock function

Press and hold the [temperature] or the [rotation] push-and-rotary control for at least two seconds until the displayed value (temperature or speed) is highlighted with a white frame.

In this state, the push-and-rotary control is locked to prevent unintended operation.

Press and hold the [temperature] or the [rotation] push-and-rotary control again for at least two seconds until the white marker frame disappears. In this state, the push-and-rotary control is unlocked.

The two push-and-rotary controls [temperature] and [rotation] can be locked at the same time as described.

### Set the heating temperature

#### **WARNING**

The heating plate is designed for a maximum operating temperature of 300 °C. Contact with surfaces above 50 °C can lead to serious injury.

Always pay attention to the heating temperature indicator of the device.

Avoid direct skin contact with the heating plate during operation and also observe the residual heat indicator after switching off!

Do not place heat-sensitive objects on the heating plate.

Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.

- Turn on the device as described in section "Switching the device on/off" on page 34
  - The display shows the current temperature of the heating plate.
- Turn the [Temperature] push-and-rotary control clockwise or counter-clockwise to set the desired temperature (setting range: 20 – 300 °C).
  - The setpoint can be adjusted at any time during operation.
- Press the [Temperature] push-and-rotary control to activate the heating function.
  - When the heating function is activated, the LED ring of the [Temperature] control lights up orange.
- Press the [Temperature] push-and-rotary control again to turn off the heating function.

### Residual heat indicator

When the device is switched off, the display shows the current temperature of the heating plate. At the same time, the LED ring of the [Temperature] control flashes.

The optical indicators (display and LED ring) will turn off as soon as the surface temperature of the heating plate reaches a level of 50 °C with a falling tendency.

### **WARNING**

As long as the device is not switched off, the display shows the current temperature of the heating plate or the measured value from the external temperature sensor.

Note that the sensor reading may differ significantly from the current surface temperature of the heating plate!

Avoid direct skin contact with the heating plate!

Do not place heat-sensitive objects on the heating plate.

## Setting the rotation speed

### **WARNING**

When using open vessels, there is a risk that fluid will spray out.

Whenever possible, use closed vessels for processing corrosive, toxic or biohazardous substances and seal them safely.

Especially using open vessels, increase the rotation speed only gradually and observe the fluid movements.

Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.

- Turn on the device as described in section "Switching the device on/off" on page 34
  - The display shows 0 rpm.
- Turn the [Rotation] push-and-rotary control clockwise or counter-clockwise to set the desired nominal speed (setting range: 100 – 1,400 rpm).
- Press the [Rotation] control to activate the mixing function.
  - When the mixing function is activated, the LED ring of the [Rotation] control lights up white.
  - The displayed rotation speed value continuously updates until the setpoint is reached.
- The setpoint can be adjusted at any time:
  - Use the [Rotation] control to adjust the setpoint as required.
  - The rotation speed is adjusted without any delay (observe the display!).
- Press the [Rotation] control again to turn off the mixing function.

## Operation with external temperature sensor

Devices of the type Hei-PLATE Mix'n'Heat Core<sup>®</sup> can be operated with an external temperature controller (PT1000, optional accessory).

For detailed instructions on how to set up and connect a temperature controller, refer to the associated operating instructions.

When using an external temperature sensor, the current temperature is not measured on the heating plate, but on the connected sensor. The temperature indicator on the display is marked with the [ext] symbol.

Calibrate the external sensor as described in section "Sensor calibration" on page 36 to obtain correct readings and work results.

Two heating modes are available when using an external temperature sensor:

- FAST mode: In this mode, the heating plate is heated up quickly to the preset temperature (factory setting).
- PrECIS mode: In this mode, the heating plate is heated slowly and without overshooting to the preset temperature.

### **WARNING**

As long as the temperature sensor is not immersed in the sample, the temperature of the ambient air is displayed. The heating plate can therefore heat up to 300 °C unnoticed!

Always immerse the temperature sensor in the sample before switching on the device or before calibrating!

Do not place heat-sensitive objects on the heating plate.

Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.

### **Procedure**

- Ensure that the external temperature sensor is correctly connected.
- Switch on the device and set the heating temperature, see section "Set the heating temperature" on page 38.
  - The device operates in the last active heating mode (FAST or PrECIS).
- Press the [Temperature] control to change the heating mode if necessary (LED ring flashes rapidly).
  - Activate FAST mode: Press the control once and wait until the LED ring lights up continuously.
  - Activate PrECIS mode: Press the control twice and wait until the LED ring lights up continuously.
- Immerse the external temperature sensor at least 20 mm deep in the sample to obtain a stable reading.
  - The display shows the measured temperature of the sensor: Symbol [ext].
- Press the [Temperature] push-and-rotary control to activate the heating function.
  - When the heating function is activated, the LED ring of the [Temperature] control lights up orange.
- Press the [Temperature] push-and-rotary control again to turn off the heating function.

## Troubleshooting

Malfunction	Possible cause/ remedy
On/off button illumination remains off	<ul style="list-style-type: none"> <li>▪ Mains voltage not present: Check the power supply cable for damage, check the connection plug for correct seating, check the fuse of the house installation.</li> <li>▪ LED defective, contact technical service.</li> </ul>
No mixing function	<ul style="list-style-type: none"> <li>▪ No magnetic stirring bar in the vessel, insert stirring bars.</li> <li>▪ Rotation speed set to zero, increase rotation speed.</li> </ul>
No heating function	<ul style="list-style-type: none"> <li>▪ Heating temperature setpoint below current temperature, adjust setpoint.</li> <li>▪ Heating plate defective, contact technical service!</li> </ul>
Temperature display not plausible with temperature sensor connected	<ul style="list-style-type: none"> <li>▪ Check temperature sensor connection.</li> <li>▪ The immersion depth is too low, note the minimum immersion depth of 20 mm.</li> </ul>
Deviation setpoint/ current value of sample	<ul style="list-style-type: none"> <li>▪ Setpoint too low, take into account heat losses.</li> <li>▪ Heating plate defective, contact technical service!</li> </ul>

If a fault cannot be rectified with the described suggestions, please contact an authorized sales representative or our technical service (see section "Contact details" on page 45).

## Error codes

E11	Sensor break, internal sensor
E12	Sensor break or short circuit, internal hardware
E13	Internal sensor difference < 15 °C
E14	Internal safety temperature (T heating plate > Tset + 25 °C)
E21	External temperature sensor: media contact lost in the current process
E22	External temperature sensor without media contact when switching on
E23	External safety temperature exceeded
E36	IO expander: Loss of communication
E41	Drive failure
E51	External sensor break
E52	External sensor difference > 15 °C
E53	External temperature sensor disconnected or plugged in when the heating is switched on

## Technical specifications

### General device data

Model	<b>Hei-PLATE Mix'n'Heat Core+</b>
Dimensions (W × H × D)	168 × 101 × 299 mm
Usable surface heating plate	Ø 135 mm bzw. Ø 145 mm (depending on model, see ordering data)
Weight	approx. 3 kg
Maximum permissible load	25 kg
Drive	EC motor, left-turning
Speed range	100 – 1,400 rpm
Adjustment accuracy	5 rpm
Heating power	800 W at 230 V (EU) 600 W at 115 V (US)
Heating temperature range	20 – 300 °C
Indication/measurement range	15 – 325°C (with temperature sensor Pt1000)
Heating control	PID
Adjustment accuracy	1 K
Measurement accuracy (DIN IEC 751 Class a)	±0.2 K, plus tolerance PT1000
Measurement resolution	1 K
Display	LCD display
Protection class (EN 60529)	IP42
Acoustic pressure	< 50 db(A)
<b>Electrical data</b>	
Rated voltage	230 V, 50/60 Hz (EU) 115 V, 50/60 Hz (US)
Connection	L+N+PE
Protection class	I
Oversupply category	II
Degree of pollution	2
Power input	Normal operation 230 V: 825 W (EU) Normal operation 115 V: 625 W (US) Standby mode: 1.7 W
EMC class	B, Group 1
<b>Permissible ambient conditions</b>	
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 mamsl

### Scope of delivery

Item	Quantity	Product no.
Hei-PLATE Mix'n'Heat Core+	1	506-21100-00
Safety Guide		01-005-006-68
Warranty registration	1	01-006-002-78

### Accessories

Further information on available accessories can be found on our website at [www.heidolph.com](http://www.heidolph.com).

### Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

---

#### DANGER

##### Electric shock

Live components are installed inside the device. When opening the device, there is a risk of touching live components.

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Penetrating liquid poses the risk of an electric shock.

When cleaning, avoid the penetration of liquids.

---

### General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

---

#### CAUTION

##### Property damages

Improper cleaning can damage the surfaces of the device. Penetrating liquid can damage the electronic components inside the device.

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

---

### Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
  - service@heidolph.de.
- In your message, please provide us with the following information in addition to an error description:
  - Item number
  - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 46.

### Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example) contact our technical service.

### Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

### Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG)

Heidolph Scientific Products GmbH, headquartered in 91126 Schwabach/Germany, address Walpersdorfer Str. 12, is registered under the number DE 50705753 in the electronic waste equipment register of the foundation ear ([www.ear-system.de](http://www.ear-system.de)).

### **Warranty statement**

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage [www.heidolph.com](http://www.heidolph.com).

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

### **Contact details**

#### **Heidolph Scientific Products GmbH**

Technical service

Walpersdorfer Str. 12  
D-91126 Schwabach/Germany  
Email: [service@heidolph.de](mailto:service@heidolph.de)

#### **Representations**

To find your local Heidolph distributor please visit [www.heidolph.com](http://www.heidolph.com)

### Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

### CERTIFICATE OF DECONTAMINATION IN CASE OF RETURNS

Please fill in the required fields.

**Note: The sender must package the goods  
properly and appropriately for transport.**

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12  
91126 Schwabach

Phone: +49 (0) 9122 9920-380  
**Fax: +49 (0) 9122 9920-19**  
E-Mail: service@heidolph.de

#### SENDER

Name \_\_\_\_\_

First name \_\_\_\_\_

Company/institution \_\_\_\_\_

Department \_\_\_\_\_

Address \_\_\_\_\_

Workgroup \_\_\_\_\_

ZC/City \_\_\_\_\_

Phone \_\_\_\_\_

Country \_\_\_\_\_

Email \_\_\_\_\_

#### DEVICE DETAILS

Article number \_\_\_\_\_

Serial no. \_\_\_\_\_

Ticket number \_\_\_\_\_

Reason for sending in \_\_\_\_\_

**Has the device been cleaned, decontaminated/disinfected?**

Yes

No

(Please mark as applicable)

If yes, which measures were carried out?

**Does this device pose a risk to people and/or the environment due  
to the processing of substances that are hazardous to health,  
the environment and/or are biohazardous?**

Yes

No

(Please mark as applicable)

If yes, with which substances did the device come into contact?

#### LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200456156) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Hei-Mix S	503-02000-xx
Hei-PLATE Mix'n'Heat Core	506-1xxxx-xx
Hei-PLATE Mix'n'Heat Core+	506-2xxxx-xx
Hei-PLATE Mix'n'Heat Expert	506-3xxxx-xx
Hei-PLATE Mix'n'Heat Ultimate	506-4xxxx-xx
Hei-PLATE Mix 20 I	506-5xxxx-xx

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**UK  
CA**

## Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,  
Heidolph Scientific Products GmbH  
Walpersdorfer Straße 12  
91126 Schwabach / Germany

Product: Laboratory magnetic stirrer

*Model:*

Hei-PLATE Mix'n'Heat Core	506-1xxxx-xx
Hei-PLATE Mix'n'Heat Core+	506-2xxxx-xx
Hei-PLATE Mix'n'Heat Expert	506-3xxxx-xx
Hei-PLATE Mix'n'Heat Ultimate	506-4xxxx-xx
Hei-PLATE Mix 20 I	506-5xxxx-xx
Radleys HP Core+	RR 91213-xxx

**Description:**

Magnetid stirrer with or without heating plate and human machine interface

The object of the declaration described above is in conformity with the relevant  
UK Statutory Instruments (and their amendments):

2008 No. 1597	<i>The Supply of Machinery (Safety) Regulations</i>
2008 2016 No. 1091	<i>The Electromagnetic Compatibility Regulations 2016</i>
2012 No. 3032	<i>The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012</i>

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP ( UK ) Ltd.  
8. Northumberland Av.  
London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH  
Walpersdorfer Straße 12, 91126 Schwabach / Germany

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04.September.2024

Andreas Hahn  
Managing Director

Jörg Ziel  
Quality Manager

## China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

### Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions.

During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

MATERIAL CONTENT DECLARATION FOR Heidolph Scientific Products GmbH PRODUCTS							
部件名称 Part name	有毒有害物质或元素 Hazardous substances						
	铅 Pb	汞 Hg	铬 Cd	六价铬 Cr(VI)	多溴联苯 PBB	多溴二苯醚 PBDE	环保期限 标识 EFUP
包装 Packaging	O	O	O	O	O	O	
塑料外壳/组件 Plastic housing / parts	O	O	O	O	O	O	
电池 Battery	O	O	O	O	O	O	
玻璃 Glass	O	O	O	O	O	O	
电子电气组件 Electrical and electronic parts	X	X	X	O	O	O	
控制器/测量设备 Controller / measuring device	X	O	X	O	O	O	
金属外壳/组件 Metal housing /parts	X	O	O	O	O	O	
电机 Motor	X	O	O	O	O	O	
配件 Accessories	X	O	O	O	O	O	

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**Note:** Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

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O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.

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X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances .

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

Schwabach, 04.September.2024

Andreas Hahn  
Managing Director

Jörg Ziel  
Quality Manager

© Heidolph Scientific Products GmbH

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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