Betriebsanleitung
Operating instructions
Notice d'instructions

Hei-PLATE Core⁺ RS232

 **heidolph**
research made easy

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Seite 4 – 28

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Datum _____

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Firmenstempel _____

Translation of the original instructions
Page 34 – 58

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Contents

Introduction

About this document.....	34
Typographic conventions.....	34
Copyright protection.....	34
Safety instructions in the official languages of the European Union	34

General notes

Directives applied, product certification.....	35
Copyright notice	35
Residual risk	35
Intended use.....	35
Compliant use.....	35
Reasonably foreseeable misuse	36
Transportation.....	36
Storage.....	36
Acclimatization	36
Permissible ambient conditions	36
Processing liquids at extreme cold temperatures	36
Energy saving potential and cost efficiency	37

Safety

General safety instructions.....	38
Electrical safety.....	38
Interfaces	38
Data security	38
Operational safety	39
Occupational safety	39
Electrostatic charging of liquid or powdered media.....	39
Personal protective equipment (PPE).....	39
Environmental protection.....	40
Biohazard.....	40
Other regulations.....	40

Device description

Overview magnetic stirrer Hei-PLATE Core* RS232.....	41
--	----

Commissioning

Positioning the device	42
Power supply	42
Switching the device on/off	42
Settings	42
Open the settings menu	42
Automatic restart	43
Soft start	43
Auto standby function	44
Safety function.....	44
Sensor calibration	44

Operation

General notes	45
Push-and-rotary control	45
Adjust values	45
Switch heating/mix function on/off.....	45
Display set values.....	46
Lock function	46
Residual heat indicator.....	46
Operating modes.....	47
Continuous operation	47
Operation with external temperature sensor	48
Remote mode.....	49
Troubleshooting.....	50
Error codes	50
Interface commands	51

Appendix

Technical specifications	54
Scope of delivery	55
Accessories	55
Device service.....	55
General cleaning instructions.....	55
Repairs - Return of equipment.....	56
Maintenance.....	56
Disposal.....	56
Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG).....	56
Warranty statement.....	57
Contact details	57
Certificate of decontamination.....	58

About this document

These operating instructions describe all the functions and the operation of the magnetic stirrers, type Hei-PLATE Mix'n'Heat Core® RS232. The operating instructions are an integral part of the described device!

Typographic conventions

This document uses the following symbols, signal words, and highlights:

	Description
	Warning symbols in combination with a signal word indicate dangers: DANGER Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury. WARNING Indicates a potential danger. Failure to respect the indications will result in serious injuries. CAUTION Indicates a potential hazard. If not avoided, damage to property and minor to moderate injuries can occur.
	Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.
[...]	Square brackets identify labels of control elements on the device as well as labels and entries in software masks and user interfaces.
→	The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

Safety instructions in the official languages of the European Union

A summary of all safety instructions in this document, in all official languages of the European Union, can be found in our Safety Guide for the Hei-PLATE product group (Ref. 01-005-006-68). This document is available for download on our homepage in the most up-to-date version.

Directives applied, product certification

CE marking

The device meets all requirements of the following directives:

- Machinery Directive 2006/42/EC
- EMC Directive 2004/108/EC
- RoHS 2 directive 2011/65/EU, including amendment 2015/863

Copyright notice

The software implemented in this product is protected by copyright laws. The rights holder is Heidolph Scientific Products GmbH, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website www.heidolph.com.

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The device has been developed and designed by the manufacturer for the following laboratory activities:

- Heating
- Mixing
- Stirring
- Titrating

Any other use of this device is not considered as intended!

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Compliant use

The compliance of each individual application must always be evaluated by the user. Any additional measures necessary to ensure compliance are always the responsibility of the user.

Reasonably foreseeable misuse

Additional measures may be necessary, and/or specific directives and safety regulations may have to be observed for use under conditions or for purposes deviating from the intended use. Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

Transportation

During transportation, avoid severe shocks and mechanical stresses that can cause damage to the device.

Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g., high temperature difference between inside and outside), allow the device to acclimatize at room temperature for at least two hours to prevent possible damage due to condensation before commissioning it in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections (power supply, tubing) only after the device has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

Processing liquids at extreme cold temperatures

The device is suitable for the processing of liquids at extreme cold temperatures within the stated permissible ambient conditions and observing the regulations for proper use.

Always use a suitable insulating vessel (e.g., a Deware flask) to process liquids above a temperature of $< -15\text{ }^{\circ}\text{C}$ to avoid damages to the device due to condensation and frost.

Energy saving potential and cost efficiency

Even in standby mode, small amounts of energy are used to supply individual device functions. However, especially over longer periods of nonuse, this resource consumption adds up and causes avoidable costs.

Therefore, always switch off the device when not in use, especially for long periods of time (devices with main switch: Main switch in position [O]; devices without main switch: Unplug) to effectively avoid unnecessary energy consumption and the resulting costs.

General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device.

Prior to use, check that the device and the power supply cord are free of visible damage.

Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

Interfaces

Safely separate extra low voltage inputs and outputs rated 25 V AC or 60 V DC or above according to EN 61140, or use double or reinforced insulation according to EN 60730-1 or DIN 60950-1.

Use only shielded connection cables. Connect the shield to the connector housing. Always cover unused interfaces (front and back) to prevent liquids from entering the device.

Data security

The user is responsible for ensuring data security when transferring data between the described device and other devices.

Only use secure networks for the data transfer and avoid use of critical infrastructure.

Only use high-quality shielded data cables for the data transfer.

For data transfer via USB, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

Operational safety

Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately.

Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Occupational safety

Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.

Do not operate any other devices in the immediate vicinity of the device ...

- which can generate electromagnetic fields in the frequency range between 9×10^3 Hz to 3×10^{11} Hz,
- which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0.1 μm),
- which generate ultrasonic or ionizing waves.

Do not process any substances that could release energy in an uncontrolled manner (e.g. self-ignition).

Do not process substances in which the energy input through mixing poses a danger.

Wipe off any liquid that may have spilled on the device immediately.

Before each transfer or transport (e.g., removal from the fume hood, transfer between laboratory benches, temporary storage in the laboratory area, etc.), always remove all attachments from the shaking platform of the device.

Electrostatic charging of liquid or powdered media

Liquid or powdered media as well as the laboratory vessels used can become dangerously electrostatically charged due to the shaking movement or when emptying or cleaning the vessel. The intensity of the charge generated depends on the properties of the medium, the flow rate, the process parameters as well as the size, geometry and material of the laboratory vessels (source: TRGS 727, Section 4).

The evaluation of a potential explosion hazard due to electrostatic charging and the definition of suitable safety measures are always the responsibility of the operator and are always based on the substances to be processed and the individual application.

Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The required instruction of the personnel is solely within the operating company's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operating company is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Overview magnetic stirrer Hei-PLATE Core⁺ RS232

- 1 Heating plate
- 2 On/off switch
- 3 Display
- 4 Rotary control [temperature]
- 5 Rotary control [rotation]

- 6 Appliance inlet
- 7 RS232 interface
- 8 Six-pole connection socket for external temperature controller EKT Hei-Con

The connections at the rear of the device are fitted with a protective cover. No tools are needed to open the protective cover.

Positioning the device

CAUTION

- To use the device, place it on a clean, stable, level, horizontal, and heat-resistant surface.
- Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.
- Before switching on, make sure that the device is sufficiently stable.
- Before switching on, always make sure that all setups are centered as far as possible and are fixed on the base unit with sufficient stability.
- Keep all support and contact surfaces clean and dry.
- During operation, make sure that the required minimum distance of 15 cm between the magnetic stirrer and other devices and structures is maintained.

Power supply

DANGER

Observe the notes and instructions given in section “Electrical safety” on page 38.

Switching the device on/off

To turn on the device, press the on/off button. When the device is switched on, the on/off button lights up green.

To turn off the device, press the on/off button again. When the device is turned off, the on/off button lights up white.

Unplug the device to de-energize it!

Settings

The following adjustments can be made in the device settings:

- Enable/disable auto standby
- Activating/deactivating the [automatic restart] function
- Activating/deactivating the [soft start] function
- Calibration of external temperature sensor (PT1000)
- Activating/deactivating the safety function

Open the settings menu

- Turn off the device (on/off button lights up white).
- Press simultaneously and hold both the push-and-rotary controls [temperature] and [rotation].
- Press the on/off button until AR off (first setting option) appears on the display.
- Release the push-and-rotary controls and the on/off button.

- Use the [temperature] control to select one after another the following options:
 - [AR]: Activate or deactivate automatic restart
 - [SoS]: Activate or deactivate soft start
 - [ASd]: Activate or deactivate auto standby
 - [SAF]: Activate or deactivate the security function
 - [CAL]: Calibrate the connected temperature sensor
- Use the [rotation] control to select the desired mode:
 - [on] or [off] for automatic restart, soft start, auto standby and safety function
 - Temperature deviation in [°C] for sensor calibration

If no selection is made within about five seconds, the menu is automatically closed and the device is in the on state!

Automatic restart

If the [automatic restart] function is activated, the device switches on again automatically with the set parameters after a power interruption when the power is restored.

By default, the [automatic restart] function is deactivated. The function can be activated or deactivated at any time via the device settings:

- Open device settings
- Set [AR] to [on] to activate the function

Operation behavior when the mains is restored after a power interruption

- On/off button flashes white/green.
- Heating plate reheats to the set temperature value.
- Speed is ramped up to the set value.
- Once the set operating values have been reached, the on/off button will light up green again.

- Set [AR] to [off] to disable the function

Soft start

When the [soft start] function is activated, the speed is gradually increased to the set value after the rotation is switched on.

By default, the [soft start] function is activated. The function can be activated or deactivated at any time via the device settings:

- Open device settings
- Set [SoS] to [off] to disable the function
- Set [SoS] to [on] to enable the function

Auto standby function

When the [auto standby] function is activated, the device is automatically switched off from idle status (no heating function, no agitation function) if no operation is performed for a period of three minutes.

The device must then be switched on again using the on/off button.

By default, the [auto standby] function is deactivated. The function can be activated or deactivated at any time via the device settings:

- Open device settings
- Set [ASd] to [on] to enable the function
- Set [ASd] to [off] to disable the function

Safety function

(Only relevant and available when operating the device with an external temperature sensor, see section "Operation with external temperature sensor" on page 48.)

If the safety function is activated, the heating function is automatically switched off when the sensor detects a sudden drop in temperature or when the measured temperature rises too slow.

By default, the safety function is activated. The function can be activated or deactivated via the device settings:

- Open device settings
- Set [SAF] to [off] to disable the function
- Set [SAF] to [on] to activate the function

Sensor calibration

(Only relevant and available when operating the device with an external temperature sensor, see section "Operation with external temperature sensor" on page 48.)

The calibration function can be used to set the connected temperature sensor according to the specific ambient conditions in your laboratory.

The sensor calibration function is only available when a compatible temperature sensor is connected. To calibrate the connected temperature sensor, follow these steps:

- Immerse the sensor in the reference sample
- Open device settings
- Select parameter [CAL]
- Adjust the displayed temperature value as desired

The value is retained until a new calibration is performed.

General notes

CAUTION

At high speeds in combination with a large load and/or a tall setup on the heating plate, there is a risk of the device swinging up and falling.

- Only increase the speed to the required level gradually under the named conditions, and pay attention to the stability of the setup.
- Lower the speed or reduce the total load on the heating plate if the device starts to move independently during operation.

If the distribution of the load on the heating plate is uneven, there is a risk of the device swinging up and falling.

- Always ensure that the vessels are evenly distributed on the heating plate.
- Pay particular attention to even distribution of the weight load when vessels of different size and/or differently filled vessels are placed on the heating plate at the same time!

High speeds carry the risk of improperly fastened vessels falling off the heating plate and contaminating the environment by leaking chemicals.

- Before switching on the device, make sure that all vessels are properly fixed on the heating plate.
- Clean contaminated surfaces immediately.

Always place individual vessels in the center of the heating plate.

Use suitable attachments for processing multiple samples/vessels simultaneously.

Do not switch on the appliance until all the vessels are securely positioned on the heating plate.

Always wear the specified personal protective equipment (PPE).

Push-and-rotary control

The device is equipped with two push-and-rotary controls for setting [temperature] and [rotation], for adjusting various device parameters, and for switching the heating and mix functions on and off (see section "Overview magnetic stirrer Hei-PLATE Core+ RS232" on page 41).

The operation of the push-and-rotary controls is described in detail in the following sections:

Adjust values

Quickly turn the two knobs clockwise or counter-clockwise to increase or decrease the displayed value (temperature, speed, setting value).

Switch heating/mix function on/off

Press the [temperature] control to turn the heating function on and off.

Press the [rotation] control to turn the mixing function on and off.

Display set values

During normal operation, move the two controls one detent position clockwise or counter-clockwise to display the set point (temperature, speed).

Lock function

Press and hold the [temperature] or the [rotation] push-and-rotary control for at least two seconds until the displayed value (temperature or speed) is highlighted with a white frame.

In this state, the rotary function of the controller is disabled to prevent unintentional change of setting values. The press function to switch off the heating or rotation function of the appliance remains active!

Press and hold the [temperature] or the [rotation] push-and-rotary control again for at least two seconds until the white marker frame disappears. In this state, the push-and-rotary control is unlocked.

The two push-and-rotary controls [temperature] and [rotation] can be locked at the same time as described.

Residual heat indicator

When the device is switched off, the display shows the current temperature of the heating plate. At the same time, the LED ring of the [Temperature] control flashes.

The optical indicators (display and LED ring) will turn off as soon as the surface temperature of the heating plate reaches a level of 50 °C with a falling tendency.

WARNING

As long as the device is not switched off, the display shows the current temperature of the heating plate or the measured value from the external temperature sensor.

Note that the sensor reading may differ significantly from the current surface temperature of the heating plate!

Operating modes

The device can be used in continuous operation, with external temperature sensor or remote controlled. Regardless of the selected operating mode, observe the following safety instructions:

WARNING

The heating plate is designed for a maximum operating temperature of 300°C. Contact with surfaces above 50 °C can lead to serious injury.

Always pay attention to the heating temperature indicator of the device.

Whenever possible, use closed vessels for processing corrosive, toxic or biohazardous substances and seal them safely.

Avoid direct skin contact with the heating plate during operation and also observe the residual heat indicator after switching off!

Do not place heat-sensitive objects on the heating plate.

Always use the necessary personal protective equipment (heat-resistant gloves, eye protection, safety clothing) to process samples from a temperature of 50 °C.

Continuous operation

In continuous operation, the heating and rotation functions are set, started and stopped manually, see the following sections.

Set the heating temperature

- Switch on the device. The display shows the current temperature of the heating plate.
- Turn the [Temperature] push-and-rotary control clockwise or counter-clockwise to set the desired temperature (setting range: 20 – 300 °C). The setpoint can be adjusted at any time during operation.
- Press the [Temperature] push-and-rotary control to activate the heating function. When the heating function is activated, the LED ring of the [Temperature] control lights up orange.
- Press the [Temperature] push-and-rotary control again to turn off the heating function.

Setting the rotation speed

- Switch on the device. The display shows 0 rpm.
- Turn the [Rotation] push-and-rotary control clockwise or counter-clockwise to set the desired nominal speed (setting range: 100 – 1,400 rpm).
- Press the [Rotation] control to activate the mixing function.
 - When the mixing function is activated, the LED ring of the [Rotation] control lights up white.
 - The displayed rotation speed value continuously updates until the setpoint is reached.
- The setpoint can be adjusted at any time:
 - Use the [Rotation] control to adjust the setpoint as required.
 - The rotation speed is adjusted without any delay (observe the display!).
- Press the [Rotation] control again to turn off the mixing function.

Operation with external temperature sensor

Devices of the type Hei-PLATE Mix'n'Heat Core⁺ can be operated with an external temperature controller (PT1000, optional accessories, detailed instructions for setting up and connecting a temperature controller can be found in the associated operating instructions).

When using an external temperature sensor, the current temperature is not measured on the heating plate, but on the connected sensor. The temperature indicator on the display is marked with the [ext] symbol.

Calibrate the external sensor as described in section "Sensor calibration" on page 44 to obtain correct readings and work results.

Two heating modes are available when using an external temperature sensor:

- **FAST mode:** In this mode, the heating plate is heated up quickly to the preset temperature (factory setting).
- **PrECIS mode:** In this mode, the heating plate is heated slowly and without overshooting to the preset temperature.

WARNING

As long as the temperature sensor is not immersed in the sample, the temperature of the ambient air is displayed. The heating plate can therefore heat up to 300 °C unnoticed!

Always immerse the temperature sensor in the sample before switching on the device or before calibrating!

Procedure

- Ensure that the external temperature sensor is correctly connected.
- Switch on the device and set the heating temperature.
 - The device operates in the last active heating mode (FAST or PrECIS).
- Press the [Temperature] control to change the heating mode if necessary (LED ring flashes rapidly).
 - **Activate FAST mode:** Press the control once and wait until the LED ring lights up continuously.
 - **Activate PrECIS mode:** Press the control twice and wait until the LED ring lights up continuously.
- Immerse the external temperature sensor at least 20 mm deep in the sample to obtain a stable reading.
 - The display shows the measured temperature of the sensor: Symbol [ext].
- Press the [Temperature] push-and-rotary control to activate the heating function.
 - When the heating function is activated, the LED ring of the [Temperature] control lights up orange.
- Press the [Temperature] push-and-rotary control again to turn off the heating function.

Remote mode

WARNING

Secure the device in remote mode with a clearly visible warning sign and, if necessary, take further ambient-specific protective measures that protect against damage to property and injuries in the event of unexpected/unmonitored starting of the device.

To activate the remote access of the device, a message of type "SetTargetSpeed" or "SetMotorState" has to be sent. Afterwards, the remote symbol appears on the display and remote commands are executed normally. See section "Schnittstellenbefehle" on page 22 for detailed information.

Device activation via PC

The device can be addressed via the integrated RS232 interface via an external control. The data interface enables activation of the device and recording of the process data in remote mode. A prerequisite for this is a suitable terminal program or software solution such as the Hei-PROCESSING Solutions from Heidolph Scientific Products GmbH

Device activation via server-based software

The device can be addressed via the integrated RS232 interface via a server-based software. If required, contact an authorized dealer or our technical service, see section "Kontaktdaten" on page 28.

Operation behavior in remote operation

In remote mode, neither the rotation speed on the device can be adjusted nor the timer function can be activated. The [Timer] and [Fast Run] buttons and the push-and-rotary control are not working.

If necessary, the remote operation can be terminated by pressing the push-and-rotary control [Rotation] or the device can be switched off via the on/off switch.

Troubleshooting

Malfunction	Possible cause/ remedy
On/off button illumination remains off	- Mains voltage not present: Check the power supply cable for damage, check the connection plug for correct seating, check the fuse of the house installation. - LED defective, contact technical service.
No mixing function	- No magnetic stirring bar in the vessel, insert stirring bars. - Rotation speed set to zero, increase rotation speed.
No heating function	- Heating temperature setpoint below current temperature, adjust setpoint. - Heating plate defective, contact technical service!
Temperature display not plausible with temperature sensor connected	- Check temperature sensor connection. - The immersion depth is too low, note the minimum immersion depth of 20 mm.
Deviation setpoint/ current value of sample	- Setpoint too low, take into account heat losses. - Heating plate defective, contact technical service!

If a fault cannot be rectified with the described suggestions, please contact an authorized sales representative or our technical service (see section "Contact details" on page 57).

Error codes

E11	Sensor break, internal sensor
E12	Sensor break or short circuit, internal hardware
E13	Internal sensor difference < 15 °C
E14	Internal safety temperature (T heating plate > Tset + 25 °C)
E21	External temperature sensor: media contact lost in the current process
E22	External temperature sensor without media contact when switching on
E23	External safety temperature exceeded
E36	IO expander: Loss of communication
E41	Drive failure
E51	External sensor break
E52	External sensor difference > 15 °C
E53	External temperature sensor disconnected or plugged in when the heating is switched on

Interface commands

Query actual value

Command	Response	Description
IN_PV_1\X\n	IN_PV_1 X\X\n	X = actual value, sample temperature sensor (°C)
IN_PV_2\X\n	IN_PV_2 X\X\n	X = actual value, sample safety temperature (°C)
IN_PV_3\X\n	IN_PV_3 X\X\n	X = actual heating plate temperature value (°C)
IN_PV_4\X\n	IN_PV_4 X\X\n	X = actual heating plate safety temperature value (°C)
IN_PV_5\X\n	IN_PV_5 X\X\n	X = actual motor speed value (rpm)

Query setpoint value

Command	Response	Description
IN_SP_1\X\n	IN_SP_1 X\X\n	X = setpoint temperature, sample/heating plate (°C)
IN_SP_2\X\n	IN_SP_2 X\X\n	X = setpoint safety temperature delta (°C)
IN_SP_3\X\n	IN_SP_3 X\X\n	X = set motor speed value (rpm)

Query mode

Command	Response	Description
IN_MODE_1\Y\n	IN_MODE_1 Y\X\n	Query ext. Sensor Y = 0: No ext. sensor (heating plate temperature control) Y = 1: external medium temperature sensor (2 Pt1000)
IN_MODE_2\Y\n	IN_MODE_2 Y\X\n	Behavior after power failure Y = 0: All OFF after power restored, Y = 1: Heating/motor ON after power restored
IN_MODE_4\Y\n	IN_MODE_4 Y\X\n	Query heating control 0: Heating control in PRECISE mode 1: Heating control in FAST mode,
IN_MODE_6\Y\n	IN_MODE_6 Y\X\n	Query motor acceleration level 1: medium 2: fast
IN_MODE_7\Y\n	IN_MODE_7 Y\X\n	Query ext. Sensor monitoring 1: on 0: off

Query status, type, version

Command	Response	Description
STATUS\X\n	STATUS Y\X\n	Status query Y = 0: Manual operation at the device Y = 1: Remote operation START_1/START_2 Y = 2: Remote operation STOP_1/STOP_2 Y < 0: Error code Y = -1: Remote blocked (device stopped manually)
S_OP\X\n	S_OP Y\X\n	Query current control status Y = 0: All off; Y = 1: Heating on; Y = 2: Motor on; Y = 3: Heating and motor on
SW_VERS\X\n	SW_VERS V Y. Y. Y\X\n	Y.Y.Y = HMI version (e.g.: 1.2.345)
M\X\n	M\X\n	Identify device: Display of the tripped device begins to flash
T\X\n	T MR_CORE_PLUS\X\n	Query device type
QR_INFO\X\n	S#S#S\X\n	Barcode string, format: Material number: max. 18 digits alphanumeric. Serial number: max. 12 digits alphanumeric. Production order: max. 8 digits numerical Separator: #

Setting set values

Command	Response	Description
OUT_SP_1 Y\X\n	OUT_SP_1 X\X\n	Set setpoint temperature, sample/heating plate (°C)
OUT_SP_2 Y\X\n	OUT_SP_2 X\X\n	Set setpoint value, safety deltaT (°C)
OUT_SP_3 Y\X\n	OUT_SP_3 X\X\n	Set setpoint speed (rpm)

Setting modes

Command	Response	Description
OUT_MODE_2 Y\X\n	OUT_MODE_2 Y\X\n	Set behavior after power failure Y = 0: All OFF after power restored Y = 1: Heating/motor ON after power restored
OUT_MODE_4 Y\X\n	OUT_MODE_4 Y\X\n	Set heating control 0 = in Precise mode 1 = in Fast mode
OUT_MODE_6 Y\X\n	OUT_MODE_6 Y\X\n	Set motor acceleration level 1 = medium 2 = fast
OUT_MODE_7 Y\X\n	OUT_MODE_7 Y\X\n	Set ext. monitoring Sensor 1 = on 0 = off

Device control

Command	Response	Description
START_1\X\n	START_1\X\n	Start heating
START_2\X\n	START_2\X\n	Start rotation
STOP_1\X\n	STOP_1\X\n	Switch off heating
STOP_2\X\n	STOP_2\X\n	Switch off rotation
RESET\X\n	RESET\X\n	Reset all, heating and rotation off, deactivate remotely
PA_NEW\X\n	PA_NEW\X\n	Switch over to enhanced interface protocol
PA_OLD\X\n	PA_OLD\X\n	Switch over to old interface protocol
CC_ON\X\n	CC_ON\X\n	Connection test on
CC_OFF\X\n	CC_OFF\X\n	Connection test off

**RS-232 interface parameters**

- 9600 Baud
- Parity: even
- Data: 7 Bit
- Stop: 1 Bit

IN commands

- X = numerical value with one place after the decimal
- Y = Decimal number
- Z = multiple numerical values, separated by semicolon
- S = ASCII string
- If no external temperature sensor is connected, all relevant actual values are sent in **999.9** format.

Command termination

All ASCII strings must end with [...\X\n].

For further information on interface handling, please contact our technical service.

Technical specifications

Model	Hei-PLATE Mix'n'Heat Core+ RS232
Dimensions (W × H × D)	168 × 101 × 299 mm
Weight	approx. 3 kg
Usable surface heating plate	Ø 135 mm bzw. Ø 145 mm (depending on model, see ordering data)
Max. load	25 kg
Max. stirring volumen (H ₂ O)	20 L
Display	LCD display
Speed range	100 – 1,400 rpm
Adjustment accuracy	5 rpm
Drive	EC motor, left-turning
Operating mode	Continuous operation, remote operation
Heating power	800 W at 230 V. 600 W at 115 V.
Heating control	PID
Heating temperature range	20 – 300 °C
Adjustment accuracy	1 K
Measurement accuracy (DIN IEC 751 Class a)	±0.2 K
Measurement resolution	1 K
Protection class	IP42 (acc. to IEC 60529)
Protection class	1 (according to IEC 61140)
Overvoltage category	2
Degree of pollution	2
Acoustic pressure	< 50 (dB(A)) (in accordance with IEC 61010)
Electrical connection	1/N/PE 230 V AC 50/60 Hz 1/N/PE 115 V AC 50/60 Hz
Power input	630 W at 115 V 830 W at 230 V 1.7 W in standby mode
EMC class	B, Group 1
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 mamsl
Interfaces	1 × RS232 1 × PT1000

Scope of delivery

Item	Quantity	Product no.
Hei-PLATE Mix'n'Heat Core+ RS232	1	506-21218-00-0
Safety Guide		01-005-006-68
Warranty registration	1	01-006-002-78

Accessories

Further information on available accessories can be found on our website at www.heidolph.com.

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

DANGER

Electric shock

Live components are installed inside the device. When opening the device, there is a risk of touching live components.

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

Penetrating liquid poses the risk of an electric shock.

When cleaning, avoid the penetration of liquids.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION

Property damages

Improper cleaning can damage the surfaces of the device. Penetrating liquid can damage the electronic components inside the device.

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
 - service@heidolph.de
- In your message, please provide us with the following information in addition to an error description:
 - Item number
 - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 58.

Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example) contact our technical service.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG)

Heidolph Scientific Products GmbH, headquartered in 91126 Schwabach/Germany, address Walpersdorfer Str. 12, is registered under the number DE 50705753 in the electronic waste equipment register of the foundation ear (www.ear-system.de).

Warranty statement

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Contact details

Heidolph Scientific Products GmbH

Technical service

Walpersdorfer Str. 12

D-91126 Schwabach/Germany

Email: service@heidolph.de

Representations

To find your local Heidolph distributor please visit www.heidolph.com

Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

CERTIFICATE OF DECONTAMINATION IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12
91126 Schwabach

Phone: +49 (0) 9122 9920-380

Fax: +49 (0) 9122 9920-19

E-Mail: service@heidolph.de

SENDER

Name _____

First name _____

Company/institution _____

Department _____

Workgroup _____

Address _____

ZC/City _____

Country _____ Phone _____

Email _____

DEVICE DETAILS

Article number _____

Serial no. _____

Ticket number _____

Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected?

Yes

No

(Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous?

Yes

No

(Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200456156) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Hei-Mix S	503-02000-xx
Hei-PLATE Mix'n'Heat Core	506-1xxxx-xx
Hei-PLATE Mix'n'Heat Core+ /RS232	506-2xxxx-xx
Hei-PLATE Mix'n'Heat Expert	506-3xxxx-xx
Hei-PLATE Mix'n'Heat Ultimate	506-4xxxx-xx
Hei-PLATE Mix 20 I	506-5xxxx-xx

Maschinenrichtlinie / Machinery Directive 2006/42/EG
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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Jürgen Pirthauer
Managing Director

Jörg Ziel

Jörg Ziel
Director of QM and Technical Service

Digital signed von Jörg Ziel
DN: c=Germany, o=Heidolph Scientific Products GmbH, CN=Jörg Ziel, email=jorg.ziel@heidolph.de
I am the author of this document
Datum: 2026.09.21 10:46:09+02'00'
Foxit PDF Reader Version: 2024.3.0

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,
Heidolph Scientific Products GmbH
Walpersdorfer Straße 12
91126 Schwabach / Germany

Product: Laboratory magnetic stirrer

Model:

Hei-PLATE Mix'n'Heat Core
Hei-PLATE Mix'n'Heat Core+ /RS232
Hei-PLATE Mix'n'Heat Expert
Hei-PLATE Mix'n'Heat Ultimate
Hei-PLATE Mix 20 I
Radleys HP Core+

506-1xxxx-xx
506-2xxxx-xx
506-3xxxx-xx
506-4xxxx-xx
506-5xxxx-xx
RR 91213-xxx

Description:

Magnetid stirrer with or without heating plate and human machine interface

The object of the declaration described above is in conformity with the relevant
UK Statutory Instruments (and their amendments):

2008 No. 1597

The Supply of Machinery (Safety) Regulations

2008 2016 No. 1091

The Electromagnetic Compatibility Regulations 2016

2012 No. 3032

*The Restriction of the Use of Certain Hazardous Substances in
Electrical and Electronic Equipment Regulations 2012*

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP (UK) Ltd.
8. Northumberland Av.
London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH
Walpersdorfer Straße 12, 91126 Schwabach / Germany

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

© Heidolph Scientific Products GmbH

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Scientific Products GmbH PRODUCTS