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Operating instructions

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Hei-TORQUE Core RS232

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Seite 4 – 26

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Datum _____

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Firmenstempel _____

Translation of the original instructions

Page 30 – 52

Traduction de la notice originale

Page 56 – 78

Zertifikate/Certifications

Contents

Introduction

About this document.....	30
Typographic conventions.....	30
Copyright protection.....	30
Safety instructions in the official languages of the European Union	30

General notes

Directives applied, product certification.....	31
Copyright notice	31
California Residents	31
Residual risk	31
Intended use	31
Compliant use	31
Reasonably foreseeable misuse	32
Transportation	32
Storage	32
Acclimatization	32
Permissible ambient conditions	32
Energy saving potential and cost efficiency	32

Safety

General safety instructions.....	33
Electrical safety.....	33
Interfaces	33
Data security	33
Operational safety	34
Occupational safety	34
Electrostatic charging of liquid or powdered media.....	34
Personal protective equipment (PPE).....	34
Environmental protection	35
Biohazard	35
Other regulations.....	35

Device description

Overview overhead stirrer Hei-TORQUE Core RS232	36
---	----

Commissioning

Mounting the support rod	37
Positioning the device	37
Chuck for stirrer tools	38
Clamping/removing stirrer tools	38
Power supply	39
Switching the device on/off	39
Automatic restart	40

Operation

Operating modes.....	41
Continuous operation	41
Timer mode.....	41
Fast run mode.....	42
Remote mode.....	43
Setting the rotation speed.....	44
Overload protection	44

Troubleshooting

Troubleshooting.....	45
----------------------	----

Appendix

Technical specifications	46
Recommended speeds	47
Performance range	47
RS232 interface assignment.....	47
Interface commands	48
Scope of delivery.....	49
Device service.....	49
General cleaning instructions.....	49
Repairs - Return of equipment.....	50
Maintenance.....	50
Disposal.....	50
Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG).....	50
Warranty statement.....	51
Contact details	51
Certificate of decontamination.....	52

About this document

These operating instructions describe all functions and operation of overhead stirrers of the type Hei-TORQUE Core RS232. The operating instructions are an integral part of the delivery.

Typographic conventions

This document uses the following symbols, signal words, and highlights:

	Description
	Warning symbols in combination with a signal word indicate dangers: DANGER Indicates an immediate dangerous situation. Failure to respect the indications will result in death or serious injury. WARNING Indicates a potential danger. Failure to respect the indications will result in serious injuries. CAUTION Indicates a potential hazard. If not avoided, damage to property and minor to moderate injuries can occur.
	Mandatory signs are used to indicate important and useful information on handling a product. This information is used to ensure operational safety and to maintain the value of the product.
[...]	Square brackets identify labels of control elements on the device as well as labels and entries in software masks and user interfaces.
→	The arrow indicates specific instructions to be followed to ensure operational safety when handling the product.

Copyright protection

This document is protected by copyright and is intended for use by the purchaser of the product only.

No transfer to third parties, reproduction in any form, including excerpts, and by any means, as well as utilization and/or disclosure of the contents is permitted without the prior written consent of Heidolph Scientific Products GmbH. Any violation is subject to compensation for damage.

Safety instructions in the official languages of the European Union

A summary of all safety instructions in this document, in all official languages of the European Union, can be found in our Safety Guide for the Hei-TORQUE product group (Ref. 01-005-007-28). This document is available for download on our homepage in the most up-to-date version.

Directives applied, product certification

CE marking

The product complies with the following standards:

- Machinery Directive 2006/42/EE
- EMC Directive 2014/30/EU

Copyright notice

The software implemented in this device is protected by copyright laws. The rights holder is Heidolph Scientific Products GmbH, Walpersdorfer Straße 12, 91126 Schwabach/Germany. Any open source software components in this software are excluded from our copyright. Further information is available in the service area on our website www.heidolph.com.

Violations of copyright (e.g. unauthorized use or modification of the software) may result in civil claims (e.g. omission, compensation for damages) and/or entail criminal penalties within the meaning of the legislation of the destination country.

California Residents

Important information for California residents regarding Prop 65. Please visit www.P65Warnings.ca.gov for more information.

Residual risk

The device has been designed and manufactured in accordance with the state-of-the-art standards at the time of development and the recognized safety regulations. During mounting and use, as well as during maintenance, repair and cleaning work, there are nevertheless certain residual risks associated with the described device.

These are identified and described at the appropriate points in this document.

Intended use

The described device has been specially developed for the following tasks: Stirring, mixing, gassing and degassing, emulsifying and suspending.

The areas of application of the described device include chemical, biological and environmental-analytical laboratory and research applications, applications of basic research and comparable facilities.

Due to its design, the device in its delivery condition may only be used in analytical processes or in laboratory-like conditions in the food, cosmetics, and pharmaceutical industries as well as other comparable industries that manufacture products intended for consumption by humans or animals, or for use on humans or animals.

Any other use of this device is not considered as intended!

Compliant use

The compliance of each individual application must always be evaluated by the user. Any additional measures necessary to ensure compliance are always the responsibility of the user.

Reasonably foreseeable misuse

Additional measures may be necessary for use under conditions or for purposes deviating from the intended use, and/or specific guidelines and safety regulations must be observed. Corresponding requirements must be evaluated and implemented by the operator in each individual case.

Compliance with and implementation of all relevant directives and safety measures for the respective field of application is the sole responsibility of the operator.

All risks resulting from improper use are borne solely by the operator.

The device may only be operated by authorized and instructed personnel. Training and qualification of the operating personnel as well as ensuring that the device is handled responsibly are the sole responsibility of the operator!

Transportation

During transportation, avoid strong vibrations and mechanical stresses that can cause damage to the product. Keep the original packaging in a dry and protected place for later use.

Storage

Always store the device in its original packaging. To protect against damage and excessive material aging, store the device in an environment that is as dry, temperature-stable and dust-free as possible.

Acclimatization

After each transportation and after storage under critical climatic conditions (e.g. high temperature difference between inside and outside), allow the product to acclimatize at room temperature for a minimum of two hours to prevent possible damage from condensation before putting it into operation in the place of use. If necessary, extend the acclimatization phase if the temperature differences are very high.

Make all supply connections only after the product has been acclimatized!

Permissible ambient conditions

The device is designed for indoor use only. The device is **NOT** suitable for outdoor use! The device is **NOT** suitable for use in potentially explosive areas!

When used in corrosive atmospheres, the service life of the device may be reduced depending on the concentration, duration and frequency of exposure.

Energy saving potential and cost efficiency

Even in standby mode, small amounts of energy are used to supply individual device functions. However, especially over longer periods of nonuse, this resource consumption adds up and causes avoidable costs.

Therefore, always switch off the device when not in use, especially for long periods of time (devices with main switch: Main switch in position [O]; devices without main switch: Unplug) to effectively avoid unnecessary energy consumption and the resulting costs.

General safety instructions

Before commissioning and using the device, familiarize yourself with all the safety regulations and occupational safety guidelines applicable at the place of use and observe them at all times.

Only operate the device if it is in faultless technical condition. In particular, ensure that there is no visible damage on the device itself and, where applicable, on connected devices or the supply connections.

If there is missing or misleading information on the device or regarding occupational safety, contact the responsible safety specialist or our technical service.

Only use the device in accordance with the regulations for intended use.

Electrical safety

Before connecting the device to the power supply, ensure that the voltage indicated on the rating plate matches the specifications of the local power utility company.

Ensure that the power supply circuit provided is protected by means of a residual-current device (RCD).

Always use the power supply cord provided with the device.

Prior to use, check that the device and the power supply cord are free of visible damage.

Make sure that the power plug can be reached directly at any time, in order to disconnect the device from the power supply at any time without delay.

Have repairs and/or maintenance work on the device carried out exclusively by an authorized electrician or by the technical service department of Heidolph Scientific Products GmbH.

Always disconnect the device from the power supply system before carrying out any maintenance, cleaning or repair work.

Interfaces

Safely separate extra low voltage inputs and outputs rated 25 V AC or 60 V DC or above according to EN 61140, or use double or reinforced insulation according to EN 60730-1 or DIN 60950-1.

Use only shielded connection cables. Connect the shield to the connector housing. Always cover unused interfaces (front and back) to prevent liquids from entering the device.

Data security

The user is responsible for ensuring data security when transferring data between the described device and other devices.

Only use secure networks for the data transfer and avoid use of critical infrastructure.

Only use high-quality shielded data cables for the data transfer.

For data transfer via USB, an industrial standard USB hub should be preferably used to ensure the most stable connection possible.

Operational safety

Operate the device under a closed ventilated fume hood when working with potentially hazardous substances (see EN 14175 and DIN 12924).

Do not make any unauthorized changes or modifications to the device!

Only use genuine spare parts and accessories, or those expressly approved by the manufacturer!

Rectify malfunctions or faults on the device immediately.

Switch off and disconnect the device from the power supply, preventing reconnection, if it is not possible to eliminate the malfunction or rectify the fault immediately.

Observe all other applicable regulations such as laboratory and workplace guidelines, recognized safety technology rules and special local regulations.

Occupational safety

Always use the prescribed personal protective equipment (PPE) such as protective clothing, safety goggles, protective gloves, safety shoes, etc.

Do not operate any other devices in the immediate vicinity of the device ...

- which can generate electromagnetic fields in the frequency range between 9×10^3 Hz to 3×10^{11} Hz,
- which generate emission or radiation sources in the frequency range 3×10^{11} Hz to 3×10^{15} Hz (in the optical spectral range wavelengths from 1,000 μm to 0.1 μm),
- which generate ultrasonic or ionizing waves.

Do not process any substances that could release energy in an uncontrolled manner (e.g. self-ignition).

Do not process substances in which the energy input through mixing poses a danger.

Wipe off any liquid that may have spilled on the device immediately.

Before each transfer or transport (e.g., removal from the fume hood, transfer between laboratory benches, temporary storage in the laboratory area, etc.), always remove all attachments from the shaking platform of the device.

Electrostatic charging of liquid or powdered media

Liquid or powdered media as well as the laboratory vessels used can become dangerously electrostatically charged due to the shaking movement or when emptying or cleaning the vessel. The intensity of the charge generated depends on the properties of the medium, the flow rate, the process parameters as well as the size, geometry and material of the laboratory vessels (source: TRGS 727, Section 4).

The evaluation of a potential explosion hazard due to electrostatic charging and the definition of suitable safety measures are always the responsibility of the operator and are always based on the substances to be processed and the individual application.

Personal protective equipment (PPE)

The operating company must determine and provide the necessary PPE, depending on the respective application and the media and chemicals used.

The required instruction of the personnel is solely within the operating company's responsibility.

Environmental protection

When processing environmentally hazardous substances, take appropriate measures to avoid hazards to the environment.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Biohazard

When processing biohazardous substances, take appropriate measures to prevent hazards to persons and the environment, including:

- Instruction of the personnel regarding the necessary safety measures.
- Provision of personal protective equipment (PPE) and instruction of the personnel in its use.
- Marking of the device with the biohazard warning symbol.

The evaluation of appropriate measures such as the marking of a hazardous area, their implementation, and the training of the relevant personnel is the sole responsibility of the operator!

Other regulations

In addition to the notes and instructions in this document, observe all other applicable regulations such as laboratory and workplace guidelines, hazardous substances ordinances, recognized rules of safety engineering and occupational medicine as well as particular local regulations!

Noncompliance will invalidate any warranty claims against Heidolph Scientific Products GmbH.

The operating company is solely liable for all damage resulting from unauthorized changes or modifications to the device, from the use of unapproved or non-genuine spare parts and accessories, or from disregarding the safety instructions and hazard warnings or the manufacturer's instructions!

Overview overhead stirrer Hei-TORQUE Core RS232

1 Display

2 Push-and-rotary control, used to adjust the rotation speed and timer function


3 [Fast mode] button

4 [Timer] button

5 On/off switch

6 RS232 interface (rear)

7 Status indicator [Timer] - is displayed when the timer function is activated

8 Status display [automatic restart] - is displayed when the automatic restart function is activated

9 [External Control Incoming] indicator - flashes while the unit is receiving external control commands

10 Status indicator [remote control] - appears when the device is connected to a remote controller via the RS 232 interface

11 Overload indicator

12 Rotation indicator

13 Setpoint speed indicator


14 On time indicator

Mounting the support rod

The supplied support rod is used to fit the stirrer on a stand or a wall grid.

Insert the support rod in the receiving hole at the back of the stirrer so that the screwed-in grub screw presses on the flat area of the support rod and thus fixes it securely:

Turn the grub screw clockwise to fix the support rod.

Tighten the grub screw by hand only and with the help of the hex wrench supplied. Make sure that you do not overtighten the grub screw, as this damages the thread.

Turn the grub screw counterclockwise to loosen and dismantle the support rod. Then fix the fitted support rod to an adequately load-bearing wall grid or a stand with the help of suitable clamps (for suitable accessories, see www.heidolph.com).

Positioning the device

The described overhead stirrers can be mounted on a stand or on a wall grid. Any necessary mounting materials are available as accessories, see www.heidolph.com. Follow the safety instructions below:

WARNING

If assembled improperly, the device can also fall or tip over when unloaded!


If the overhead stirrer is mounted on a wall grid, ensure that the wall grid and all mounting elements have sufficient self-protection and load capacity!

When mounting the overhead stirrer on a stand, make sure that the setup is sufficiently stable.

- The stand may only be set up on a stable surface.
- When mounting on a stand, follow the specific setup instructions for the stand.
 - Always align the overhead stirrer so that the quick-action chuck vertically faces the notional centerline between the **long** stand base bars, see following figure:

After assembly of the overhead stirrer, check all screw connections / clamping screws for tight fit.

Chuck for stirrer tools

WARNING

Before mounting/dismantling the quick-action chuck, switch off the device and disconnect it from the power supply.

The chuck is fixed on the stirrer shaft by a grub screw so that it cannot rotate:

- Push the chuck onto the stirrer shaft so that the screwed-in grub screw engages in the assembly hole of the stirrer shaft and the chuck can thus be fixed securely (above figure, right).
- Use the hex wrench supplied to turn the grub screw clockwise, in order to fix the chuck.

Tighten the grub screw by hand only and with the help of the hex wrench supplied.

Make sure that you do not overtighten the grub screw, as this damages the thread.

Clamping/removing stirrer tools

WARNING

If necessary, stop the rotational movement and before clamping/removing stirrer tools, switch off the device to prevent it from starting up unexpectedly!

- Turn the chuck clockwise to open it until the required stirrer tool can be inserted.
- Insert the stirrer tool into the chuck. Make sure that the working height is as required!
- Hold the stirrer tool as close as possible to the middle, and close the chuck by carefully turning it counterclockwise until the stirrer tool is tightly clamped.

The quick-action chuck can be overturned by excessive force applied on opening! In this case, the quick-action chuck can no longer be closed normally and a latching noise can be heard on turning it in the closing direction.

Remedial action: Turn the quick-action chuck again in the opening direction, beyond the latching point.

The quick-action chuck can then be operated normally again.

Power supply

DANGER

Observe the instructions in section "General safety instructions" on page 33, in particular Electrical safety.

Power supply cord with locking system

Only use the supplied power supply cord with locking system!

Connect the power supply cord to the IEC-appliance inlet on the back of the device. Make sure that the locking system engages correctly.

Connect the power supply cord to a properly secured power socket.

Disconnecting the power supply cord

Disconnect the power supply cord from the mains socket-outlet.

Press in the locking lever on the cable coupling and disconnect the power supply cord from the IEC appliance plug on the back of the device.

Switching the device on/off

To switch on and off, use the main switch on the front of of the device..

Automatic restart

This function can be used to set the device so that the stirrer is restarted automatically after being switched off or after a power outage, and the speed is run up to the last set value.

CAUTION

Always activate or deactivate the automatic restart function as required before the actual process start!

Pay attention to the safety symbols on the device for signaling the functional status for all process sequences.

If the automatic restart function is active, secure the device with an information sign if necessary.

Proceed as follows to activate/deactivate the automatic restart function:

- Switch on the device and use the speed controller to set the rotational speed value to [213 rpm]: Value has white background field.
- Wait for approx. six (6) seconds until the value is accepted: Value is displayed without background field.
- Then use the speed controller to set the rotational speed values in chronological order to [214 rpm], [215 rpm], [216 rpm] and [217 rpm] and wait each time until the value is accepted.
- As soon as the last value [217 rpm] has been accepted, the status change is completed:
 - On activating the function, the [restart on] notice is displayed briefly.
 - In operation with active automatic restart, a small white caution symbol is shown continuously at the top edge of the screen.
 - On deactivating the function, the [restart off] notice is displayed briefly.
 - As soon as the automatic restart function is deactivated, the warning symbol at the top edge of the screen goes out.

On restarting after a power interruption, a [Rotation] warning triangle appears on the display briefly, the white caution symbol at the top edge of the screen then flashes until the rotational movement is stopped normally.

A power interruption that occurs during unsupervised operation can be identified by this flashing function.

Operating modes

The device can be used in continuous operation, timer controlled or remote controlled. Regardless of the selected operating mode, observe the following safety instructions:

WARNING

Risk of injury due to glass break and leaking medium

- Select the diameter of the laboratory vessel so that the stirrer tool turns inside the vessel without touching it.
- After clamping the stirrer tool and before the process starts, check the stirrer tool and the chuck for smooth running.
- Never start up stirrer tools or chucks that have an unbalance, and replace them immediately.

Risk of injury on rotating parts

- Never grip turning parts.
- Always wear close-fitting clothing when working with the device.
- Do not wear any accessories such as ties, scarves, shawls, pieces of jewelry, etc. when working with the device.
- When working with the device, tie back long hair or wear a cap or hair net.

Damage to property, loss of production

- Adjust the start-up behavior of the device to the specific properties of the sample to be processed to avoid overflowing/ejection of the sample or overloading of the stirrer tool (safety settings!).
- Follow the specific safety instructions for the substances to be processed (refer to the safety data sheets if applicable).
- Always wear the required individual personal protective equipment (PPE).

Continuous operation

In continuous operation, the rotational movement is started and stopped manually. Start and stop the rotational movement by pushing the push-and-rotary control.

Timer mode

In timer mode, a desired stirring time is preselected. Activate and set the timer as follows:

- Stop the rotational movement. The timer cannot be activated/set while the device is running.
- Press the [Timer] button: The [Timer] button lights up, the [Timer] setting menu opens.
- Use the speed controller to select the [Timer set] menu option.
- Press the speed controller to confirm the selection.

- Select the [hh], [mm] and [ss] options one after the other, and set a required stirring duration in hours, minutes and seconds (from one second up to 99 hours, 59 minutes and 59 seconds).
 - Confirm each individual value by pressing the speed controller.
 - Turn the speed controller to mark (select) and adjust the next value.
 - After confirming the value for seconds, the [ON] option is selected automatically: Press the speed controller to confirm and activate the timer.
- The timer is started on starting the rotational movement.

For physical reasons, the maximum speed can only be reached after a certain delay. This factor must be taken into consideration, in particular if, in timer mode, high speeds are to be reached in a very short time.

The effective delay until the required rotation speed is reached depends on the prevailing ambient conditions, on the stirrer tool used, and the substance processed.

- If the rotational movement is stopped before the timer has expired, it is stopped and continues running when the rotational movement restarts.
- If the timer time expires normally (no interruption due to manual stopping of the rotational movement), the rotational movement is stopped automatically and the timer is reset to the set time, e.g. 60 s.
- If the rotational movement is restarted after the timer expired normally, the timer is restarted with the set time.

While the timer is running, the remaining time is shown on the display to the nearest second. At the same time, the timer symbol appears on the left at the top of the screen.

To deactivate the timer, it must be switched to [OFF] via the Settings menu.

Fast run mode

With this function, the device can be set for a certain time to the maximum speed of 2,000 rpm with the press of a button:

- Press and hold the [Fast Run] button on the front of the device: The rotation speed is increased and maintained to 2,000 rpm.
- Let go of the [Fast mode] button to reduce the rotation speed back to the previously set value.

The fast mode function can be activated at any time and independently of the operating state. I.e., the stirrer tool starts, from a standstill if applicable, as soon as the [Fast mode] button is pressed.

Fast mode phases are not recorded in timer mode.

Remote mode

WARNING

Secure the device in remote mode with a clearly visible warning sign and, if necessary, take further ambient-specific protective measures that protect against damage to property and injuries in the event of unexpected/unmonitored starting of the device.

To activate the remote access of the device, a message of type "SetTargetSpeed" or "SetMotorState" has to be sent. Afterwards, the remote symbol appears on the display and remote commands are executed normally.

The transmission of the messages is made via ASCII strings. The protocol is active after the start. By a defined ADIP "Connect" command, you can switch to the ADIP mode.

RS232 interface parameters

- 9600 Baud
- Parity: none
- Data: 8 Bit
- Stop: 1 Bit

Command termination

All ASCII strings have to be finished with `\r\n` (0x0D 0x0A):

- ...`\r\n`

For further information on interface handling, please contact our technical service, see "Contact details" on page 51.

Device activation via PC

The device can be addressed via the integrated RS232 interface via an external control. The data interface enables activation of the device and recording of the process data in remote mode. A prerequisite for this is a suitable terminal program or software solution such as the Hei-PROCESSING Solutions from Heidolph Scientific Products GmbH

Device activation via server-based software

The device can be addressed via the integrated RS232 interface via a server-based software. If required, contact an authorized dealer or our technical service, see section "Contact details" on page 51.

Operation behavior in remote operation

In remote mode, neither the rotation speed on the device can be adjusted nor the timer function can be activated. The [Timer] and [Fast Run] buttons and the push-and-rotary control are not working.

If necessary, the remote operation can be terminated by pressing the push-and-rotary control by the operator or the device can be switched off via the on/off switch.

Setting the rotation speed

A rotational speed in the range of 20 – 2,000 rpm can be set. To set the required speed, do the following:

- Turn the speed controller clockwise to increase the (rotational) speed.
- Turn the speed controller counterclockwise to reduce the (rotational) speed.

When setting the required speed, monitor the set value on the display:

- In the [Rotation off] state the display value has a white background when changed. After six seconds, the new set value is adopted automatically and the background field is hidden.
- In the [Rotation on] state, the rotation speed is adjusted immediately according to the new setting.

The actual display is shown on the display as the actual value.

Overload protection

The device is equipped with an overload protection function: The performance indicator on the display visualizes the effective loading in the form of a progress bar.

In case of an overload (e.g. if speed is too high in highly viscous substances), the outermost right-hand segment of the performance indicator flashes. In this case, reduce the rotation speed until the performance indicator returns to normal mode.

In the event of persistent overload and in case of a blockage, the rotational movement is stopped. In this state, two restart attempts are started automatically. If the load does not reduce sufficiently, the rotational movement is switched off.

In this case, check the stirrer tool and the device for blockages and adjust the rotation speed to the substance to be processed if necessary. The rotational movement can then be started again normally.

Troubleshooting

The following table shows possible faults and solutions:

Error (message)	Possible cause, remedy
Display remains dark after switching on	- Power supply failure: check that the power plug is correctly connected - Electronic defect: contact technical service!
Rotation cannot be started, rotation stops unexpectedly, [Overload] error message	- Motor or power supply unit is overloaded: torque too high, stirrer blades blocked - Rotation stops completely after two restart attempts: Eliminate the cause of the blockage, switch on rotation again
[Overheating] error message	Motor/Electronics too hot: Overheat protection activated: Automatic shutdown. Allow the motor to cool and switch the device back on
Controls on the device do not function	Device is operating remotely, if necessary quit remote mode

In case of recurring errors please contact the responsible sales department or our technical service (see section "Contact details" on page 51).


Technical specifications

Model	Hei-TORQUE Core RS232
Dimensions (W × H × D)	70 × 282 × 195 mm
Weight	approx. 2.5 kg
Interfaces	1 × RS232 (Sub-D, 9 pole)
Stirrer shaft receptacle	Ø 0,5 – 10,5 mm
Control panel	TFT display, monochrome, 2 inch
Speed range	20 – 2,000 rpm
Speed accuracy	±1 %
Speed control	electronic
Drive	EC motor
Motor protection	Software monitoring with error message
Overload protection	Automatic cut-out
Permissible duty cycle	suitable for continuous operation
Volume H2O	25 l max.
Viscosity sample	max. 10,000
Torque	Nominal: 40 Ncm max. 55 Ncm
Protection class	IP42 (IEC 60529)
Protection class	I  (IEC 61140)
Overvoltage category	II
Degree of pollution	2
Sound pressure level (dB(A))	< 55 (based on IEC 61010)
Electromagnetic compatibility (EMC) acc. to EN 61326-1:2013	• Immunity: industrial electromagnetic environment • Emission: Class B, Group 1
Rated voltage	1/N/PE 100 – 240 V, AC 50/60 Hz
Power input	105 W
Operating and storage temperature	5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
Installation altitude	up to 2,000 m asl

Recommended speeds

Stirrer tool	Maximum rpm
Blade impeller: BR 10, BR 11, BR 12, BR 14	≤ 2,000
Blade impeller: BR 13	≤ 800
Half-moon impeller: H 18	≤ 800
Pitched-blade impeller: PR 39, PR 33	≤ 800
Pitched-blade impeller: PR 30, PR 31, PR 32	≤ 2,000
Radial-flow impeller: TR 20, TR 21	≤ 2,000
Anchor-type impeller: AR 19	≤ 800
VISCO JET® stirring system: VISCO JET® 60-120, VISCO JET® CRACK 80-120	≤ 500

Performance range

RS232 interface assignment

SUB-D9 con- nector	→	SUB-D9 socket
RXD (2)	→	(2) RXD
TXD (3)	→	(3) TXD
GND (5)	→	(5) GND

Interface commands

Command	Feedback	Meaning
r\x\n	RPM: xxxx\x\n	Display actual speed (rpm)
s\x\n	SET: xxxx\x\n	Display set speed (rpm)
m\x\n	NCM: xxxx\x\n	Display torque (Nm)
f\x\n	FLT: No Error!\x\n	No fault
	FLT: overload!\x\n	Motor overloaded
	FLT: Stopped Manually!\x\n	Stirrer stopped manually
	Note If the device is manually stopped in remote mode or an overload occurs, it automatically switches to manual mode. The device can then no longer be addressed by the remote controller. Feedback in both cases: FLT: Stopped Manually!\x\n]	
z\x\n	RPM: xxxx\x\n ZACT: HHMMSS\x\n	Display set rotation period and process times
t\x\n	TMP: xxx\x\n	Display PCB temperature
Rxxxx\x\n	Rxxxx\ nSET:xxxx\x\n	Apply set speed xxxx (rpm)
	Note If Interface Mode is set to [Safe] in the device menu and the motor has not been started manually, this command will switch the device to remote mode. To switch back to manual mode, the D\x\n command must be sent to the device or the device must be switched off and on again.	
	C\x\n	Clear Error\x\n
D\x\n	Direct control\x\n	Delete "overload" error status
M\x\n	M\x\n	Switch off connection:
		- Remote mode will be disabled - Device can again be operated manually
M\x\n	M\x\n	Identify device:
T\x\n	HT: 40\x\n	Display of the tripped device begins to flash
T\x\n	HT: 40\x\n	Display device type
F\x\n	F\x\n	Check connection between PC and stirrer (only if Interface Mode [Safe]. After 60 seconds without command input, the stirrer is switched off)
e\x\n	e\x\n	Set interface mode to [Expert]:
s\x\n	s\x\n	The device can be controlled remotely via RS232 command, while all local controls remain active
		Set interface mode to [Safe]:
		The device can either be controlled remotely via RS232 command or only manually.
Q\x\n	SN\x\n	View serial number

Scope of delivery

Item	Quantity	Product number
Hei-TORQUE Core RS232	1	501-60431-00
Support rod	1	22-02-14-01-41
Warranty registration	1	01-006-002-78
Power supply cord	1	country-dependent

Further information, particularly on the available accessories can be found on our website at www.heidolph.com!

Device service

When carrying out service work on the device (cleaning, maintenance, repair), observe the general instructions and safety information described in this section.

DANGER

Switch the device's main switch off and disconnect it from the power supply before carrying out maintenance work, cleaning, or repairs.

When cleaning, avoid the penetration of liquids.

Before replacing the fuses, switch off the device and disconnect the power supply cord.

Always replace the two device fuses in pairs with original manufacturer fuses. Further information on available accessories can be found on our website at www.heidolph.com.

After fuse replacement, check the device for a safe condition according to IEC 61010-1.

General cleaning instructions

Wipe all surfaces and the control panel with a damp cloth if necessary. Persistent contamination can be removed with mild soapy water.

CAUTION

Clean the device's surfaces with a soft, lint-free and only slightly moistened cloth.

Never use any aggressive or abrasive cleaning agents or aids.

Repairs - Return of equipment

Repairs to the device may only be carried out by authorized experts! Unauthorized repairs during the warranty period will result in the loss of the warranty claim. Regardless of the warranty claim, the owner is solely liable for damage caused by unauthorized repairs.

- In case of repair and before returning your device, contact our technical service at the following e-mail address:
 - service@heidolph.de.
- In your message, please provide us with the following information in addition to an error description:
 - Item number
 - Serial no.

The required data can be found on the rating plate of the device.

A service representative will contact you as soon as possible to agree on the next steps.

Include the completed certificate of decontamination with every device return, see "Certificate of decontamination" on page 52.

Maintenance

The device contains no user-serviceable components. If necessary (in the event of abnormal operating behavior such as excessive noise or heat generation, for example) contact our technical service.

Disposal

- When disposing of the device, observe the provisions of the WEEE Directive 2012/19/EU and its transposition into national law in the country of use.
- When disposing of portable batteries, observe the provisions of the European Battery Directive 2013/56/EU and their implementation in national law in the country of use.
- Check the device and all components for residues of substances that are hazardous to health, the environment and biohazardous before disposal.
- Properly remove and dispose of residues of substances that are hazardous to health, the environment and biohazardous!

Registration in accordance with the Electrical and Electronic Equipment Act (ElektroG)

Heidolph Scientific Products GmbH, headquartered in 91126 Schwabach/Germany, address Walpersdorfer Str. 12, is registered under the number DE 50705753 in the electronic waste equipment register of the foundation ear (www.ear-system.de).

Warranty statement

Heidolph Scientific Products GmbH provides a three-year warranty against material and manufacturing defects.

Glass and wear parts, transportation damage, and damage resulting from improper handling or non-intended use of the product are excluded from the warranty.

The warranty period for registered products begins on the date of purchase. Register the product with the enclosed warranty card or on our homepage www.heidolph.com.

For non-registered products, the warranty period begins with the date of the serial production (to be determined by the serial number).

In the event of material or manufacturing defects, the product will either be repaired or replaced free of charge within the warranty period.

Contact details

Heidolph Scientific Products GmbH

Technical service

Walpersdorfer Str. 12

D-91126 Schwabach/Germany

Email: service@heidolph.de

Representations

To find your local Heidolph distributor please visit www.heidolph.com

Certificate of decontamination

Enclose the certificate of decontamination, duly completed, with your device return. Submissions without a certificate of decontamination cannot be processed!

CERTIFICATE OF DECONTAMINATION IN CASE OF RETURNS

Please fill in the required fields.

Note: The sender must package the goods properly and appropriately for transport.

Heidolph Scientific Products GmbH

Walpersdorfer Straße 12
91126 Schwabach

Phone: +49 (0) 9122 9920-380

Fax: +49 (0) 9122 9920-19

E-Mail: service@heidolph.de

SENDER

Name _____

First name _____

Company/institution _____

Department _____

Workgroup _____

Address _____

ZC/City _____

Country _____

Phone _____

Email _____

DEVICE DETAILS

Article number _____

Serial no. _____

Ticket number _____

Reason for sending in _____

Has the device been cleaned, decontaminated/disinfected?

Yes

No

(Please mark as applicable)

If yes, which measures were carried out?

Does this device pose a risk to people and/or the environment due to the processing of substances that are hazardous to health, the environment and/or are biohazardous?

Yes

No

(Please mark as applicable)

If yes, with which substances did the device come into contact?

LEGALLY BINDING DECLARATION

The principal/consignor is aware that they are liable to the agent/consignee for losses or damage incurred due to incomplete and incorrect information.

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hereby declare, that the product (from serial number on 200067433) designated below is in compliance with the basic requirements of all applicable EU-directives stated below with regard to design, type of model sold and manufactured by us. This certificate will be invalid if the product is modified without the prior written consent and agreement of the manufacturer.

Hei-TORQUE Expert 100	501-61011-xx
Hei-TORQUE Expert 200	501-62011-xx
Hei-TORQUE Expert 400	501-64011-xx
Hei-TORQUE Ultimate 100	501-61021-xx
Hei-TORQUE Ultimate 200	501-62021-xx
Hei-TORQUE Ultimate 400	501-64021-xx
Hei-TORQUE Core	501-60411-xx
Hei-TORQUE Core RS232	501-60430-xx
Hei-TORQUE Safety Core Paket	501-60417-xx
Hei-TORQUE Safety Ultimate 100 Paket	501-61027-xx
Hei-TORQUE Safety Ultimate 200 Paket	501-62027-xx
Hei-TORQUE Safety Ultimate 400 Paket	501-64027-xx

Maschinenrichtlinie / Machinery Directive 2006/42/EG
EMV-Richtlinie / Electromagnetic Compatibility Directive 2014/30/EU
Delegierte (RoHS-) Richtlinie / Delegated (RoHS) Directive 2015/863/EU
Angewandte (harmonisierte) Normen / (Harmonized) Standards applied:
EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

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Jörg Ziel
Director of QM and Technical Service

Declaration of Conformity

In accordance with UK Government guidance

This declaration is issued under the sole responsibility of the manufacturer,

Heidolph Scientific Products GmbH
Walpersdorfer Straße 12
91126 Schwabach / Germany

Product: Laboratory stirring device

Model:

Hei TORQUE Expert 100	501-61011-xx
Hei TORQUE Expert 200	501-62011-xx
Hei TORQUE Expert 400	501-64011-xx
Hei TORQUE Ultimate 100	501-61021-xx
Hei TORQUE Ultimate 200	501-62021-xx
Hei TORQUE Ultimate 400	501-64021-xx
Hei TORQUE Core	501-60411-xx
Hei-TORQUE Core RS232	501-60430-xx

Description:

Overhead stirrer to work with stirring tools eg. impeller.

The object of the declaration described above is in conformity with the relevant UK Statutory Instruments (and their amendments):

2008 No. 1597
2008 2016 No. 1091
2012 No. 3032

The Supply of Machinery (Safety) Regulations
The Electromagnetic Compatibility Regulations 2016
The Restriction of the Use of Certain Hazardous Substances in
Electrical and Electronic Equipment Regulations 2012

and complies with the following technical standards :

EN ISO 12100:2010, EN 61326-1:2013, EN 61010-1:2010, EN IEC 63000:2018

UK Authorised Representative (for authorities only):

ProductIP (UK) Ltd.
8. Northumberland Av.
London WC2N 5BY

Signed for and on behalf of Heidolph Scientific Products GmbH
Walpersdorfer Straße 12, 91126 Schwabach / Germany

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With this declaration, we confirm (according to current knowledge) that all sold laboratory devices and components by Heidolph Scientific Products GmbH fulfill the requirements of the EU directive 2011/65/EU (RoHS) and the extension 2015/863. All devices are compatible with the requirements of the RoHS Directive for the following materials:

Max. 0,01% of weight in homogenous material for cadmium and max. 0,1% of the weight in homogeneous material for Lead, Mercury, Hexavalent chromium, Polybrominated biphenyls (PBB), Polybrominated diphenyl ethers (PBDE), Bis(2-ethylhexyl) phthalate (DEHP), Butyl benzyl phthalate (BBP), Dibutyl phthalate (DBP), Diisobutyl phthalate (DIBP).

In the case of individual assemblies, maximum concentrations maybe exceeded within the permissible exceptions of the Directive.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

China RoHS DECLARATION OF CONFORMITY

Heidolph Scientific Products GmbH has made reasonable efforts to ensure that hazardous materials and substances may not be used in its products.

In order to determine the concentration of hazardous substances in all homogeneous materials of the subassemblies, a "Product Conformity Assessment" (PCA) procedure was performed. As defined in GB/T 26572 the "Maximum Concentration Value" limits (MCV) apply to these restricted substances:

- Lead (Pb): 0.1%
- Mercury (Hg): 0.1%
- Cadmium (Cd): 0.01%
- Hexavalent chromium (Cr(VI)): 0.1%
- Polybrominated biphenyls (PBB): 0.1%
- Polybrominated diphenyl ether (PBDE): 0.1%

Environmental Friendly Use Period (EFUP)

EFUP defines the period in years during which the hazardous substances contained in electrical and electronic products will not leak or mutate under normal operating conditions. During normal use by the user such electrical and electronic products will not result in serious environmental pollution, cause serious bodily injury or damage to the user's assets.

The Environmental Friendly Use Period for Heidolph Scientific Products GmbH products is 25 years.

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This table is created according to SJ/T 11364-2014

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Note: Table applies to all products. Some of the components or parts listed above may not be part of the enclosed product.

- O: 表示该有毒有害物质在该部件所有均质材料中的含量均在GB/T 26572规定的限量要求以下。
O: Indicates that the above mentioned hazardous substance contained in all homogeneous materials of the part is below the required limit as defined in GB/T 26572.
- X: 表示该有毒有害物质至少在该部件某一均质材料中的含量超出GB/T 26572规定的限量要求。
X: Indicates that the above mentioned hazardous substance contained in at least one of the homogeneous materials of this part is above the required limit as defined in GB/T 26572.

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Apart from the disclosures in the above table, the subassemblies are not intentionally manufactured or formulated with lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), and polybrominated diphenyl ethers (PBDE).

Products manufactured by Heidolph Scientific Products GmbH may enter into further devices or can be used together with other appliances.

With these products and appliances in particular, Heidolph Scientific Products GmbH will not take responsibility for the EFUP of those products and appliances.

Schwabach, 04.September.2024

Andreas Hahn
Managing Director

Jörg Ziel
Quality Manager

© Heidolph Scientific Products GmbH

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Subject to change without notice. The printed version of this document is not regularly updated. The latest issue of this document can be found by visiting our homepage.

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MATERIAL CONTENT DECLARATION FOR Heidolph Scientific Products GmbH PRODUCTS