

Hei-FLOW Peristaltic Pumps Tubing

TUBING SPECIFICATIONS							
Tubing Type	Application	Features	Color and Material	Standard Compliant	Temp. Range	Sterilization	Restrictions
Tygon® standard 	General applications in the laboratory	Non-toxic, non-oxidizing Good resistance to acids, lyes and inorganic media Very low gas permeability, long service life	Transparent Thermoplastic soft PVC	FDA (21 CFR 177.2601) and GLP	-50 to +75 °C	✓ Can be autoclaved for 30 min at 1 bar and 120 °C (material may change color) or with ethylene oxide	Release of plasticizers possible
Tygon® 2001 for food 	Ideal for products with a high fat content	Extremely chemical-resistant, e.g. appropriate for the use of polar solvents Plasticizer- and oil-free Particularly long service life Translucent to aid visual inspection Outstanding flexibility	Transparent Thermoplastic	USP Class VI, FDA (21 CFR 177.2600) and GLP	-78 °C to +71 °C	✓ Can be autoclaved for 30 min at 1 bar, sterilized by irradiation or with ethylene oxide	-
Tygon® for hydrocarbons 	Especially for hydrocarbons, petroleum products and distillates	Ideal for petrol, kerosene, heating oil, cutting solutions and coolants on a glycol base Resistant to ozone and UV	Translucent Yellow Thermoplastic soft PVC	GLP	-40 to +75 °C	-	Not suitable for strong lyes and acids as well as food and pharmaceuticals.
PharMed® 	Ideal for medical, laboratory and research applications	High flexural fatigue strength Non-toxic, biocompatible Very low gas permeability Well suited for acids and lyes	Opaque Beige Polypropylene-based thermoplastic elastomer with plasticizers	USP Class VI, GLP, USP and Ph. Eur.	-51 to +135 °C	✓ Can be autoclaved or sterilized with ethylene oxide or sterilized by irradiation	Release of additives possible
Silicone 	For use in pharmacy and biology	Extremely smooth interior prevents bacterial growth Biocompatible, minimal adsorption and absorption Best flow properties, high temperature stability Absolutely inert, plasticizer-free	Translucent White Polydimethylsiloxane with silica earth and silicone additives, excellent contact pressure resistance	USP Class VI, meets GLP and NSF	-80 to +200 °C	✓ Can be autoclaved for 30 min at 1 bar or sterilized by irradiation	Unsuitable for concentrated solvents, oils, acids or diluted caustic soda, relatively high gas permeability
Viton® 	Excellent acid resistance - at high temperatures	Low gas permeability Resistant to solvents and corrosive media	Opaque Black Fluorocarbon rubber, thermoformed Viton B (67% fluorinated)	GLP	-30 to +205 °C	-	Limited service life

COMPATIBILITY						
Acids	Lyes	Solvents	Pressure	Vacuum	Viscous media	Sterile media
B Good	B Good	D Unsuitable	B Good	B Good	A Excellent	C Conditional
A Excellent	A Excellent	B Good	B Good	B Good	B Good	B Good
B Good	B Good	C Conditional	B Good	B Good	A Excellent	C Conditional
B Good	B Good	D Unsuitable	C Good	A Excellent	B Good	A Excellent
C Conditional	C Conditional	D Unsuitable	C Conditional	B Good	C Conditional	A Excellent
A Excellent	A Excellent	- Varying Tests recommended	B Good	B Good	B Good	C Conditional

Tubing Resistance Index

KEY

Tubing

P = PharMed®
 S = Silicone
 T = Tygon® standard
 TU = Tygon® for hydrocarbons
 TK = Tygon® 2001 for food
 V = Viton®

Resistance

A = excellent
 B = good
 C = conditional
 D = unsuitable
 - = not tested

Chemical	P	S	T	TU	TK	V
Butanol (butyl alcohol)	D	B	D	D	A	A
Butyl acetate	B	D	D	D	D	D
Butyric acid	B	D	D	C	D	-

Chemical	P	S	T	TU	TK	V
C Calcium oxide	A	A	A	A	A	-
Carbon disulfide	D	D	D	D	D	-
Carbon tetrachloride	D	D	D	D	D	A
Chlorine gas, wet	D	D	B	B	C	B
Chloroacetic acid, 20 % in W.	B	A	A	D	A	D
Chlorobenzene (Mono, Di, Tri)	D	D	D	D	C	A
Chloroform	D	D	D	D	C	A
Chlorobromomethane	B	D	D	D	-	A
Chromic acid, 20 % in W.	A	D	B	C	B	A
Chromic acid, 50 % in W.	C	D	C	D	-	-
Chlorosulfonic acid	D	D	D	D	D	D
Cupric Chloride 40% in W.	A	A	A	A	A	-
Cyclohexane	D	D	D	C	D	A
Cyclohexanone	D	D	D	D	C	D

Chemical	P	S	T	TU	TK	V
D Diesel	D	D	D	B	-	-
Dimethylformamide	B	B	D	D	A	D

Chemical	P	S	T	TU	TK	V
E Ethanol (ethyl alcohol)	A	B	D	B	A	A
Ether	C	D	D	C	D	-
Ethyl acetate	B	D	D	D	D	D
Ethylamine	D	C	D	D	B	-
Ethyl bromide	D	D	D	D	C	-
Ethyl chloride	C	D	D	D	D	A
Ethylene chlorohydrin	A	B	D	B	A	A
Ethylenedichloride	C	D	D	D	D	B
Ethylene glycol	A	A	A	A	A	A
Ethylene oxide	A	D	A	A	A	D

Chemical	P	S	T	TU	TK	V
F Fatty acids	C	B	B	C	C	C
Fluoroborate salts	A	-	A	A	A	-
Formaldehyde, 37 % in W.	D	C	D	D	C	D
Formic acid, 25 % in W.	A	A	A	C	A	D
Freon 11	A	A	A	A	-	-
Fruit juice	A	A	A	A	A	A

Chemical	P	S	T	TU	TK	V
G Glycerin	A	A	A	A	A	A

Chemical	P	S	T	TU	TK	V
H Hydrobromic acid, 20 – 50 %	D	D	A	A	A	A
Hydrochloric acid, 10 % in W.	A	D	A	A	A	A
Hydrochloric acid, 37 % in W.	B	D	A	D	A	B
Hydrofluoric acid, 10 % in W.	D	D	C	A	A	B
Hydrofluoric acid, 50 %	D	D	D	D	A	A

Chemical	P	S	T	TU	TK	V
Hydrogen cyanide	A	A	A	A	A	A
Hydrogen peroxide, 10 % in W.	A	A	A	A	A	A
Hydrogen peroxide, 90 % in W.	B	C	D	D	B	-
Hydroiodic acid	B	B	A	A	A	-
Hypochlorous acid, 25 % in W.	A	A	A	A	A	A

Chemical	P	S	T	TU	TK	V
I Iodine solutions	A	C	A	A	A	-
Iron II chloride 40 % in W.	A	A	A	A	A	B
Iron II sulfate 5 % in W.	A	A	A	A	A	A
Iron III chloride 43 % in W.	A	A	A	A	A	-
Iron III sulfate 5 % in W.	A	A	A	A	A	-

Chemical	P	S	T	TU	TK	V
K Kerosene	D	D	D	D	B	A
Ketones	D	D	D	D	C	-

Chemical	P	S	T	TU	TK	V
L Lactic acid, 10 % in W.	A	A	A	A	A	-
Lactic acid, 85 % in W.	B	D	D	D	-	-
Lead acetate, 35 % in W.	A	A	A	A	A	-

Chemical	P	S	T	TU	TK	V
M Magnesium chloride, 35 % in W.	A	A	A	A	A	A
Magnesium sulfate, 25 % in W.	A	A	A	A	A	-
Manganese salts	A	A	A	A	A	-
Mercury salts	A	A	A	A	A	-
Methane	A	-	A	A	A	A
Methanol	A	B	D	B	A	D
Methyl Ethyl Ketone	D	D	D	D	C	D
Monoethanolamine	C	D	D	D	D	D

Chemical	P	S	T	TU	TK	V
N Naphthalene	D	D	D	D	D	A
Nickel salts	A	A	A	A	A	-
Nitric acid, 10 % in W.	A	C	A	D	A	A
Nitric acid, 35 % in W.	A	D	A	D	A	A
Nitric acid, 68 – 71 % in W.	D	D	D	D	D	-
Nitrobenzene	D	D	D	D	C	-
Nitrous acid, 10 % in W.	A	B	A	C	A	-

Chemical	P	S	T	TU	TK	V
O Oils, animal	C	A	D	A	B	-
Oils, mineral	D	D	C	A	D	A
Oleic acid	C	B	D	B	D	B

Chemical	P	S	T	TU	TK	V
P Perchloroethylene	C	D	D	D	D	A
Perchloric acid, 67 % in W.	A	D	C	D	A	A
Phenol	A	D	D	C	A	-
Phosphoric acid 25 % in W.	A	D	A	A	A	A
Phthalic acid, 9 % in Alc.	A	B	D	C	B	-
Potassium carbonate, 55 % in W.	A	A	A	A	A	-
Potassium cyanide, 33 % in W.	A	A	A	A	-	-

Chemical	P	S	T	TU	TK	V
Potassium hydroxide, <10 % in W.	A	A	A	D	-	B
Potassium iodide, 56 % in W.	A	A	A	A	A	-
Propanol (propyl alcohol)	C	A	D	D	A	B
Pyridine	C	D	D	D	C	D

Chemical	P	S	T	TU	TK	V
S Silicone oils	C	D	B	A	B	A
Silver nitrate, 55 % in W.	A	A	A	A	A	A
Soap solutions	B	A	A	A	A	A
Sodium bicarbonate, 7 % in W.	A	A	A	A	A	A
Sodium bisulfate	A	-	A	A	A	-
Sodium borate	A	A	A	A	A	A
Sodium carbonate	A	A	A	A	A	B
Sodium dithionite	A	-	A	A	-	-
Sodium ferrocyanide	A	A	A	D	-	-
Sodium hydroxide, 10 – 15 % in W.	A	A	A	D	A	B
Sodium hydroxide, 30 – 40 % in W.	A	C	C	D	A	B
Sodium nitrate, 3.5 % in W.	A	A	A	A	A	-
Sodium sulfate, 3.6 % in W.	A	A	A	A	-	A
Sodium sulfide, 13 % in W.	A	A	A	A	A	-
Stearic acid, 5 % in Alc.	C	D	D	B	B	-
Sulfuric acid, 10 % in W.	A	A	A	B	A	A
Sulfuric acid, 30 % in W.	A	B	A	B	A	A
Sulfuric acid, 95 – 98 % in W.	D	D	D	D	C	A
Sulfurous acid	A	A	A	A	A	A

Chemical	P	S	T	TU	TK	V
T Tannic acid, 75 % in W.	B	A	B	A	D	-
Tartaric acid, 56 % in W.	A	A	A	A	A	A
Tin salts	A	A	A	A	A	-
Toluene	D	D	D	D	C	A
Trichloroacetic acid, 90 % in W.	B	D	A	D	A	C
Trichlorethylene	C	D	D	D	C	A
Trisodium phosphate	A	A	A	A	A	A
Turpentine	D	D	D	B	A	A

Chemical	P	S	T	TU	TK	V
U Urea, 20 % in W.	A	A	A	A	A	-
Uric acid	A	A	A	C	A	-

Chemical	P	S	T	TU	TK	V
X Xylene	D	D	D	D	C	B

Chemical	P	S	T	TU	TK	V
Z Zinc chloride, 80 % in W.	A	A	A	A	A	A

Please note:
 All information provided here is not guaranteed to be correct.
 Recommended testing of tubing prior to application use.