




America

CERTIFICATE

No. U8 069850 0013 Rev. 00

Holder of Certificate: **Heidolph Instruments GmbH & Co. KG**
Walpersdorfer Str. 12
91126 Schwabach
GERMANY

Certification Mark:



Product: **Laboratory Equipment**

This product was voluntarily tested to the relevant safety requirements referenced on this certificate. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited Certification body.

Test report no.: 713223865

Date, 2022-01-24

(Benedikt Pulver)



America

CERTIFICATE

No. U8 069850 0013 Rev. 00

Model(s): Hei-PLATE Mix'n'Heat Core
 Hei-PLATE Mix'n'Heat Core+
 Hei-PLATE Mix 20L
 Hei-PLATE Mix'n'Heat Expert
 Hei-PLATE Mix'n'Heat Ultimate

Brand Name: Heidolph

Tested according to: UL 61010-1:2012/R:2019-07
 UL 61010-2-010:2019
 UL 61010-2-051:2019
 CAN/CSA-C22.2 No. 61010-1:2012/A1:2018-11
 CSA C22.2 No. 61010-2-010:2019
 CSA C22.2 No. 61010-2-051:2019

Parameters:

Model:	Hei-PLATE Mix 20L		Hei-PLATE Mix'n'Heat Core Hei-PLATE Mix'n'Heat Core+ Hei-PLATE Mix'n'Heat Expert		Hei-PLATE Mix'n'Heat Ultimate	
Voltage	115V	230-240V	115V	230-240V	115V	230-240V
Frequency	50/60Hz	50/60Hz	50/60Hz	50/60Hz	50/60Hz	50/60Hz
Power	25W	25W	625W	825W	660W	860W