

2540 EPV SPECIFICATIONS

NEW ELECTRONIC BENCHTOP AUTOCLAVE

ITEMS TO STERILIZE



Glass



Tools



Wrapped
Tools



2540 EPV

PN

023210751



Chamber Volume (L)

23

Max Available
Height (in)

17

Chamber Dimensions
(in) (Øxd)

10 x 18

External Dimensions
(in) (w x h x d)

19 x 19.7 x 21.1

Power

120 V, 50-60 Hz

Connectivity

USB port (2), Wifi compatibility,
Ethernet port (1)

Number of Trays

5

Wastewater tank (L)



3-YEAR WARRANTY

PLUG & PLAY

Includes built-in reservoir.
No need to connect to
building utilities

T-CONNECT

Compatible with
T-Connect cloud software

COMPLETE CLOSED DOOR DRYING

Vacuum pump
completely dries load

TOUCHPAD

Touchpad interface for
easy operation

FAST HEAT-UP TIMES

Keeps chamber warm
between cycles to reduce
cycle times

