

HEIDOLPH NORTH AMERICA'S

DYNAMIC DUO



HEI-DRO CLEAN 165UC GLASSWASHER

Under Counter Laboratory Glassware Washer

3850ELP OR 3870ELP

Lab-Line Benchtop Autoclave



HEPA Filtered
Forced Air Drying



Built-in Water
Softener



Data logging



Direct Injection
Spindles



20 Preprogrammed
and 20 Custom
Cycles



Spray Arm and
Conductivity
Monitoring



No Utility
Connections
Required



2 Built-in Temperature
Probes for Complete
Liquid Sterilization



Create up to 16
Custom Cycles



4 Validated Cycles:
Glass, Plastic, Liquid A,
and Liquid B-Waste



Recycles Water



Free, built-in printer

Sterilize your samples with this clean combination for cost savings!

DISINFECTION AND STERILIZATION SYSTEM

A single package for all your laboratory disinfection and sterilization needs

Sterilization System Part number	Description
023250041	Hei-Dro Clean 165UC Washer** + 3850ELP Lab Line Benchtop Autoclave (Includes: Detergent, Neutralizer, Water Softener Salt, & Sterilizer Stand)
023250042	Hei-Dro Clean 165UC Washer** + 3870ELP Lab Line Benchtop Autoclave (Includes: Detergent, Neutralizer, Water Softener Salt, & Sterilizer Stand)

**All Hei-Dro Clean Washers come with upper standard rack & lower spindle rack

2 YEAR WARRANTY

HEI-DRO CLEAN 165UC WASHER



Chamber volume (L)	165
Chamber Dimensions (in) (w x h x d)	21.7 x 19.7 x 23.6
External Dimensions (in) (w x h x d)	23.6 x 25.6 x 33.3
Power	220 - 240 V, 3-phase NEMA L15-30 Ⓢ
Utilities	Cold Water, Hot Water, RO/DI Water drain
PN	023250040

INCLUDED WITH WASHER:

- Detergent: 023213791
- Neutralizer: 023213793
- Water Softener Salt



RO Water System *	Description
023213850	Hei-Dro Clean 165UC Washer RO water filtration system with 9.2 gal. expansion tank

*RO Water System sold separately

AUTOCLAVES



MODEL	3850ELP	3870ELP
Chamber volume (L)	65	85
Max available height (in)	12.1	12.1
Chamber Dimensions (in) (Øxd)	15 x 19.7	15 x 27.2
External Dimensions (in) (w x h x d)	28.3 x 21.3 x 30.1	28.3 x 21.3 x 37
Power	230V, 1-phase (60Hz) NEMA 6-20P Ⓢ	230V, 3-phase (60Hz) NEMA L15-30 Ⓢ
Utilities	No Plumbing Required	
Autoclave Stand	Included with Autoclave (023212609)	
PN	023210681	023210672

INCLUDED WITH AUTOCLAVE:

- Autoclave Stand: 023212609
- Chamber Brite: 023212803
- Printer Paper (Pkg. of 5): 023212200

