


Flexible Sample Formats

Genevac accessories are thoughtfully designed to flex with your laboratory's evolving workflow needs, helping to eliminate unnecessary steps and streamline your processes. Explore the systems at-a-glance:



EZ-2 4.0
(Standard, Plus, Elite)



HT Series 3i
(HT 6, HT 12)



miVac 2
(DNA, Duo, Quattro)

Sample Carriers / Holders	2 Sample Carriers	HT 6: Up to 6 Sample Carriers HT 12: Up to 12 Sample Carriers	Fixed-angle or Swing Rotors
Sample Formats	Microplates, Vials, Tubes, Bottles, Flasks	Microplates, Vials, Tubes, Bottles, Flasks	Tubes, Microplates, Vials, Flasks
Capacity	Up to 8 Microplates, 48 Tubes, 80 Vials	HT 6: Up to 24 microplates, 288 tubes, 480 vials HT 12: Up to 48 microplates, 576 tubes, 960 vials	DNA: 48 tubes or up to 6 microplates Duo: up to 6 shallow plates Quattro: up to 20 shallow or 8 deep-well plates
Solvent Compatibility	Standard: low BP inc, water, acetonitrile, DCM, methanol, TFA, organic solvents Plus: Medium BP inc. DMF Elite: High BP, Including DMSO *HCl and flammable solvent compatible models available	Entire range of solvents in previous column	Volatile organic solvents, water, medium boiling point solvents DNA: methanol, ethanol, water