

COMBINE CLIMATE-FRIENDLY
HEI-CHILL PELTIER

MADE IN
GERMANY

The average operational lifespan of 10 years is backed by a 3 year warranty and makes your purchase a worthwhile investment.

Thermostat Hei-CHILL Peltier

Chilling - Heating - Process integration

Environmentally friendly as no internal refrigerant is used.



The thermostat has a temperature range of $-10\text{ }^{\circ}\text{C}$ to $70\text{ }^{\circ}\text{C}$ and can also be used for smaller reactors with a working volume of 5 liters.

The power of 600 W is sufficient, also for chilling tasks when using rotary evaporators

Thanks to the plug-in system, the tubes are quickly connected to each other - thus, the connection of the device is ensured by simple steps.

The compact size of $27 \times 43 \times 43\text{ cm}$ enables the use on laboratory benches

Different interfaces: Ethernet, RS 232 and RS 485 allow the integration of the thermostat in processes



www.heidolph.com

For use in research purposes only – not for use in the medical area or pharmaceutical production!
California Residents: Important information for California residents regarding Prop 65.
Please visit www.P65Warnings.ca.gov for more information.

 **heidolph**
research made easy