

# PRECISE HEATING & STIRRING HEI-PLATE MIX 'N' HEAT MAGNETIC STIRRER

MADE IN  
GERMANY

HIGH-RESISTANT HEATING PLATE

REFRACTORY HOUSING

800 W

3-year warranty on  
all devices and an op-  
erational lifespan of  
more than 10 years

## Hei-PLATE Lift Package Ultimate

More comfort – More efficiency – More safety

### Interface for more precision

For the Pt1000 or the double sensor Pt1000 to measure the temperature exactly until inside the liquid

RS 232, LAN, USB A and USB B interfaces for datalogging as well as the remote control and monitoring via PC

USB C: Controllability of the electric Heidolph lab jack via magnetic stirrer and PC

### Well-designed features for a carefree working

Different heating modes (FAST, PRECISE, PRECISE+) and safety features handle your samples optimally

Ramp programming and creating profiles for a precise and quick reproducibility of your results



### Combination with e-Lift for more safety and comfort

Automatic lowering in case of an error or by timer setting (end of the process)

Ergonomic test set-up possible, regardless of the stand grids

Fine adjustment of the height via knob

Adjusted footprint

Robust mechanic made of V4A