

# SHAKING



## Hei-SHAKE Orbital Core

### THE ALL-ROUNDER UP TO 15 KG



#### Intuitive handling

Standby button in view to switch on and off the device and one knob to set values

Display and LED rings for an easy control and identification of the active functions

RS232 for remote control and recording data via software

Innovative body design for a secure handling and reliable trapping of the liquids from the device

#### Limitless possibilities

Broad rpm range (1 – 530)

Made for various vessels

Flexible attachments (i.e. changeable angle) for tool-free assembly

Platforms 1000 and 2000

#### High robustness and efficiency

Excellent load to footprint ratio (max. 15 kg load)

Up to 50% lower energy consumption when compared to other similar shakers

Low temperature rise on platform

For more information visit [our website](#)



Energy Efficiency



Safety




Sustainability



### MADE IN GERMANY

3-year warranty on all devices and an operational lifespan of more than 10 years

### RECOMMENDED ACCESSORIES



Perforated platform  
1000 and 2000



Separatory  
funnel attachment



Attachment for Falcon tubes,  
PCR tubes and Eppis



Multi-tier  
attachment