

SHAKING



Hei-SHAKE Wave Core

THE BIOSPECIALIST FOR
CONTINUOUS RINSING



Energy Efficiency



Safety



Sustainability



MADE IN GERMANY

3-year warranty on all devices and an
operational lifespan of more than 10 years

Smooth and efficient shaking

Broad rpm range (1 – 80)

Available in 10° angle to avoid foaming

Available in 5° angle to avoid foaming,
breakage and spillage

Intuitive handling

Standby button in view and one knob with
LED ring to set values and identify the
active functions

Display and timer for an easy control of the
process parameters

RS232 for remote control and data
recording via software

High robustness and safety

Excellent load to footprint ratio
(max. 15 kg load)

Up to 50% lower energy consumption
when compared to other similar shakers

Innovative body design for a
secure handling

Low temperature rise on platform

For more information visit [our website](#)

RECOMMENDED ACCESSORIES



Rubber Mat 1000 and 2000



Set for upgrade to platform 2000



Multi-tier attachment