

SHAKING



Hei-SHAKE Orbital Expert

THE ONE FOR AUTOMATISATION



Energy Efficiency*



Safety



Limitless possibilities

Broad rpm range (1 – 530)

Programmable alignment of the platform position for integration into automated systems

Made for various vessels and flexible attachments for tool-free assembly

Great ease of use

Touch display and timer for an easy control of the process parameters

Saving profiles and favorites for quick and easy repetition of routine applications

Wide range of interfaces (USB, RS232) for data recording and remote control

High robustness and efficiency

Excellent load to footprint ratio (max. 15 kg load)


Up to 50% lower energy consumption when compared to similar shakers

Low temperature rise on platform

For more information visit [our website](#)

*Claims based on internal tests. Not made by third party.

RECOMMENDED ACCESSORIES



Set for upgrade to platform 2000



Separatory funnel attachment



Attachment for test tubes, Falcon tubes, and Eppis



Multi-tier attachment