

SHAKING



Hei-SHAKE Reciprocating Core

THE MOST STURDY WITH
36 MM STROKE



Energy Efficiency



Safety



Sustainability



Limitless possibilities

Stroke of 36 mm & broad rpm range (1 – 300)

Made for various vessels

Flexible attachments (i.e. changeable angle) for tool-free assembly Platforms 1000 and 2000

Intuitive handling

Standby button in view and one knob with LED ring to set values and identify the active functions

Display and timer for an easy control of the process parameters

RS232 for remote control and data recording via software

Innovative body design for a secure handling

High robustness and efficiency

Excellent load to footprint ratio (max. 15 kg load)

Up to 50% lower energy consumption when compared to other similar shakers

Low temperature rise on platform

For more information visit [our website](#)



MADE IN GERMANY

3-year warranty on all devices and an operational lifespan of more than 10 years

RECOMMENDED ACCESSORIES



Flat rubber mat /
Set for upgrade 2000



Separatory
funnel attachment



Attch. for test tubes & Falcons
(centrifuge tubes) (45-90°)



Multi-tier
attachment