

SHAKING



Hei-SHAKE Rock Core

THE BIOSPECIALIST WITH SMOOTH GROOVE



Energy Efficiency



Safety



Sustainability



Smooth and efficient shaking

Broad rpm range (1 – 80)

Available in 10° angle to avoid foaming

Available in 5° angle to avoid foaming, breakage and spillage

Intuitive handling

Standby button in view and one knob with LED ring to set values and identify the active functions

Display and timer for an easy control of the process parameters

RS232 for remote control and data recording via software

High robustness and safety

Excellent load to footprint ratio (max. 15 kg load)

Up to 50% lower energy consumption when compared to other similar shakers

Innovative body design for a secure handling

Low temperature rise on platform

For more information visit [our website](#)



MADE IN GERMANY

3-year warranty on all devices and an operational lifespan of more than 10 years

RECOMMENDED ACCESSORIES



Set for upgrade to platform 2000



Multi-tier attachment