CONTINUOUS DISTILLATION: HEI-VOLUME AUTOMATIC MODULE

MADE IN GERMANY

The average operational lifespan of 10 years is backed by a 3 year warranty and makes your purchase a worthwhile investment.

Hei-VOLUME Distimatic Pro Benchtop with Hei-VAP

Saves time – increases throughput – reduces costs

Workload reductions at a glance:

No manual filling or draining

No handling of heavy glass flasks

Automatic backfeed of the heating bath

Draining of viscous residues

Processing of changing media by calibration-free sensors

Replacement of storage and collection vessels during the ongoing process

System shutdown including peripherals after the end of the process and in case of fault

Automatic cleaning program

Manual mode for processing even small quantities

Operation 24 hours a day, 7 days a week without monitoring

Enabling and maintaining constant conditions for higher throughput heidolp Educt sensor vessel Collector Control box with compressor

Process flow display with pictograms and colors for clear identification of process progress and multi-lingual user interface

