

CONTINUOUS DISTILLATION: HEI-VOLUME AUTOMATIC MODULE

MADE IN
GERMANY

The average operational lifespan of 10 years is backed by a 3 year warranty and makes your purchase a worthwhile investment.


Hei-VOLUME Distimatic Pro Industrial with Hei-VAP Industrial

Saves time – increases throughput – reduces costs

Workload reductions at a glance:

- No manual filling or draining
- No handling of heavy glass flasks
- Automatic backfeed of the heating bath
- Draining of viscous residues
- Processing of changing media by calibration-free sensors
- Replacement of storage and collection vessels during the ongoing process
- System shutdown including peripherals after the end of the process and in case of fault
- Automatic cleaning program
- Manual mode for processing even small quantities
- Operation 24 hours a day, 7 days a week without monitoring

Enabling and maintaining constant conditions for higher throughput



www.heidolph.com

For use in research purposes only – not for use in the medical area or pharmaceutical production!
California Residents: Important information for California residents regarding Prop 65.
Please visit www.P65Warnings.ca.gov for more information.

 **heidolph**
research made easy