

# EVAPORATE **THE OLD,** DISCOVER **THE NEW**

Trade in your old Heidolph rotary evaporator for the latest model and experience better performance, reliability, and efficiency. We'll handle disposal or recycling, so upgrading is seamless.

## Not sure if you need an upgrade?

Schedule an evaluation of your current equipment or easily modernize your Hei-VAP Generation 1 with an upgrade kit—featuring a new panel and all the latest features.



### HEIDOLPH'S LR4000

PTFE Vacuum Seal

Flask ejector presses off  
vapor tube

2 Independent power switches for  
control panel & heating bath

Adjustable bath kit to accommodate  
different flask sizes and/or bump trap  
with the most stability that supports  
the entire bath.



### THE HEI-VAP SERIES GENERATION 1

PTFE Vacuum Seal

Easy-Clip for releasing vacuum  
pressure and securing flask, no metal  
connection

NRTL certified for North American  
electrical safety compliance

Built-in extension cord to provide full  
control of your system from outside of  
the fume hood



### THE HEI-VAP SERIES GENERATION 2

Leading system tightness  
with FFKM vacuum seal to maximize  
life-span

One size that fits all applications from  
small scale up to a 5L flask

Dynamic AutoAccurate feature allows  
you to completely automate your  
evaporation process

\*Available on only Hei-VAP Expert Control and Hei-VAP  
Ultimate Control

Grease free system with zero glass on  
glass connections. An easier process  
to maintain your new screw threaded  
"stopcock"