

Hei-REACT Lab Reactors

THE SYSTEM WITH UNIQUE HANDLING COMFORT



Safety



CE Certified



Unique Handling Comfort

Quick and flexible handling of stirrer and lid via swiveling telescopic stand(s) with just one hand

Fast and intuitive set-up. Easily reproducible set-up alignment

Tool-free PEEK hose couplings and draindown-manifold kit for fast and mess free vessel exchange

Transportable tripod, safe storage on lab benches

Modular System

Broad selection of system components

Heidolph's modular system adapts to customer needs with preconfigured packages

Hei-REACT Base Configuration showcases unlimited modularity

Scale Without Compromise

Radley's Premium Glas Quality

Optimized geometry and height to width ratio (1.25 to 1) ensure easy scale up



MADE IN GERMANY

3-year warranty on all devices and an operational lifespan of more than 10 years

For more information visit [our website](#)

EXCERPT FROM CONFIGURATION LIST



Jacketed Reactor Vessels

Single or Vacuum Jacketed, Sizes ranging from 100ml to 5,000ml



5 Neck Lid

1xB19, 2xB24, 1xB29, 1xB34



Overhead Stirrers

RS232 interface for networking and external control using our laboratory automation software.



Thermoregulators

Open or closed systems for precise, dynamic temperature control over a wide temperature range.