


# UNLIMITED COMBINATIONS: HEI-SHAKE SHAKER

MADE IN  
GERMANY

The average operational lifespan of 10 years is backed by a 3 year warranty and makes your purchase a worthwhile investment.

## Unimax 1010

Shaking temperature controlled - Process integration being possible - Compact size



Interface RS 232 for integration into semi-automatic laboratory processes with Hei-PROCESS

Usable area of 290 x 258 mm provides sufficient space for different sample vessels

Time display of the timer to the nearest second provide full control of experiments

Automatic restart support the smooth laboratory work flow during unplanned interruptions

The overheat protection switches off the device in an emergency - this is important for temporarily unlimited continuous use

Integration in the modular incubator system 1000 ensures flexible applications

### ACCESSORIES



Modular incubator system 1000



Perforated platform 1000



Erlenmeyer clamps  
(25 - 2000 mL)



Tension roller  
attachment