

PRECISE HEATING & STIRRING HEI-PLATE MIX 'N' HEAT MAGNETIC STIRRER

MADE IN
GERMANY

HIGH-RESISTANT HEATING PLATE

REFRACTORY HOUSING

800 W

3-year warranty on
all devices and an op-
erational lifespan of
more than 10 years

Hei-PLATE Mix 'n' Heat Core⁺

Interface for more precision

For the Pt1000 to measure the temperature in the setup exactly

Two heating modes (FAST and PRE-CISE) handle your samples optimally according to the setup

Robust and durable for all challenges in the lab

Kera-Disk® plate – scratch-resistant, refractory, chemical-resistant and easy to clean

Refractory die-cast aluminum housing

New innovative body design offers additional protection against burns, secure handling and reliable trapping of liquids from the device

Intuitive handling for a quick start

Standby button in view to switch the device easily on and off

Display to set the set values and have an eye on the actual values at anytime

The locking function of the push-and-rotary controls protects against accidental changing of the values.

LED rings for the easy identification of the active functions and the residual heat



ACCESSORIES RECOMMENDATION



Silicone protective cover



Temperature sensor



Heating and bath attachments



High-performance condenser