

Hei-VAP Industrial with Hei-VOLUME Distimatic Pro

THE LARGE SCALE PACKAGE FOR VOLUMES UP TO 1,000 LITERS & MORE



CE Certified



Energy Efficiency



Safety



Sustainability



6 to 20 liter flask volume



MADE IN GERMANY

3-year warranty on all devices and an operational lifespan of more than 10 years

Automation to increase performance and reproducibility

Various modes for carrying out complex applications

Simple creation of profiles

Production of defined distillate quantities

Dispensing of viscous concentrates from the rotary flask

Manual mode for processing small quantities

Automation for continuous distillation

Calibration-free sensors monitor the process

Automatic filling and emptying (no need to change rotary flask). 24/7 operation - without supervision

No pressure fluctuations during the entire process cycle

System shutdown including peripherals at the end of the process

High level of user-friendliness


Recording and exporting of all process data possible

All components in contact with solvents are made from FDA-compliant materials


All process parameters at a glance thanks to 7" touch display

For more information visit [our website](#)

RECOMMENDED ACCESSORIES



Glassware A1 for higher performance distillation



Glassware R for high boiler distillation



Hei-VAC Valve Industrial



Base Cart



Base Cart Chiller 5000