

MULTIFARIOUS POSSIBILITIES: HEI-SHAKE SHAKERS

MADE IN
GERMANY

The average operational lifespan of 10 years is backed by a 3 year warranty and makes your purchase a worthwhile investment.

Multi Reax

Mixing powerfully – increasing sample throughput – saving time



The 3 mm shaking orbit in combination with a speed of 150 – 2,000 rpm guarantees the best mixing results

Two different adapters can be attached to the vortexer: for 2 tubes with 16 – 32 mm Ø or 6 tubes with 10 – 16 mm Ø

Comfortable inlets guarantee a secure hold when transporting the device

The speed can be adjusted via the speed controller. The parameters are shown precisely on the digital display

The timer allows the accomplishment of process steps efficiently