

3870ELP - 230V Technical datasheet



| | |
|---|----------------------------------|
| Height | 12.1" |
| Power | 230V, 3 phase (60Hz) NEMA L15-30 |
| Chamber Volume | 85 L |
| Chamber Dimensions | 15" x 27.2" |
| External Dimensions | 28.3" x 21.3" x 37" |
| Utilities | No plumbing required |
| Autoclave Stand | N/A |
| Autoclave Stand with Castors | 023212601 |
| Table - SS Top (36" x 48") Castors & Lower Shelf | 023212602 |
| Warranty | 3 Year |
| Certifications | UL, cUL, and ASME |