

Hei-Dro Clean 165UC Washer

Technical datasheet



| | |
|--|-------------------------------------|
| Recommended Flow Pressure | 44 psi (Min. 29 psi to Max. 73 psi) |
| Cold Water Temperature | Min. 5°C - Max. 20°C |
| Hot Water Temperature | Min. 45°C - 70°C |
| RO/DI Water Temperature | Min. 5°C - Max. 20°C |
| Connection Thread (with flask gasket) | G 3/4" (BSPP 3/4") |
| Length of Water Inlet Hose | 1.85 m/6.07' |
| Max Drain Water Temperature | 93°C/199°F |
| Max Pressure | 2.9 psi |
| Max Drain Pump Head Height from Floor Level | 1 m/3.28' |