

Hei-SHAKE Rock Core (5°) Technical datasheet



| | |
|--|---|
| Dimensions (W × H × D) | 356 × 199 × 403 mm |
| Weight | approx. 14 kg |
| Maximum permissible load | max. 15 kg |
| Drive | EC motor |
| Speed range | 0 – 80 rpm |
| Acoustic pressure | < 85 dB (A), in accordance with IEC 61010 |
| Rated voltage | 1/N/PE 100-240V AC 50/60 Hz |
| Protection Class | IP32 (acc. to EN 60529) |
| Power input | 90 W |
| Usable area | aprox. 300 × 270 mm |
| Movement | Rocking |
| Interfaces | 1 × RS232 (Sub-D, 9 pole) |
| Tilt Angle | 5° |
| Operating and Storage Temperature | 5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| Installed Altitude | up to 2,000 m asl |