

Model 5075 ELVCP - 230V  
**Technical datasheet**



<b>Height</b>	
<b>Power</b>	230V, 3 phase (60Hz) NEMA L15-30
<b>Chamber Volume</b>	160 L
<b>Chamber Dimensions</b>	19.7" x 29.5"
<b>External Dimensions</b>	34.3" x 42.9" x 30.3"
<b>Utilities</b>	Treated Water, Tap Water, Compressed Air, Floor Sink
<b>Wastewater Tank</b>	N/A
<b>Cooling Water Tank</b>	N/A
<b>Quiet Air Compressor for Models with Cooling</b>	023250023