


Multi Reax 115V

Technical datasheet



| | |
|------------------------------------|--|
| Dimensions (W × H × D) | 260 × 168 × 381 mm |
| Weight | approx. 10 kg |
| Maximum permissible load | 1.5 kg |
| Drive | Asynchronous motor |
| Speed range | 150 – 2,000 rpm; Indicator: digital |
| Protection class (EN 60529) | IP30 |
| Acoustic pressure | < 45 dB (A), in accordance with IEC 61010 |
| Rated voltage | 230 V / 50 Hz, 110 V / 60 Hz |
| Protection Class | I |
| Overvoltage category | n.a. |
| Degree of pollution | n.a. |
| Power input | 50 W |
| Operating temperature | 5 °C – 31 °C at up to 80 % rel. humidity 32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| Installation altitude | up to 2,000 m asl |
| Overheat protection | self-resetting |
| Speed control | n.a. |
| Usable area | n.a. |
| Movement | circular vibrating |
| Orbit | 3 mm |

| | |
|---|------|
| Fuse | n.a. |
| Permissible mains voltage fluctuations | n.a. |
| Bottle dimensions | n.a. |
| Interfaces | |