



## Reactor-Ready Pilot Core System (5 L to 35 L)

### Technical datasheet

<b>Frame Height</b>	1900mm with wheels (1800mm with feet)
<b>Frame Height - 015950404 200mm Height Reduction Kit</b>	1700mm with wheels (1600mm with feet)
<b>Frame Height - 015950403 300mm Height Reduction Kit</b>	1600mm with wheels (1500mm with feet)
<b>Frame Width</b>	540mm (approx. 800mm with hoses)
<b>Frame Depth</b>	650mm
<b>Vessel Sizes</b>	5L, 10L, 15L, 20L
<b>Thermofluid Temperature Range</b>	-70°C to +230°C
<b>Hose Size and Bore</b>	M30, 20mm ID
<b>Flange Size</b>	DN200 flange
<b>Lid Ports</b>	6 or 7 ports
<b>Max. Vacuum</b>	10mbar
<b>Max. Pressure</b>	0.5 Bar