

## System 27 Technical datasheet



<b>Height</b>	7 in/180 mm
<b>Reaction Positions</b>	24 (MTO format)
<b>Stirring Speed</b>	0 rpm to 1,100 rpm
<b>Colume Range</b>	0.5 ml to 7 ml
<b>Temperature Range, Cooling</b>	Ambient down to -7°C
<b>Temperature Range, Heating</b>	Ambient up to 150°C
<b>Time to Max. Heat Up</b>	30 min
<b>Temperature Stability</b>	+/- 0.5°C
<b>Base Diameter</b>	6.9 in/176 mm
<b>Top Diameter</b>	7.5 in/190 mm
<b>Well Diameter</b>	0.49 in/12.5 mm
<b>Operating Pressure</b>	3 psi