

Unimax 1010

Technical datasheet



| | |
|---|--|
| Dimensions (W × H × D) | 320 × 125 × 375 mm |
| Weight | approx. 8 kg |
| Maximum permissible load | 5 kg |
| Drive | EC motor |
| Speed range | 30 – 500 rpm |
| Acoustic pressure | < 50 dB (A), in accordance with IEC 61010 |
| Rated voltage | 100 – 240 V, 50/60 Hz |
| Overvoltage category | n.a. |
| Degree of pollution | n.a. |
| Power input | 50 W |
| Operating temperature | 5 °C – 31 °C up to 80 % rel. humidity 32 °C – 40 °C up to 50 % rel. humidity (decreasing linearly) |
| Installation altitude | up to 2,000 m asl |
| Overheat protection | n.a. |
| Speed control | n.a. |
| Usable area | 290 × 258 mm |
| Movement | orbital |
| Orbit | 10 mm |
| Fuse | M 1,25 A (230 V) or T 2,0 A (115 V) |
| Permissible mains voltage fluctuations | n.a. |
| Bottle dimensions | n.a. |

Interfaces

1 × RS232 (Sub-D, 9 pole)

Protection Class