

Technical datasheet



| | |
|---|----------------------|
| Height | 5.25" |
| Power | 120V |
| Chamber Volume | 7.5L |
| Chamber Dimensions | 7" x 12" |
| External Dimensions | 17" x 12.5" x 16" |
| Utilities | No Plumbing Required |
| Autoclave Stand | N/A |
| Autoclave Stand with Castors | N/A |
| Table - SS Top (36" x 48") Castors & Lower Shelf | N/A |
| Warranty | 3 Year |
| Certifications | UL, cUL, and ASME |