

## Technical datasheet



<b>Temperature range</b>	Ambient to 150°C
<b>Max. Pressure</b>	3 psi
<b>Reaction Positions</b>	24 (MTP format)
<b>Stirring Speed</b>	0 rpm to 1,100 rpm
<b>Temperature Accuracy</b>	Better than 0.5°C
<b>Inert Gas</b>	Yes
<b>Reflux</b>	Yes
<b>Time to Max. Heat Up</b>	30 min
<b>Well Diameter</b>	16 mm
<b>Working Volume</b>	0.5 ml - 7.0 ml