

## Hei-Dro Clean 165UC Washer

### Technical datasheet



<b>Recommended Flow Pressure</b>	44 psi (Min. 29 psi to Max. 73 psi)
<b>Cold Water Temperature</b>	Min. 5°C - Max. 20°C
<b>Hot Water Temperature</b>	Min. 45°C - 70°C
<b>RO/DI Water Temperature</b>	Min. 5°C - Max. 20°C
<b>Connection Thread (with flask gasket)</b>	G 3/4" (BSPP 3/4")
<b>Length of Water Inlet Hose</b>	1.85 m/6.07'
<b>Max Drain Water Temperature</b>	93°C/199°F
<b>Max Pressure</b>	2.9 psi
<b>Max Drain Pump Head Height from Floor Level</b>	1 m/3.28'