

Hei-SHAKE Reciprocating Core 2000

Technical datasheet



Dimensions (W × H × D)	356 x 129 x 403 mm
Weight	approx. 14 kg
Maximum permissible load	15 kg
Drive	EC Motor
Speed range	1 - 300 rpm
Acoustic pressure	< 85 dB (A), in accordance with IEC 61010
Rated voltage	1/N/PE 100-240 V, AC 50/60 Hz
Overvoltage category	N/A
Degree of pollution	N/A
Power input	90 W
Operating temperature	5 °C – 31 °C at up to 80 % rel. humidity
Installation altitude	up to 2,000 m asl
Overheat protection	N/A
Speed control	N/A
Usable area	(W x D) 390 × 340 mm
Movement	reciprocating (horizontal shaking motion)
Orbit	N/A
Fuse	N/A
Permissible mains voltage fluctuations	N/A
Bottle dimensions	N/A
Interfaces	N/A

Protection Class

IP32 (acc. to EN 60529)