

# Hei-VOLUME Distimatic Pro 24/7 Industrial con adaptador para vidrios de la serie A

## Technical datasheet



<b>Dimensions (W × H × D)</b>	186 × 429 × 521 mm, including control panel
<b>Weight</b>	approx. 17 kg
<b>Protection class (EN 60529)</b>	IP42
<b>Rated voltage</b>	100 – 240 V AC; 50/60 Hz
<b>Connection</b>	L+N+PE
<b>Protection class</b>	I
<b>Overvoltage category</b>	II
<b>Degree of pollution</b>	2
<b>Power input</b>	1,500 W max
<b>Operating temperature</b>	5 °C – 31 °C at up to 80 % rel. humidity   32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly)
<b>Installation altitude</b>	up to 2,000 m asl
<b>Fuse</b>	15 A