

## Radleys Tornado Overhead Stirring System Technical datasheet



|                                      |                       |
|--------------------------------------|-----------------------|
| <b>Weight</b>                        | 17.6 lb/8 kg          |
| <b>Height</b>                        | 23.1 in/587 mm        |
| <b>Temperature range</b>             | Ambient up to 180°C   |
| <b>Stirring Speed</b>                | 0 rpm to 1,000 rpm    |
| <b>Inert Gas</b>                     | Yes                   |
| <b>Reflux</b>                        | Yes                   |
| <b>Volume Range</b>                  | 20 ml to 150 ml       |
| <b>Gear Selections</b>               | 1:1 and 2:1           |
| <b>Stirring Seal Relief Pressure</b> | 1 bar                 |
| <b>Maximum Viscosity</b>             | 10,000 cps at 500 rpm |