



Reactor-Ready Pilot Core System

Technical datasheet

| | |
|--|---------------------------------------|
| Frame Height | 1900mm with wheels (1800mm with feet) |
| Frame Height - 015950404 200mm Height Reduction Kit | 1700mm with wheels (1600mm with feet) |
| Frame Height - 015950403 300mm Height Reduction Kit | 1600mm with wheels (1500mm with feet) |
| Frame Width | 540mm (approx. 800mm with hoses) |
| Frame Depth | 650mm |
| Vessel Sizes | 5L, 10L, 15L, 20L |
| Thermofluid Temperature Range | -70°C to +230°C |
| Hose Size and Bore | M30, 20mm ID |
| Flange Size | DN200 flange |
| Lid Ports | 6 or 7 ports |
| Max. Vacuum | 10mbar |
| Max. Pressure | 0.5 Bar |