

## Technical datasheet



<b>Hauteur</b>	5.25"
<b>Puissance</b>	120V
<b>Chamber Volume</b>	7.5L
<b>Chamber Dimensions</b>	7" x 12"
<b>External Dimensions</b>	17" x 12.5" x 16"
<b>Utilities</b>	No Plumbing Required
<b>Autoclave Stand</b>	N/A
<b>Autoclave Stand with Castors</b>	N/A
<b>Table - SS Top (36" x 48") Castors &amp; Lower Shelf</b>	N/A
<b>Warranty</b>	3 Year
<b>Certifications</b>	UL, cUL, and ASME