

Technical datasheet



| | |
|-----------------------------|--------------------|
| Temperature range | Ambient to 150°C |
| Max. Pressure | 3 psi |
| Reaction Positions | 24 (MTP format) |
| Stirring Speed | 0 rpm to 1,100 rpm |
| Temperature Accuracy | Better than 0.5°C |
| Inert Gas | Yes |
| Reflux | Yes |
| Time to Max. Heat Up | 30 min |
| Well Diameter | 16 mm |
| Working Volume | 0.5 ml - 7.0 ml |