

# Hei-CHILL 35

## Technical datasheet



<b>Dimensions (l × h × p)</b>	22.8x26.0x32.3 (in)
<b>Cooling</b>	Power at -10°C: 2   at -20°C: 1   at 0°C: 2.5   at 10°C: 3.5   at 15°C: 4
<b>Alimentation électrique</b>	208V 2~ 60 Hz
<b>Bypass for adjusting flow pressure</b>	x
<b>Connections all G thread male fittings</b>	1/2"
<b>Potential free contact for alarm</b>	Optional
<b>Pump flow max.</b>	30
<b>Pump max. (bar)</b>	4
<b>Remote Controller incl. 40m cable</b>	not possible
<b>RS232</b>	x
<b>Tank size</b>	15
<b>Temperature range</b>	15 to -20°C
<b>Thermofluid Water/Monoethylene Glycol 50:50</b>	x