

Hei-CHILL 35

Technical datasheet



Dimensions (l × h × p)	22.8x26.0x32.3 (in)
Cooling	Power at -10°C: 2 at -20°C: 1 at 0°C: 2.5 at 10°C: 3.5 at 15°C: 4
Alimentation électrique	208V 2~ 60 Hz
Bypass for adjusting flow pressure	x
Connections all G thread male fittings	1/2"
Potential free contact for alarm	Optional
Pump flow max.	30
Pump max. (bar)	4
Remote Controller incl. 40m cable	not possible
RS232	x
Tank size	15
Temperature range	15 to -20°C
Thermofluid Water/Monoethylene Glycol 50:50	x