

Hei-CHILL 5000 Base Cart

Technical datasheet



| | |
|--|--|
| Dimensions (l × h × p) | External dimensions [WxDxH] 1100 x 730 x 730 mm Installation dimensions [WxD] 1100 x 950 mm |
| Poids | 145 kg |
| Classe de protection (EN 60529) | IP 20 |
| Niveau de pression acoustique | 61 dB(A) [distance 1m] |
| Type de raccordement | Electrical: 208 V 3~/N/PE 50 Hz |
| Précision de réglage | at -10°C: +/- 1,5 K |
| Working temperature range | -10°C to +40°C |
| Display resolution | 0,1 K |
| Cooling capacity | at 20°C: 5000 W at 10°C: 3500 W at 0°C: 2500 W at -10°C: 1100 W |
| Circulation pump | Flow rate [max.] 40 l/min Flow pressure [max.] 3,0 bar |
| Hose connections | 3/4" female thread |
| Coolant tank | 28 l |
| Permitted ambient temperature range | 5°C to 40°C |
| Max. temperature of coolant | 60°C |
| Current [max.] | 16,0 A |
| Class. acc. DIN 12876 | I / NFL |
| Max. operating pressure | 28 bar |

| | |
|---------------------------|-----------------|
| Refrigerant | R513a (GWP 631) |
| CO2 equivalent (t) | 0,57 t |
| Alimentation | 4000 W |
| Filling quantity | 0,900 kg |