

Radleys Reactor-Ready Pilot Core System (5 L to 50 L)

Technical datasheet



Frame Height	1900mm with wheels (1800mm with feet)
Frame Height - 015950404 200mm Height Reduction Kit	1700mm with wheels (1600mm with feet)
Frame Height - 015950403 300mm Height Reduction Kit	1600mm with wheels (1500mm with feet)
Frame Width	540mm (approx. 800mm with hoses)
Frame Depth	650mm
Vessel Sizes	5L, 10L, 15L, 20L
Thermofluid Temperature Range	-70°C to +230°C
Hose Size and Bore	M30, 20mm ID
Flange Size	DN200 flange
Lid Ports	6 or 7 ports
Max. Vacuum	10mbar
Max. Pressure	0.5 Bar