

## Tornado Overhead Stirring System Technical datasheet



<b>Poids</b>	17.6 lb/8 kg
<b>Hauteur</b>	23.1 in/587 kg
<b>Temperature range</b>	Ambient up to 180°C
<b>Stirring Speed</b>	0 rpm to 1,000 rpm
<b>Inert Gas</b>	Yes
<b>Reflux</b>	Yes
<b>Volume Range</b>	20 ml to 150 ml
<b>Gear Selections</b>	1:1 and 2:1
<b>Stirring Seal Relief Pressure</b>	1 bar
<b>Maximum Viscosity</b>	10,000 cps at 500 rpm