

Reax top Technical datasheet



| | |
|---|--|
| Dimensions (W × H × D) | 134 x 105 x 172 mm |
| Weight | 2.8 kg |
| Maximum permissible load | n.a. |
| Drive | shaded pole motor |
| Speed range | 100 – 2,500 rpm; 200 – 3,200 rpm |
| Protection class (EN 60529) | IP 22 |
| Acoustic pressure | < 50 up to 2,500 rpm; < 58 up to 3,200 rpm |
| Rated voltage | 230 V 50 Hz; 230 V 60 Hz; 115 V 60 Hz |
| Protection class | 1 |
| Overvoltage category | n.a. |
| Degree of pollution | n.a. |
| Power input | 51 W |
| Operating temperature | n.a. |
| Installation altitude | n.a. |
| Overheat protection | self-resetting |
| Speed control | scale with reference points |
| Usable area | n.a. |
| Movement | circular vibrating |
| Orbit | 5 mm |
| Fuse | n.a. |
| Permissible mains voltage fluctuations | n.a. |

Bottle dimensions

n.a.