

## Technical datasheet



|  |   |
|--|---|
| <b>Dimensions (W × H × D)</b>            | 356 × 199 × 403 mm  |
| <b>Weight</b>                            | approx. 14 kg   |
| <b>Maximum permissible load</b>          | max. 15 kg  |
| <b>Drive</b>                             | EC motor  |
| <b>Speed range</b>                       | 0 – 80 rpm  |
| <b>Acoustic pressure</b>                 | < 85 dB (A), in accordance with IEC 61010   |
| <b>Rated voltage</b>                     | 1/N/PE 100-240V AC 50/60 Hz   |
| <b>Protection Class</b>                  | IP32 (acc. to EN 60529)   |
| <b>Power input</b>                       | 90 W  |
| <b>Usable area</b>                       | aprox. 300 × 270 mm   |
| <b>Movement</b>                          | Rocking   |
| <b>Interfaces</b>                        | 1 × RS232 (Sub-D, 9 pole)   |
| <b>Tilt Angle</b>                        | 5°  |
| <b>Operating and Storage Temperature</b> | 5 °C – 31 °C at up to 80 % rel. humidity<br>32 °C – 40 °C at up to 50 % rel. humidity (decreasing linearly) |
| <b>Installed Altitude</b>                | up to 2,000 m asl   |